## IC10/20

Rev. 21 / 02.2020

<b>INTERNATIONAL SYMBOL USAGE</b> You may see one or more of these symbols on the labeling/packaging of this product: <b>Key guide to symbols (Guía de símbolos, Guia de símbolo, Guide des symbols, Guida ai simboli, Sembollerin aina rehben).</b>	
	Atención! Superficie caliente / Caution! Hot surface / Atenção! Superficie quente / Attention! Surface chaude / Attenzione! Superficie calda / Dikkat! Sıcak yüzey
	Corriente continua / Direct current / Corrente contínua / Courant continu / Corrente continua / Doğru akım
	Atención / Caution / Atenção / Attention / Attenzione / Uyarı
	Código de lote / Batch code / Código do lote / Code du lot / Codice di lotto / Lot numarası
	Humedad Relativa de funcionamiento / Operating Relative Humidity / Umidade relativa de funcionamento / Humidité relative de fonctionnement / Umidità relativa di funzionamento / Bağıl alıřma nemi
	Temperatura ambiente de funcionamiento / Operating temperature / Temperatura ambiente de funcionamento / Température de fonctionnement / Temperatura ambiente di funzionamento / Ortam alıřma sıcaklıđı
	Mantener alejado de la luz solar / Keep away from sunlight / Manter afastado da luz solar / Garder à l'abri de la lumière du soleil / Mantener lontano dalla luce solare / Güneř ışıđından uzak tutun
	Recolección independiente para el desecho de equipos eléctricos y electrónicos / Separate collection for waste of electrical and electronic equipment / Colheita independente para a eliminaão de equipamentos elétricos e eletrônicos / Collecte sélective pour les équipements électriques et électroniques / Raccolta differenziata per lo smaltimento dei dispositivi elettronici ed elettronici / Elektrikli ve elektronik cihaz atıkları için ayrı toplama

## ES Incubadora Dual

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## EN Dual Incubator

for Biological Indicators and Hygiene Monitoring Systems

### Product description

Bionova® IC10/20 incubator has been designed for a wide range of products. It has a heating block with:  
-26 positions for conventional and Rapid Readout Colorimetric Self-contained Biological Indicators.  
-10 positions for Self-contained Ampoules, Culture Medium and Protein Pen Detection System.  
-A special hole for external thermometer.  
-An ampoule crusher for IBs.  
The incubator is pre-set at 37 or 60±2 °C. These two temperatures cannot be modified.

### Operating conditions

Bionova® IC10/20 incubator works at room temperature of 10-30 °C, Relative Humidity of 30-80% and a supply voltage of 12 V Direct Current. Altitude 3500 Meters (máx).

### Power supply specifications:

-Input parameters: 110/220 V, 0,2 A, 50/60 Hz.  
-Output parameters: 12 V, 4 A.

### WARNING

Do not pour any liquid inside.  
Do not immerse into any liquid.  
The incubator's metal surfaces may be hot during its functioning. Avoid contact with them.  
Only for indoor use.  
The incubator should be sent to the manufacturer for repair.  
Make sure that the incubator is connected to a properly rated power cord.  
Do not use this product in a manner not specified by Terragene® S.A., otherwise the protection provided by the product might be affected.  
Do not replace the power supply cord for other not provided within the product.  
Only personnel authorized by Terragene® S.A can access the inside of the auto-reader and its components. Parts or components inside the machine should not be manipulated by the user.

### Directions for use

- Place the auto-reader on a firm surface, free from vibrations, away from sunlight, currents of hot or cold air, chemicals and corrosive or flammable substances. Do not place the equipment in a way that disconnection of the plug from the power supply could be difficult. Leave a gap of at least 10 cm from the wall. **NOTE:** The proper switching-on of the device may be verified by the blinking of the temperature blue light.
- Turn on the device by first connecting the power supply AC outlet to the power grid and then connecting the other end of the power supply (plug) to the back of the unit.
- Select the incubation temperature using the switch button located at the base of the incubator. During heating, the blue light will remain blinking until it reaches the selected temperature.  
**IMPORTANT:** Check that the incubation temperature is the correct one for the product to be incubated.
- Incubate at 37 or 60±2 °C.  
Color change of the growth indicating medium (in conventional and Rapid Readout Colorimetric Self-contained Biological Indicators, Self-contained Ampoules or Culture Medium) shows a failure in the sterilization process.  
A color change in the Protein Pen Detection System indicates the presence of detectable protein and/or allergens levels.  
Used IBs or Protein Pen Detection Systems cannot be reused.  
**NOTE:** For the incubation and detailed analysis of the results, see specifications contained in each product's directions for use.

### Calibration

Bionova® IC10/20 incubator allows to calibrate the incubation temperature with the Bionova® TB-IC1020 thermometer. It is advised that this procedure be carried out only by qualified technicians. To perform the calibration, follow this steps:  
1. Turn on the Incubator and select the incubation temperature (37 °C or 60 °C). Wait until the selected temperature stabilizes; the blue light will stop blinking.  
2. Locate the calibration hole at the bottom of the incubator:

- Insert a precision flat-head screwdriver of 1,8 mm in the hole.
- To rise the incubation temperature, slightly rotate the screwdriver in a clockwise direction. Wait 2 minutes, after the blue light becomes stable again, and verify the temperature increase.
- To lower the incubation temperature, rotate the screwdriver in a counterclockwise direction. Wait 2 minutes, after the blue light becomes stable again, and verify the temperature variation.  
**NOTE:** During the calibration process, the light might become intermittent.  
**IMPORTANT:** Due to the high sensitive of the calibration hole, it's recommended to verify the temperature increase every time the screwdriver is rotated ¼ of a turn.

### Cleaning

Unplug the device from the power supply. If the device is hot, wait until it has cooled down before handling. Clean the outside of the device with a damp cloth and a small amount of detergent. Do not allow any liquid to run inside the device during its cleaning. Do not use abrasive or corrosive cleaners or disinfectants.

### Warranty

Terragene® S.A. guarantees both the quality of the product material components and the quality of its manufacturing process. Should any material or manufacturing faults be detected within the warranty period (1 year from the time of purchase), the only obligation of Terragene® S.A. will be product repairing or substitution.

### Limitation of liability

Terragene® S.A. shall not be liable for any loss or damage that may result from the unsuitable use of the equipment, negligence or user's full responsibility.

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Terragene S.A