QUICKSTART GUIDE

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


BIONOVA®
Hyper

Terragene®

Rev. 7 | Jun 2024

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EN Instructions for use

Note: Remember to set-up the Auto-reader local date and time before starting incubations.

Setting Up your Auto-reader

Place the Auto-reader on a firm surface, free from vibrations, away from direct sunlight, currents of hot or cold air, chemical and corrosive or flammable substances. Do not place the Auto-reader in a way that disconnection of the power supply's AC plug from mains could be difficult. Leave a space of at least 10 cm from the Auto-reader to the closest wall. Do not move the Auto-reader periodically or during its use. Connect the Auto-reader to a secure and stable electrical mains outlet socket. Do not wet or heat the Auto-reader. If liquid is spilled on the device, disconnect it and dry it immediately. Before turning on, verify that all reading positions are empty.

Using your Auto-reader

- 1| Power on your Auto-reader by connecting the power supply's AC plug to mains and then connect the plug at the other end of power supply to the rear of the Auto-reader. (See image 1)
- 2| Set-up the Auto-reader local date and time running the Bionova® Wireless Assistant on a compatible mobile device.
- 3| Wait until the Temperature progress & stability Indicator light stops blinking. Place a SCBI in a readout position. The incubation will start automatically. (See image 2)
- 4| Wait until the incubation is complete. The position status indicator light will turn green to indicate a negative result, or red to indicate a positive result. If a positive result is found, press the Configuration & alarm cancellation button to silence the alarm. Do not remove the SCBI until a result is informed. (See image 3)
- 5| Run Bionova® Cloud Web Application or the Bionova® Wireless Assistant mobile app to check the incubation result e-ticket.

Disposal: Discard the SCBIs according to your country's sanitary regulations.

Connectivity

To connect your Auto-reader to a PC by USB, connect one end of the included USB cable to the rear of the Auto-reader and the other end to an USB port on your PC. Connecting the Auto-reader to a PC allows using Bionova® Cloud Web Application. For more information, please refer to the user manual. To access the extended features of your Auto-reader using Wi-Fi or Bluetooth®, please download the Bionova® Wireless Assistant Application to your preferred mobile device and follow the steps detailed on the Application's tutorial.

Note: This Auto-reader is only compatible with 2.4 GHz WiFi networks.

Note: Due cybersecurity reasons, the Auto-reader will only be able to connect to WPA2-PSK WiFi networks.

Note: Bluetooth® connection mode is not available for iOS mobile devices.

For a detailed overview on the connection process, please refer to the user manual.

Connectivity credentials

During the connection process, the following credentials can come useful:

Device Internal Wi-Fi network credentials

Wi-Fi network name: HYPER_[batch][serial number]

Wi-Fi network password: abcd1234

Bluetooth® credentials

Device Bluetooth® ID: HYPER_[batch][serial number]

Product Compatibility and Incubation Programs

Auto-reader Incubation Programs

The following incubation programs are available for your Auto-reader.

Temperature selection	Time selection	Incubation program
60 °C	Hyper (5 min.)	Hyper at 60 °C

Note: Follow each Indicator's Instructions for use before its incubation.

Compatible indicators

We develop new indicators for extending the features of our Auto-readers regularly. Please check out your Auto-reader's compatibility matrix to find out more about all compatible indicators for your device.

Troubleshooting

Refer to the user manual for a detailed troubleshooting routine. Please visit www.terrargene.com/BIONOVA-HYPER or scan the QR code to access the Auto-reader user manual.

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