


QUICKSTART GUIDE

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BIONOVA®
IC10/20FR

Terragene®
Rev. 0 | June 2021

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EN Instructions for use

Setting Up your Auto-reader

Place the Bionova® IC10/20FR Auto-reader on a firm surface, free from vibrations, away from direct sunlight, currents of hot or cold air, chemical and corrosive or flammable substances. Do not place the Auto-reader in a way that disconnection of the power supply's AC plug from mains could be difficult. Leave a space of at least 10 cm from the Auto-reader to the closest wall. Do not move the Auto-reader periodically or during its use. Connect the Auto-reader to a secure and stable electrical mains outlet socket. Do not wet or heat the Auto-reader. If liquid is spilled on the device, disconnect it and dry it immediately. Before turning on, verify that all reading positions are empty.

Using your Auto-reader

- 1 | Turn on your IC10/20FR Auto-Reader by connecting the power supply's AC plug to mains and then connect the plug at the other end of power supply to the rear of the Auto-reader. (See image 1)
- 2 | Press and hold the temperature selection button for 3 seconds to enter the temperature selection mode. Press again to select the desired temperature. After 4 seconds the temperature modification will be accepted. (See image 2)
- 3 | Wait until the temperature stability indicator stops blinking.
- 4.A | SCBI incubation procedure:
4.A.1 | Select the appropriate incubation time for the SCBI you are planning to incubate by pressing the incubation time selection button. (See image 3.A)
- 4.A.2 | Place the SCBI in the chosen position. The incubation will start automatically. (See image 3.B)
- 4.B | Chemdy® PRO1 MICRO Indicator Incubation procedure
4.B.1 | Press the Pro button for 1 second to enter the reading mode. (See image 4.A)
- 4.B.2 | Wait until the qualitative indicator of protein concentration corresponding to 1 µg stops blinking. Place the protein pen in the corresponding position. Press the Pro button again for 2 seconds to start the incubation. (See image 4.B)
- 5 | Wait until the incubation is completed. The position status indicator light will turn green to indicate a negative result or red to indicate a positive result. For the PRO1 MICRO indicator, the concentration of the protein detected will also be indicated. If a positive result is detected, press the alarm button to cancel the alarm. Do not remove the indicator until the ticket with the result of the incubation is printed.
Note: To modify the protein concentration threshold to determine a positive result for the PRO1 MICRO Indicator, please refer the user's manual.

Disposal: Dispose of the indicators in accordance with your country's waste-disposal guidelines.

Connectivity

To connect your IC10/20FR Auto-reader to a PC by USB, connect one end of the included USB cable to the rear of the incubator and the other end to a USB port on your PC. Connecting the Auto-reader to a PC allows the use of the Bionova® Traceability Software. Please see the User's Manual for further information.

Product Compatibility and Incubation Programs

Bionova® IC10/20FR Auto-reader Incubation Programs
The following incubation programs are available for your Bionova® IC10/20FR Auto-reader.

Temperature	Time	Incubation Program
60 °C	20 min. (Ultra)	20 min. at 60 °C
	30 min.	30 min. at 60 °C
	1 hs	1 hour at 60 °C
	2 hs	2 hours at 60 °C
	3 hs	3 hours at 60 °C
37 °C	PRO (15 min.)	PRO at 60 °C
	4 hs	4 hours at 37 °C

⚠ Follow each Indicator's Instructions for use before its incubation.

Bionova® IC10/20FR Auto-reader compatible indicators

We develop new indicators for extending the features of our Auto-readers regularly. Please visit www.terrargene.com/IC1020FR to access the Compatibility Matrix for Terragene® Auto-readers. The Compatibility Matrix specifies all compatible Indicators developed for your Auto-reader up to date. Stay tuned for new and exciting updates!

Troubleshooting

Refer to the user manual for a detailed troubleshooting routine. Please visit www.terrargene.com/IC1020FR or scan the QR code to access the Auto-reader user manual.

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