
# QUICKSTART GUIDE

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## EN Instructions for use

**IMPORTANT:** Remember to set-up the Auto-reader local date and time before starting incubations.

### Setting Up your Auto-reader

Place the Bionova® Photon Auto-reader incubator on a firm surface, free from vibrations, away from direct sunlight, currents of hot or cold air, chemical and corrosive or flammable substances. Do not place the Auto-reader in a way that disconnection of the power supply's AC plug from mains could be difficult. Leave a space of at least 10 cm from the Auto-reader to the closest wall. Do not move the Auto-reader periodically or during its use. Connect the Auto-reader to a secure and stable electrical mains outlet socket. Do not wet or heat the Auto-reader. If liquid is spilled on the device, disconnect it and dry it immediately. Before turning on, verify that all reading positions are empty.

### Using your Auto-reader

1) Power on your Bionova® Photon Auto-reader incubator by connecting the power supply's AC plug to mains and then connect the plug at the other end of power supply to the rear of the Auto-reader. (See image 1)  
2) After turning on the device, a two-steps initialization sequence will start. Please wait until the sequence is complete before using the Auto-reader.

■ During the first step, the Auto-reader will set and stabilize the incubation temperature. When the correct temperature is reached, the Temperature progress & stability indicator will stop blinking.  
■ During the final step, the Auto-reader will set the readout system of each incubation position. When the readout system is ready for a reading, the Terragene logo will stop blinking. Once both indicators have stopped blinking, the initialization sequence is complete and the device is ready to process a sample.

**IMPORTANT:** The initialization sequence is carried out each time the incubator is turned on and can take between 10 and 30 minutes depending on the Auto-reader's current temperature. Plan accordingly.

3) Place a SCBI in a readout position. The incubation will start automatically. (See image 2)  
4) Wait until the incubation is complete. If a negative result is found, the position status indicator will turn green and a short audible sound will be emitted. If a positive result is found, the position status indicator will turn red and an audible alarm will be triggered, press the alarm/configuration button to silence the alarm. Do not remove the SCBI until a result is informed. (See image 3)  
5) Run Bionova® Traceability Software or the Bionova® Wireless Assistant mobile app to check the incubation result e-ticket.

**Disposal:** Discard the SCBIs following your Country disposal directives.

### Connectivity

To connect your Bionova® Photon Auto-reader incubator to a PC by USB, connect one end of the included USB cable to the rear of the incubator and the other end to an USB port on your PC. Auto-reader to PC Wi-Fi connectivity is also available and can be Set-up using the Bionova® Wireless Assistant. Connecting the Auto-reader to a PC (by USB, or Wi-Fi) allows using Bionova® Traceability Software to save, manage, and print e-tickets reports. Connecting the Auto-reader to an Android device (by Wi-Fi or Bluetooth®) allows using the Bionova® Wireless Assistant APP to check the status of the device, access e-tickets, and more. Check the User Manual for more information.

### Connectivity credentials

For connecting your Auto-reader to the Bionova® Wireless assistant APP please follow the steps detailed on the APP's Interactive tutorial. During the connection process, the following credentials can come handy:

#### Device Internal Wi-Fi network credentials

Wi-Fi network name: PHOTON\_[batch][serial number]  
Wi-Fi network password: abcd1234

#### Bluetooth® credentials

Device Bluetooth ID: PHOTON\_[batch][serial number]

Please check the device user manual for a detailed overview on the connection process.

### Product Compatibility and Incubation Programs

#### Bionova® Photon Auto-reader Incubation Programs

The following incubation programs are available for your Bionova® Photon Auto-reader.

| Temperature | Time             | Incubation Program |
|-------------|------------------|--------------------|
| 60 °C       | 7 sec. (Instant) | Instant at 60 °C   |

! Follow each Indicator's instructions for use before its incubation.

#### Bionova® Photon Auto-reader compatible indicators

We develop new indicators for extending the features of our Auto-readers regularly. Please visit [www.terragene.com/PHOTON](http://www.terragene.com/PHOTON) to access the Compatibility Matrix for Terragene® Auto-readers. The Compatibility Matrix specifies all compatible Indicators developed for your Auto-reader up to date. Stay tuned for new and exciting updates!

#### Troubleshooting

Refer to the user manual for a detailed troubleshooting routine. Please visit [www.terragene.com/PHOTON](http://www.terragene.com/PHOTON) or scan the QR code to access the Auto-reader user manual.

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