
# QUICKSTART GUIDE

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

Rev. 8 | Jun 2024

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## EN Instructions for use

**Note:** Remember to set-up the Auto-reader local date and time before starting incubations.

### Setting Up your Auto-reader

Place the Auto-reader incubator on a firm surface, free from vibrations, away from direct sunlight, currents of hot or cold air, chemical and corrosive or flammable substances. Do not place the Auto-reader in a way that disconnection of the power supply's AC plug from mains could be difficult. Leave a space of at least 10 cm from the Auto-reader to the closest wall. Do not move the Auto-reader periodically or during its use. Connect the Auto-reader to a secure and stable electrical mains outlet socket. Do not wet or heat the Auto-reader. If liquid is spilled on the device, disconnect it and dry it immediately. Before turning on, verify that all reading positions are empty.

### Using your Auto-reader

- 1) Power on your Auto-reader by connecting the power supply's AC plug to mains and then connect the plug at the other end of power supply to the rear of the Auto-reader. (See image 1)
- 2) Set-up the Auto-reader local date and time running the Bionova® Wireless Assistant on a compatible mobile device.
- 3) After turning on the device, a two-steps initialization sequence will start. Please wait until the sequence is complete before using the Auto-reader.
  - During the first step, the Auto-reader will set and stabilize the incubation temperature. When the correct temperature is reached, the Temperature progress & stability indicator light will stop blinking.
  - During the final step, the Auto-reader will set the readout system of each incubation position. When the readout system is ready for a reading, the Terragene® indicator light will stop blinking.
- 4) Once both indicator lights have stopped blinking, the initialization sequence is complete and the device is ready to process a sample.
- 5) The initialization sequence is carried out each time the Auto-reader is turned on and can take between 10 and 30 minutes depending on the Auto-reader's current temperature. Plan accordingly.
- 6) Place a SCBI in a readout position. The incubation will start automatically. (See image 2)
- 7) Wait until the incubation is complete. The position status indicator light will turn green to indicate a negative result, or red to indicate a positive result. If a positive result is found, press the Configuration & alarm cancellation button to silence the alarm. Do not remove the SCBI until a result is informed. (See image 3)
- 8) Run Bionova® Traceability Software or the Bionova® Wireless Assistant mobile app to check the incubation result e-ticket.

**Disposal:** Discard the SCBIs according to your country's sanitary regulations.

### Connectivity

To connect your Auto-reader to a PC by USB, connect one end of the included USB cable to the rear of the Auto-reader and the other end to an USB port on your PC. Connecting the Auto-reader to a PC allows using Bionova® Cloud Web Application. For more information, please refer to the user manual.

To access the extended features of your Auto-reader using Wi-Fi or Bluetooth®, please download the Bionova® Wireless Assistant Application to your preferred mobile device and follow the steps detailed on the Application's tutorial.

- Note:** This Auto-reader is only compatible with 2.4 GHz WiFi networks.
- Note:** Due cybersecurity reasons, the Auto-reader will only be able to connect to WPA2-PSK WiFi networks.

**Note:** Bluetooth® connection mode is not available for iOS mobile devices.

For a detailed overview on the connection process, please refer to the user manual.

**Connectivity credentials**  
During the connection process, the following credentials can come useful:

**Device Internal Wi-Fi network credentials**  
Wi-Fi network name: PHOTON\_[batch][serial number]  
Wi-Fi network password: abcd1234

**Bluetooth® credentials**  
Device Bluetooth® ID: PHOTON\_[batch][serial number]

**Product Compatibility and Incubation Programs**  
**Auto-reader Incubation Programs**  
The following incubation programs are available for your Auto-reader.

Temperature selection	Time selection	Incubation program
60 °C	Instant (7 sec.)	Instant at 60 °C

**Note:** Follow each Indicator's Instructions for use before its incubation.

**Compatible indicators**  
We develop new indicators for extending the features of our Auto-readers regularly. Please check out your Auto-reader's compatibility matrix to find out more about all compatible indicators for your device.

**Troubleshooting**  
Refer to the user manual for a detailed troubleshooting routine.  
Please visit [www.terragene.com/PHOTON](http://www.terragene.com/PHOTON) or scan the QR code to access the Auto-reader user manual.

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