
BT10 Biological Indicator

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Quality certification Certificado de calidad Bionova® BT10

EO sterilization / Esterilización por OE
Bacillus atrophaeus ATCC 9372

LOT	
	
	
Heat Shock Population / Población _____CFU / UFC	
D - value / Valor D _____min. (64 °C, 60% RH, EO 600 mg/L)	
Survival time / Tiempo de sobrevivencia _____min. Survival time = (log ₁₀ labeled population - 2) x labeled D-value	
Kill time / Tiempo de muerte _____min. Kill time = (log ₁₀ labeled population + 4) x labeled D-value	

Parameters determined at time of manufacture according to ISO 11138-1: 2017, ISO 11138-2:2017 and IRAM 37102: 1999 (Parts 1 and 2) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant.	
ATCC is a registered trademark of American Type Culture Collection.	
	
Lic. Adrian J. Rovetto Director Técnico Technical Director	

Uso exclusivo para profesionales e Instituciones Sanitarias.			
Explanation of Symbols		Manufacture Date.	
	Product designed for use with EO sterilization cycles.		Expiration Date.
	Batch number		Manufacturer.

EN Biological Indicators For Ethylene Oxide sterilization

Indications for Use
Terragene® Bionova® BT10 Self-contained Biological Indicators (SCBIs) have been designed for quick and easy monitoring of Ethylene Oxide sterilization cycles.

Product description
Bionova® BT10 is a SCBI designed for monitoring Ethylene Oxide (EO) sterilization processes. Each tube contains a population of spores of *Bacillus atrophaeus* ATCC 9372 soaked on a carrier. It also contains blue culture medium within a glass ampoule, which is on top of the carrier, on the base of the tube. If the sterilization process is not successful, culture medium will turn from blue to yellow after incubation at 37 ± 2 °C, thus indicating the presence of living *Bacillus atrophaeus* spores.

If the sterilization process is successful, the indicator medium will remain blue after incubation. The final readout should be performed after 48 hours of incubation at 37 ± 2 °C.

Precautions
Do not use Terragene® Bionova® BT10 SCBI to control Steam, Dry Heat, Formaldehyde or other sterilization processes different from specified.

Do not re use biological indicators. Do not reuse the sterilizer until the biological indicator test results are negative.

Instructions for use
1. Identify the Bionova® BT10 Indicator tube by writing the sterilizer number (in case of having more than one sterilizer), load number, and processing date on the indicator label.
2. Pack the biological indicator along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place this package in those areas which a priori you consider most inaccessible for the sterilizing agent (EO). Generally a problematic area is the center of the load or near the door.
3. Sterilize as usual.
4. After the sterilization process has finished, you can:

A. Remove the biological indicator from the sterilization package for incubation prior to material aeration.

B. Aerate the package containing the biological indicator with the rest of the load and then remove the biological indicator for incubation.
CAUTION: Wear safety glasses and gloves when removing the Bionova® BT10 Biological Indicator from the sterilized package.
WARNING: Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.

5. Check the chemical indicator on the label of the biological indicator. A color change to green confirms that the biological indicator has been exposed to Ethylene Oxide.
IMPORTANT: This color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.

6. Crush the ampoule contained in the biological indicator with an individual ampoule crusher or with the ampoule crusher placed within the top of the incubator’s incubation area.

Then shake the tube down vigorously, with movements similar to those performed to lower the temperature of a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the biological indicator in the incubator.
IMPORTANT: Use a non-sterilized SCBI as a positive control at least once per day, when a sterilization cycle is run. The positive control ensures that correct incubation conditions were met; the capability of the medium to promote rapid growth; viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals and proper functioning of Terragene® Bionova® Auto-Reader Incubators. Both, the positive control indicator and the processed indicator, should belong to the same batch.

7. Incubate the processed biological indicator and the indicator used as positive control for a maximum of 48 hours at 37 ± 2 °C. Readout should be performed at convenient intervals of 10 hours. A color change from blue to yellow of the growth indicator medium means a sterilization process failure has occurred. If after 48 hours there is no color change in the processed indicators, a final negative result is obtained (the sterilization process was satisfactory).

The positive control indicator should show a yellow color change for the results to be valid.

NOTE: To extend the incubation time beyond 48 hours, we recommend using a humidified environment to prevent complete evaporation of the culture medium contained within the indicator. Record the positive results and discard the biological indicators immediately as described below.

Monitoring frequency

Follow the policies and procedures with the monitoring frequency specified by the professional associations and/or standards corresponding to your country. As recommended practice, and to provide optimal patient safety, Terragene® recommends that each sterilization load be monitored with the appropriate biological indicator.

Storage
Store in a dark place at temperature between 10-30 °C and 30-80 % relative humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Shelf life
Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date. Chemical Process Indicator on SCBI label have an expiration date of 2 years when used as part of SCBI.
Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal
Discard biological indicators after use according to your country’s healthcare and safety regulations. The positive biological indicator can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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■ Terragene S.A