
BT100 Biological Indicator

R10 / 02.2025

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Quality certification Certificado de calidad Bionova® BT100

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Survival time = (log₁₀ labeled population - 2) x labeled D-value

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Kill time = (log₁₀ labeled population + 4) x labeled D-value


Parameters determined at time of manufacture according to ISO 11138-1:2017, ISO 11138-5:2017 and IRAM 37102-1:1999 standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant.
ATCC® is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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For formaldehyde sterilization

Composition

Each tube contains a population of spores of *Geobacillus stearothermophilus* ATCC® 7953 soaked on a carrier. It also has a purple colored culture indicator medium contained within the glass ampoule, next to the carrier, on the base of the tube.

Product description

Bionova® BT100 Biological Indicators have been specifically designed for monitoring formaldehyde sterilization processes. If the sterilization process has not been successful, the indicator medium will turn from purple to yellow after incubation between 55-62 °C, thus indicating the presence of living *Geobacillus stearothermophilus* spores. If the sterilization process was effective, the indicator medium will remain purple after incubation. The final readout should be performed after 48 hours of incubation between 55-62 °C.

Precautions

Do not use Bionova® BT100 Biological Indicators for monitoring steam, ethylene oxide, dry heat, plasma or other sterilization processes.

Do not reuse the biological indicators.

Do not reuse the sterilizer until the biological indicator test results are negative.

Storage

Best stored in a dark place under the following conditions: Temperature between 10-30 °C, 30-80 % relative humidity.

Do not freeze.

Do not store biological indicators near sterilizing agents or other chemical products.

Instructions for use

1. Identify the Bionova® BT100 Indicator by writing the sterilizer number (in case of having more than one), load number, and processing date on the indicator label.

2. Pack the biological indicator along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place this package in those areas which a priori you consider most inaccessible for sterilizing agent. Generally, a problematic area is the center of the load and areas near the door.

3. Sterilize as usual.

4. After the sterilization process has finished, you can:

A. Remove the biological indicator from sterilization package for incubation prior to material aeration.

B. Aerate the package containing the biological indicator with the rest of the load and then remove the biological indicator for incubation. **CAUTION:** Wear safety glasses and gloves when removing the Bionova® BT100 Biological Indicator from the sterilized package.

WARNING: Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.

5. Check the chemical indicator on the label of the biological indicator. A color change to green confirms that the biological indicator has been exposed to formaldehyde.

IMPORTANT: This color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.

6. Crush the ampoule contained in the biological indicator with an individual ampoule crusher or with the ampoule crusher placed within the top of the incubator's incubation area.

Then shake the tube down vigorously, with movements similar to those performed to lower the temperature of a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the biological indicator in the incubator. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control in order to ensure that correct incubation conditions were met; capability of medium to promote rapid growth; viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals. Both the positive control indicator and the processed indicator should belong to the same batch.

7. Incubate the processed biological indicator and the indicator used as positive control for a maximum of 48 hours between 55-62 °C. Readouts should be performed at convenient intervals of 10 hours.

A color change from purple to yellow of the growth indicator medium means that a sterilization process failure has occurred. If after 48 hours there is no color change in the processed indicators, a final negative result is obtained (the sterilization process was effective). The positive control should show a color change from purple to yellow for the results to be valid. **NOTE:** To extend the incubation time to more than 48 hours, we recommend using a humidified environment to prevent the complete evaporation of the culture medium contained in the indicator.

Record the positive results and discard them immediately as shown below.

Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicators can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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