
BT102 Biological Indicator

Rev. 8 / 09.2021

Rapid Readout Fluorescence System

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ISO and USP Compliant
ATCC is a registered trademark of American Type Culture Collection
Lic. Adrián J. Rovetto Director Técnico Technical Director

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Explanation of Symbols	
FORM Product designed for use in Formaldehyde sterilization processes	Expiration Date
LOT Batch number	Manufacturer

EN Biological Indicators

For Formaldehyde sterilization

Composition
Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores soaked on a spore carrier. It also has purple colored culture medium contained within the glass ampoule, over the filter next to the carrier, on the base of the tube.

Product description
Bionova® BT102 Fluorescence Rapid Readout Biological Indicators have been designed for fast and easy monitoring of Formaldehyde sterilization cycles. The system consists of a plastic tube & a special filter, a spore carrier and a glass ampoule with culture medium. The tube has a plastic cap with holes and a barrier permeable to Formaldehyde.

Precaution

Do not use Bionova® BT102 Biological Indicators to control EO, Dry Heat, Hydrogen Peroxide or other sterilization processes different from the specified in this IFU. Do not reuse biological indicators.

Do not reuse the sterilizer until the biological indicator test result is negative.

Instructions for use
1. Identify the Bionova® BT102 Biological Indicator by writing the sterilizer number (in case of having more than one), load number and processing date on the label.
2. Pack the biological indicator along with materials to be sterilized in appropriate packages according to recommended sterilization practices. Place the package in those areas considered a priori most inaccessible for the sterilizing agent (Formaldehyde). Generally, the center of the load and areas near the door are more problematic.
3. Sterilize as usual.
4. After the sterilization process has finished, open the sterilizer door, wait five minutes and remove the biological indicator from the package. **CAUTION:** Wear safety glasses and gloves when removing the Bionova® BT102 Biological Indicator from the sterilized package. **WARNING:** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.
5. Let the biological indicator cool down until it reaches room temperature.
6. Check the chemical indicator on the label of the biological indicator. A color

change to green indicates that the biological indicator has been exposed to Formaldehyde. **IMPORTANT:** this color change does not evidence the process effectiveness to achieve sterility. If the chemical indicator color has not changed, check the sterilization process.

7. Press the cap to seal the tube. Crush the ampoule contained in the biological indicator with an individual ampoule crusher or with the ampoule crusher placed within the top of the incubator & incubate for 48 hours.
8. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature of a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the biological indicator in the incubator.

IMPORTANT: Use a non-sterilized SCBI as a positive control at least once per day, when a sterilization cycle is run. The positive control ensures that correct incubation conditions were met, capability of medium to promote rapid growth; viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals and proper functioning of Terragene® Bionova® Auto-Reader Incubators. Both, the positive control indicator and the processed indicator, should belong to the same batch.
8. Incubate the processed biological indicator and the indicator used as positive control for a maximum of 2 hours at 60±2 °C for Rapid Readout. Fluorescence readout by the reader (380 nm / r/emission 455-465 nm) means a failure in the sterilization process. If no fluorescence is detected at the end of the 2-hour incubation, then, the result is negative (the sterilization process has been effective). For results to be valid, the indicator used as positive control must be detected as such by the reader. It is a good practice to incubate a positive control for a visual color change.

Record the positive results and discard the biological indicators immediately, as shown below.

Rapid Readout: 2 hours

The Rapid Readout must be carried out in the appropriate Bionova® Auto-reader Incubators. Fluorescence is produced when the reader stimulates the carrier with UV light. Final readout of negative results is readily available after 2 hours of incubation. Fluorescence readout is an indirect measure of the ability of *Geobacillus stearothermophilus* spores to undergo sterilization process (spore survival). If there is a failure in the sterilization process can also become evident by the color change of the medium. Due to the high sensitivity of the fluorescence results at 2 hours, conventional incubation for color change of Bionova® BT102 Rapid Readout Biological Indicator is not an advantage.

Visual confirmation: 48 hours

Optional, you can confirm a visual color-change confirmation after a 48-hour incubation. If the sterilization process has not been successful, culture medium will turn to yellow during incubation at 60 °C, indicating the presence of living spores. If sterilization was successful, culture medium will remain purple after the incubation process.

A definitive negative result is obtained after the 48-hour incubation.

The positive control must show a color change from purple to yellow for results to be valid.

Visual confirmation: 7 days

A 7-day readout is optional and not intended to be routinely performed. This is an initial validation of the 2-hour readout. The 2-hour fluorescence results are compared to the 7-day visual readings. The sensitivity of the system is determined as the difference between the 7 days positive indicators and false negative indicators (negative visual results versus visually positive) in relation to the positives at 7 days (*). Based on the claimed ≥ 97% readout reliability of the 2-hour biological indicator, there is no advantage to incubate the indicator beyond 2 hours. **NOTE:** If 7-day readout is made, a humidified environment will be required to prevent medium from drying out before the seventh day.

Monitoring frequency

Follow facility policies and procedures which should specify a biological indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As the best practice and to provide optimal patient safety, Terragene® recommends that every sterilization load be monitored with an appropriate biological indicator.

Storage

Store in a dark place at temperature between 10-30 °C and 30-80 % relative humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Shelf life

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date. Chemical Process Indicator on SCBI label have an expiration date of 2 years when used as part of SCBI. Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal
Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicators can be autoclaved in a gravity air displacement sterilizer at 121 °C for 30 minutes or 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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