Spores suspension

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Spores Suspension

BT10S/2, BT10S/3, BT10S/4, BT10S/5,
BT10S/6, BT10S/7, BT10S/8, BT10S/9*Bacillus atrophaeus* ATCC 9372

Composition

Each Bionova® BT10S Spore suspension contains 10 ml of bacterial spores of *Bacillus atrophaeus* ATCC 9372 in 40% of ethanol.

Product description

Bionova® BT10S Spore Suspension is designed for direct inoculation of products or for the preparation of customized Biological Indicators used for control of Ethylene Oxide sterilization processes and/or Dry Heat sterilization processes.

WARNING: Use BT10S biological indicator only for the control of EO and/or Dry Heat sterilization processes. Do not use for other processes such as Steam, Formaldehyde and Radiation sterilization.

Storage

Store the product in a dark place, at temperature between 2-8 °C, 35-60 % Relative Humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

Directions for use

1. Perform inoculation operations in a clean area.
2. Samples to be inoculated should be representative of product being sterilized.
3. Select the spore suspension of appropriate population for the sterilization method use. For most purposes, inoculation of products with approximately one million (10⁶) spores provides a suitable challenge.

NOTE: Spore suspensions are standardized on the basis of number of spores per 0.1 ml of suspension.

4. Shake vial vigorously before each use.

5. Use a sterile pipette or syringe to accurately deliver the volume of suspension to be utilized.

6. Deposit suspension onto product. The area to be inoculated should be the most difficult for the sterilizing agent to penetrate.

7. Allow product to dry at 20-30 °C for approximately 24 hours.

8. Package inoculated product exactly like product being sterilized and identify as "Inoculated Test Sample".

9. Distribute "Inoculated Test Sample" throughout the sterilizer load.

10. After sterilization cycle is complete, test the inoculated products as soon as possible using an appropriate culture medium, and incubate at 37 °C. We suggest the use of Bionova® MC1020 or MC1030 culture media for a final result at 48 hours of incubation. For a conventional culture medium such as TSB, incubate at the optimal temperature for up to 7 days.

IMPORTANT: Each time that a processed product is incubated use an unprocessed product as a positive control. The positive control ensures that the incubation conditions were met. An uninoculated product can be used as a negative control.

A color change to yellow of the incubation medium shows a failure in the sterilization process. If after 48 hours of incubation, the purple color of the medium remains the same, the result is negative (the sterilization process was efficient). The positive control must change to yellow in order to validate the results.

The positive vials have to be registered and then have to be processed as indicated below.

WARNING: Do not use the sterilizer until the biological indicator shows a negative result. (The processed indicator remains purple after 48 hours of incubation).

Disposal

Discard Bionova® BT10S spore suspension and growth medium tubes according to your country's healthcare and safety regulations. Positive growth medium tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.

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Survival time = (log₁₀ labeled population - 2) x labeled D-value

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Kill time = (log₁₀ labeled population + 4) x labeled D-value

Valor Z / _____ °C
Z-Value

Performance values calculated using a filter paper inoculated with 0.1 ml of spore suspension as a BI.

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Parameters determined at time of manufacture according to ISO 11138-1:2006, ISO 11138-2:2006, ISO 11138-4:2006 and IRAM 37102: 1999 (Parts 1 and 2). The shown values are reproducible only under the same conditions under which they were determined.

ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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