BT202 Biological Indicator

Rev. 7 / 09.2019

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For Steam sterilization

Composition

Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores inoculated on a carrier. The tube also has a purple growth indicator medium contained within the glass ampoule.

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Bionova® BT202 Biological Indicator is specifically designed for monitoring vacuum-assisted and gravity air displacement Steam sterilization processes at 121-135 °C. If the sterilization process was not successful, the color of the culture medium will change from purple to yellow after incubation at 60 ± 2 °C, thus indicating the presence of live spores of *Geobacillus stearothermophilus*. If the sterilization process was successful, the color of culture medium will remain purple after incubation at 60±2 °C. Terragene® recommends the usage of Bionova® incubators and final readout time of 10 hours. Bionova® BT202 Biological Indicators are fully compatible with incubators currently commercialized by other brands.

Precautions

Do not use Bionova® BT202 Biological Indicator for monitoring Ethylene Oxide, Dry Heat, Chemical Vapor, Radiation or other sterilization processes.

Do not reuse the biological indicators.

Do not use the sterilizer until the biological indicator test results negative.

Shelf life: 2 years

Store in a dark place at temperatures between 10-30 °C, 30-80 % relative humidity.

Do not freeze.

Do not store biological indicators near sterilizing agents or other chemical products.

Instructions for use

1. Identify the Bionova® BT202 Indicator by writing the sterilizer number (in case of having more than one), load number and processing date on the indicator label.

2. Pack the biological indicator along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place this package in those areas which a priori are considered most inaccessible for the sterilizing agent (Steam). Generally, the center of the load and areas near the door are problematic.

3. Sterilize as usual.

4. After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the biological indicator from package. **CAUTION!** Wear safety glasses and gloves when removing the Bionova® BT202 Biological Indicator from the sterilized package. **WARNING!** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.

5. Let the biological indicator cool down until it reaches room temperature.

6. Check the chemical indicator on the label of the biological indicator. A color change to brown confirms that the biological indicator has been exposed to steam. **IMPORTANT:** this color change does not evidence the process effectiveness to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.

7. Crush the ampoule contained in the biological indicator with an individual ampoule crusher or with the ampoule crusher placed within the top of the incubator's incubation area. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature of a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the biological indicator in the incubator. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control in order to ensure that appropriate incubation conditions were met; capability of media to promote rapid growth; viability of spores have not been altered due to improper storage temperature, humidity or proximity to chemicals and proper functioning of Bionova® Incubator. Both the positive control indicator and the processed indicator should belong to the same batch.

8. Incubate the processed biological indicator and the indicator used as a positive control for a maximum of 10 hours at 60±2 °C, check results

after 3-5 hs incubation. Color change of the growth indicator medium from purple to yellow means that a sterilization process failure has occurred. If at 10 hours of incubation processed indicators show no color change, a final negative result is obtained (the sterilization process was efficient). The culture medium of the positive control indicator must show a color change from purple to yellow for the result to be valid.

NOTE: To extend the incubation time to more than 10 hours, we recommend using a humidified environment to prevent the complete evaporation of the culture medium contained in the indicator. Record the positive results and discard the biological indicators immediately as described below.

Disposal

Dispose biological indicators after use according to your country's healthcare and safety regulations. Positive biological indicators can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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ATCC is a registered trademark of American Type Culture Collection.
Terragene S.A.