Spores Suspension

BT20S/2 | BT20S/3 | BT20S/4 | BT20S/5 | BT20S/6 |
BT20S/7 | BT20S/8 | BT20S/9

Quality certification Certificado de calidad Bionova® BT20S/x

FORM STEAM
VH2O2

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Performance values for Steam and Formaldehyde were calculated using a filter paper inoculated with 0.1 ml of spore suspension as a BI. For Hydrogen Peroxide, a polypropylene tube inoculated with 0.1 ml of spore suspension was used as a BI.

Parameters determined at time of manufacture according to ISO 11138-1:2017, ISO 11138-3:2017, ISO 11138-5:2017 and IRAM 37102:1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Technical Director
Director Técnico

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Geobacillus stearothermophilus ATCC 7953

Composition

Each Bionova® BT20S Spore Suspension BT20S contains 10ml of bacterial spores of *Geobacillus stearothermophilus* ATCC 7953 in 40% of ethanol.

Product description

Bionova® BT20S Spore Suspension has been designed for direct inoculation of products or for the preparation of customized Biological Indicators used for control of Steam, Formaldehyde, Hydrogen Peroxide and/or Peracetic Acid sterilization processes. The combined use of inoculated products, pre-sterilization microorganism count determinations (bioburden), and physical process measurements and controls should be employed to validate sterilization processes.

WARNING: Use BT20S biological indicator only for the control of Steam, Formaldehyde, Hydrogen Peroxide and/or Peracetic Acid sterilization processes. Do not use for other processes such as Ethylene Oxide, Dry Heat and/or Radiation sterilization.

Storage

Store in a dark place and at temperatures between 2-8 °C, 30-80 % Relative Humidity. Do not freeze.

Do not store near sterilizing agents or other chemical products.

Directions for use

1. Perform inoculation operations in a clean area.
2. Samples to be inoculated should be representative of product being sterilized.

3. Select the spore suspension of appropriate population for the sterilization method use. For most purposes, inoculation of product with approximately one million (10⁶) spores provides a suitable challenge.

NOTE: Spore suspensions are standardized on the basis of number of spores per 0.1 ml of suspension.

4. Vortex vial before each use.

5. Use a sterile pipette or syringe to accurately deliver the volume of suspension to be utilized.

a) If syringe is used, disinfect stopper surface and pull syringe plunger halfway back. Insert needle through stopper, push plunger in, and then slowly withdraw plunger to fill syringe to desired volume.

b) If a pipette is used, remove stopper and insert pipette. Dispense the desired volume.

6. Inoculate the product appropriately. The area to be inoculated should be the most difficult for the sterilizing agent to penetrate.

7. Identify the inoculated samples.

8. Distribute "Inoculated Test Samples" throughout the sterilizer load.

9. After sterilization cycle is completed, test the inoculated products as soon as possible using an appropriate culture medium, and incubate at 60 °C. We suggest Bionova® MC20-2 or MC1020-2 media to final reading at 24 hours for Steam/Hydrogen Peroxide sterilization and at 48 hours for Formaldehyde sterilization. For conventional culture medium like TSB, incubation for 7 days at 60 °C is recommended.

IMPORTANT: Each time a processed product is incubated, use an unprocessed product sample as a positive control. The positive control ensures that the incubation conditions were met. An uninoculated product can be used as a negative control. The positive vials have to be registered and then have to be processed as indicated below.

WARNING: Do not use the sterilizer until the biological indicator shows a negative result.

Disposal

Discard Bionova® BT20S spore suspension and growth medium tubes according to your country's healthcare and safety regulations. Positive growth medium tubes can be sterilized in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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