BT221 Biological Indicator

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Quality certification Certificado de calidad Bionova® BT221

STEAM

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(* Sensitivity/ _ (N° Positives at 7 days) - (N° False negatives) Sensibilidad (N° Positives at 7 days) X 100 ≥ 97%

Parameters determined at time of manufacture according to ISO 11138-1: 2017, ISO 11138-3: 2017 and IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant. ATCC is a registered trademark of American Type Culture Collection

Terragen® and Bionova® are registered trademarks of Terragene S.A.

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Explanation of Symbols	
 Product is designed for use with	 Manufacture Date
 Steam sterilization cycles	 Expiration Date
 Batch number	 Manufacturer

Intended Use Table		
	Autoclave/Steam Cycles	
	Gravity Displacement	
	132 °C	135 °C
	3 minutes 10 minutes	3 minutes 10 minutes
BT221	Fluorescence Read Time	pH Color Change
	1 hours	48 hours

EN Biological Indicators

For gravity air displacement Steam sterilization

Indications for use

United States

Terragen® Bionova® SCBI (BT221) is a self-contained biological indicator inoculated with viable *Geobacillus stearothermophilus* bacterial spores and is intended for monitoring the efficacy of Steam sterilization processes. See Intended Use Table for reference.

Outside the United States

Terragen® Bionova® BT221 Fluorescence Rapid Readout Biological Indicators have been designed for quick and easy monitoring of Gravity Air Displacement Steam sterilization cycles at 132-135 °C.

Device description

Terragen® Bionova® BT221 Fluorescence Rapid Readout Biological Indicators are single-use Self-Contained Biological Indicators (SCBI) that consist of a polypropylene tube, a spore carrier and a glass ampoule with a culture medium, enclosed with a colored cap. Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores inoculated on a spore carrier, a plastic cap with holes and a barrier permeable to Steam. Each BT221 has a Process Indicator on label that changes from pink to brown when exposed to Steam.

Precautions

WARNING: Do not use Terragene® Bionova® BT221 SCBI to monitor 121 °C or 132 °C vacuum assisted steam sterilization cycles, 121 °C gravity steam sterilization cycles or other cycles not mentioned in the instructions for use. Do not use Terragene® Bionova® BT221 SCBI to control Hydrogen Peroxide, EO, Dry Heat, Formaldehyde or other sterilization processes different from heat/steam. Do not reuse SCBI.

WARNING: Place one or more SCBIs in sterilizing hard-to-reach areas to ensure all areas of the chamber are sterilized. Please evaluate all load configurations to ensure ALL hard-to-reach areas have been identified, and place a SCBI in each of those locations.

WARNING: Do not reuse the sterilizer until the SCBI test result is negative.

Instructions for use

1. Identify the Terragene® Bionova® BT221 SCBI by writing the sterilizer number (in case of having more than one), load number and processing date on the label. 2. Pack the SCBI along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place the package in those areas which are considered most inaccessible for the sterilizing agent (e.g., the center of the load and areas near the door). 3. Sterilize as usual. 4. After the sterilization process has finished, open the sterilizer door, wait five minutes and remove the SCBI from the package. **CAUTION:** Wear safety glasses and gloves when removing the Terragene® Bionova® BT221 SCBI from the sterilized package. **WARNING:** Do not crush or handle the SCBI excessively, since this might cause the glass ampoule to burst. Let the SCBI cool down until it reaches room temperature. 5. Check the Process Indicator on SCBI label. A color change to brown indicates that the SCBI has been exposed to Steam. **IMPORTANT:** This color change does not evidence the process effectiveness to achieve sterility. If the Process Indicator color has not changed, check the sterilization process. 6. Press the lid to seal the tube. Crush the ampoule contained in the SCBI with an individual ampoule crusher or with the ampoule crusher placed within the incubator's incubation area. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the SCBI in the incubator. **IMPORTANT:** Use a non-sterilized SCBI as a positive control in order to ensure that correct incubation conditions were met; capability of medium to promote rapid growth; viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals and proper functioning of Terragene® Bionova® Auto-Reader Incubators. Both, the positive control indicator and the processed indicator, should belong to the same batch. 7. Incubate the processed indicator and the positive control indicator in the appropriate Terragene® auto-reader incubators for a maximum of 1 hour at (60 ± 2) °C for rapid readout. **NOTE:** Holding time between sterilization and incubation should not exceed a 7-day period. Fluorescence detection by the reader (excitation 340-380 nm / emission 455-465 nm) means a failure in the sterilization process. If no fluorescence is detected at 1-hour incubation, the result is negative. The positive control must give positive fluorescence readout. It is good practice to incubate a positive control for a visual color change. Record the positive results and discard the SCBI immediately, as it is indicated below.

Rapid readout: 1 hour

The Rapid Readout must be carried out in the appropriate Terragene® Bionova® auto-reader incubators. Fluorescence is emitted after the reader stimulates the spore carrier with UV light. Final readout of negative results is readily available after 1 hour incubation. Fluorescence readout is an indirect measure of the germination and growth of *Geobacillus stearothermophilus* spores, which have survived the sterilization process (positive result). Furthermore, a failure in the sterilization process can also become evident by culture medium color change. Due to the high sensitivity of the fluorescence results at 1 hour, conventional incubation for color change of Terragene® Bionova® SCBI BT221 is not an advantage.

Visual confirmation: 48 hours

Optionally, you can perform a visual color-change confirmation after a 48-hour incubation. If the sterilization process has not been successful, culture medium will turn yellow during incubation at (60 ± 2) °C, indicating the presence of living spores. If sterilization was successful, culture medium will remain purple after incubation. The positive control must show a color change from purple to yellow for results to be valid.

Readout time: 7 days

A 7-day readout is optional and may be routinely performed. This is an initial validation of the 1-hour readout. Fluorescence results are compared to the 7-day visual readouts. **NOTE:** If 7-day readout is performed, a humidified environment will be required to prevent medium from drying out.

Monitoring frequency

Follow facility policies and procedures which should specify a biological indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As the best practice and to provide optimal patient safety, Terragene® recommends that every sterilization load be monitored with an appropriate biological indicator.

Storage

Store in a dark place at temperature between 10-30 °C and 30-80 % relative humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Model	Autoclave/Steam Cycles	
	Gravity Displacement	
	132 °C	135 °C
	3 minutes 10 minutes	3 minutes 10 minutes
BT221	Fluorescence Read Time	pH Color Change
	1 hours	48 hours

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