BT222-T

Biological Indicator

Super Rapid Readout Fluorescence System

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Parameters determined at time of manufacture according to ISO 11138-1: 2006, ISO 11138-3: 2006 and IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions at which they were determined.

ISO and USP Compliant
ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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For Steam sterilization

Composition

Bionova® BT222-T Biological Indicator consists of a plastic tube with a plastic cap with holes and a steam permeable barrier. Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores soaked in a paper strip (spores carrier) and a purple color culture medium within a glass ampoule, which is over a filter and the carrier on the base of the tube.

Product description

Bionova® BT222-T Super Rapid Readout Biological Indicators have been designed quick and easy monitoring of steam sterilization cycles at 121-135 °C.

Super rapid readout: 1 hour

The Super Rapid Readout must be carried out in Bionova® IC10/20FR Reader Incubator, which is capable of performing an optical reading of the fluorescence emission of the product resulting from breaking a specific substrate on the spores carrier. Fluorescence is produced when the reader stimulates the carrier with UV light at 360 nm. Final reading of negative results is readily available at 1-hour incubation. It is a direct measure of the germination and growth of the spores of *Geobacillus stearothermophilus* spores which have survived the sterilization process (positive result). Furthermore, a failure in the sterilization process can also become evident by the color change of the medium. Due to the high sensitivity of the fluorescence results after 1 hour, conventional incubation for color change of BT222-T Super Rapid Readout Indicator is not an advantage.

Warning!

Do not use Bionova® BT222-T Biological Indicators to monitor sterilization cycles not mentioned in instructions for use. Do not use Bionova® BT222-T Biological Indicators to control EO, Dry Heat, Formaldehyde or other sterilization processes. Do not reuse the biological indicators.

Storage

Store in a dark place, at 10-30 °C temperature, 30-80 % Relative Humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Instructions for use

1. Identify the Bionova® BT222-T Biological Indicator by writing the sterilizer number (in case of having more than one sterilizer), load number and processing date on the label.
2. Pack the biological indicator along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place the package in those areas which a priori you consider most inaccessible for the sterilizing agent (Steam). Generally, a problematic area is the center of the load and near the door.
3. Sterilize as usual.

4. After the sterilization process has finished, open the sterilizer door, wait five minutes and remove the biological indicator from the package. **CAUTION!** Wear safety glasses and gloves when removing the Bionova® BT222-T biological indicator from the sterilized package. **WARNING!** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.
5. Let the biological indicator cool until it reaches room temperature.

6. Check the chemical indicator on the label of the biological indicator. A color change to brown indicates that the biological indicator has been exposed to steam. **IMPORTANT:** This color change does not evidence the process effectiveness to achieve sterility. If the chemical indicator color has not changed, it is necessary to check the sterilization process.

7. Press the cap to seal the tube. Then, break the ampoule contained in the biological indicator. This may be done in 3 different ways:

A- With an ampoule crusher.
B- Manually, wearing protective gloves to avoid injuries. To do this, take the tube by placing your index finger and thumb on the space between the cap and the line of chemical indicator and then compress the tube.

C- With the ampoule crusher placed within the top of the IC10/20FR's incubation area. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spores carrier entirely. Finally, place the biological indicator in the incubator.

IMPORTANT: Use a non-sterilized biological indicator as a positive control in order to ensure that correct incubation conditions were met: capability of medium to promote rapid growth; viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals and proper functioning of IC10/20FR Bionova® Incubator. Both the positive control indicator and the processed indicator should belong to the same batch.

8. Incubate the processed biological indicator and the indicator used as positive control for a maximum of 1 hour at 60±2 °C for the Super Rapid readout.

Fluorescence detection by the reader (stimulation at 340-380 nm/emission at 455-465 nm) means a failure in the sterilization process. If no fluorescence is detected at the end of the 1-hour incubation, then the result is negative (sterilization process has been effective). The indicator used as positive control must be detected as such by the reader.

Record the results and discard biological indicators immediately as explained below.

WARNING! Do not reuse the sterilizer until the biological indicator test result is negative.

Visual confirmation: 48 hours

Optionally, you can perform a visual color-change confirmation after a 48-hour incubation. If the sterilization process hasn't been successful, culture medium will change to yellow during incubation at 60 °C, indicating the presence of living spores. If sterilization was successful, culture medium will remain purple after the incubation process.

A definitive negative result is obtained after the 48-hour incubation. The positive control must show a color change from purple to yellow for results to be valid.

Readout time

A 7-day readout is optional and not intended to be routinely performed. This is an initial validation of the 1-hour reading.

The 1-hour super rapid readout and the optional 48-hour visual pH color change incubation times have been correlated with a 7-day incubation period following the FDA's Reduced Incubation Time protocol. Sterilized indicators were examined at 48 hours and 7 days for detection of a visual pH color change.

The 1-hour fluorescence change readings and the 48-hour visual pH color change readings were compared to the 7-day visual pH color change readings to determine the readout time of the indicator.

Bionova® BT222-T Biological Indicators have 1-hour reduced incubation time results that correlate to the 7-day (168 hours) visual readout results in ± 97% of cases.

Bionova® BT222-T Biological Indicators have 48-hour incubation results that correlate to the 7-day (168 hours) visual readout results in ± 97% of cases.

Due to the high reliability of the 1-hour fluorescence result there is no advantage in incubating Bionova® BT222-T biological indicators beyond 1 hour.

NOTE: If 7-day readout is made, the incubator will be required to prevent media from drying out.

Monitoring frequency

Follow facility Policies and Procedures which should specify a biological indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a best practice and to provide optimal patient safety, Terragene recommends that every steam sterilization load be monitored with an appropriate biological indicator.

Disposal

Bionova® BT222-T are normally discarded after the fluorescent readout result has been recorded. Dispose biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

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Fabricado por Terragene S.A. - C. Güemes 2879 - (2000)