
BT224 Biological Indicator

Rev. 9 / 08.2023

BIONOVA®

Ultra Rapid Readout Fluorescence System

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Quality certification Certificado de calidad Bionova® BT224

Steam sterilization / Esterilización por Vapor

Geobacillus stearothermophilus ATCC® 7953

LOT	
Heat shock population / Población _____ CFU / UFC	
D - value / Valor D (121 °C) _____ min.	
Survival time / Tiempo de sobrevida _____ min.	
Survival time = (log ₁₀ labeled population - 2) x labeled D-value	
Kill time / Tiempo de muerte _____ min.	
Kill time = (log ₁₀ labeled population + 4) x labeled D-value	
D - value / Valor D (132 °C) _____ sec./seg.	
Survival time / Tiempo de sobrevida _____ min.	
Survival time = (log ₁₀ labeled population - 2) x labeled D-value	
Kill time / Tiempo de muerte _____ min.	
Kill time = (log ₁₀ labeled population + 4) x labeled D-value	
D - value / Valor D (135 °C) _____ sec./seg.	
Survival time / Tiempo de sobrevida _____ sec./seg.	
Survival time = (log ₁₀ labeled population - 2) x labeled D-value	
Kill time / Tiempo de muerte _____ min.	
Kill time = (log ₁₀ labeled population + 4) x labeled D-value	
Z-value / Valor Z _____ °C	

System's sensitivity / Sensibilidad del Sistema:

System's sensitivity is determined as the difference between the number of positive indicators after 7-day incubation and false negative indicators (negative by fluorescence readout and visually positive) in relation to the number of 7-day positive indicators. / La sensibilidad del sistema se determina como la diferencia entre el número de incubaciones positivas después de 7 días de incubación y los indicadores falsos negativos (negativos por lectura de fluorescencia y positivos visualmente) en relación al número de indicadores positivos luego de 7 días de incubación.

Sensitivity / (Nº Positives at 7 days) - (Nº False negatives)
Sensibilidad =

(

N

º

Positives at 7 days)
(

N

º

False negatives)
X 100 ≥ 97%

Parameters determined at time of manufacture according to ISO 11138-1: 2017, ISO 11138-3: 2017 and IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant
TERRAGENE® and BIONOVA® are registered trademarks of Terragene S.A. ATCC® is a registered trademark of American Type Culture Collection.

Adrián J. Rovetto
Director Técnico
Technical Director

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Explanation of Symbols

Product is designed for use with Steam sterilization cycles
 Expiration Date
 Batch number
 Manufacturer

Intended Use Table		
	Autoclave/Steam Cycles	
Model	Gravity Displacement	135 °C
	132 °C	
	10 minutes	10 minutes
	15 minutes	
	25 minutes	
BT224	Dynamic Air Removal (Vacuum Assist)	
	132 °C	135 °C
	4 minutes	3 minutes
	Fluorescence Read Time	pH Color Change
	20 minutes	48 hours

EN Biological Indicators

For Steam sterilization

Indicators by Use

United States

Terragene® Bionova® SCBI (BT224) is a self-contained biological indicator inoculated with viable 10⁸ *Geobacillus stearothermophilus* bacterial spores and is intended for monitoring the efficacy of steam sterilization processes. On each Terragene® Bionova® SCBI is a chemical process indicator that changes color from pink to brown when exposed to Steam. See Intended Use Table for reference.

Outside the United States

Terragene® Bionova® BT224 Fluorescence Ultra Rapid Readout Biological Indicators have been designed for quick and easy monitoring of vacuum assisted and gravity displacement steam sterilization cycles at 132-135 °C.

Device Description

Terragene® Bionova® BT224 Fluorescence Ultra Rapid Readout Biological Indicators are single-use Self-Contained Biological Indicators (SCBI) that consist of a polypropylene tube, a spore carrier and a glass ampoule with a culture medium, enclosed with a colored cap. Each tube contains a population of *Geobacillus stearothermophilus* ATCC® 7953 spores inoculated on a spore carrier, a plastic cap with holes and a barrier permeable to Steam. Each BT224 has a Process Indicator on label that changes from pink to brown when exposed to Steam.

Precautions

WARNING: Do not use Terragene® Bionova® BT224 SCBI to control Hydrogen peroxide, EO, Dry Heat, Formaldehyde or other sterilization processes different from the specified. Do not reuse SCBI.

WARNING: Place one or more SCBI in sterilizing hard-to-reach areas to ensure all areas of the chamber are sterilized. Please evaluate all load configurations to ensure ALL hard-to-reach areas have been identified, and place a SCBI in each of those locations.

WARNING: Do not reuse the sterilizer until the SCBI test result is negative.

Instructions for Use

1. Identify the Terragene® Bionova® BT224 SCBI by writing the sterilizer number (in case of having more than one), load number and processing date on the label. 2. Pack the SCBI along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place the package in those areas which are considered most inaccessible for the sterilizing agent (e.g., the center of the load and areas near the door).

3. Sterilize as usual.
4. After the sterilization process has finished, open the sterilizer door, wait five minutes and remove the SCBI from the package. **CAUTION:** Wear safety glasses and gloves when removing the Terragene® Bionova® BT224 SCBI from the sterilized package. **WARNING:** Do not crush or handle the SCBI excessively, since this might cause the glass ampoule to burst. Let the SCBI cool down until it reaches room temperature.

5. Check the Process Indicator on SCBI label. A color change to brown indicates that the medium reaches the base of the tube and soaks the entire area, which does not evidence the process effectiveness to achieve sterility. If the Process Indicator color has not changed, check the sterilization process.

6. Press the lid to seal the tube. Crush the ampoule contained in the SCBI with an individual ampoule crusher or with the ampoule crusher placed within the incubator's incubation area. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the SCBI in the incubator. **IMPORTANT:** Use a non-sterilized SCBI as a positive control at least once per day, when a sterilization cycle is run. The positive control ensures that correct incubation conditions were met; capability of medium to promote rapid growth; viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals and proper functioning of Terragene® Bionova® Auto-Reader incubators. Both, the positive control indicator and the processed indicator should belong to the same batch.

7. Incubate the processed indicator and the positive control indicator in the appropriate Terragene® Bionova® Auto-Reader Incubators for a maximum of 20 minutes at (60 ± 2) °C for ultra rapid readout. **NOTE:** Holding time between sterilization and incubation should not exceed a 7-day period. Fluorescence detection by the reader (excitation 340-380 nm / emission 455-465 nm) means a failure in the sterilization process. If no fluorescence is detected at 20-minute incubation, the result is negative. The positive control must give positive fluorescence readout. It is good practice to incubate a positive control for a visual color change.

Record the results and discard the SCBI immediately, as it is indicated below.

Ultra Rapid readout: 20 minutes

The Ultra Rapid Readout must be carried out in the appropriate Terragene® Bionova® Auto-Reader Incubators. Fluorescence is emitted after the reader stimulates the spore carrier with UV light. Final readout of negative results is readily available after 20 minutes of incubation. Fluorescence readout is an indirect measure of the viability of *Geobacillus stearothermophilus* spores after the sterilization process (positive results). Furthermore, a failure in the sterilization process can also become evident by culture medium color change. Due to the high sensitivity of the fluorescence results at 20 minutes, conventional incubation for color change of Terragene® Bionova® SCBI BT224 is not an advantage.

Visual confirmation: 48 hours

Optionally, you can perform a visual color-change confirmation after a 48-hour incubation. If the sterilization process has not been successful, culture medium will turn yellow during incubation at (60 ± 2) °C, indicating the presence of living spores. If sterilization was successful, culture medium will remain purple after incubation. The positive control must show a color change from purple to yellow for results to be valid.

Readout time: 7 days

A 7-day readout is optional and may be routinely performed. This is an initial validation of the 20-minute readout. Fluorescence results are compared to the 7-day visual readouts. **NOTE:** If 7-day readout is performed, a humidified environment will be required to prevent medium from drying out.

Monitoring frequency

Follow the policies and procedures with the monitoring frequency specified by the professional associations and/or standards corresponding to your country. As recommended practice, and to provide optimal patient safety, Terragene recommends that each sterilization load be monitored with the appropriate biological indicator.

Storage

Store in a dark place at temperature between 10-30 °C and 30-90% relative humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Shelf life

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date. Chemical Process Indicator on SCBI label have an expiration date of 2 years when used as part of SCBI.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicators can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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