# BT23 Biological Indicator

Rev. 10 / 09.2021

## Spore Ampoules

BT23/4 | BT23/5 | BT23/6

Purple = Sterile

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## Quality certification Certificado de calidad Bionova® BT23/

Steam sterilization / Esterilización por Vapor  
*Geobacillus stearothermophilus* ATCC 7953

LOT

Population / Población \_\_\_\_\_ CFU / UFC

D - value / Valor D (121 °C) \_\_\_\_\_ min.

Survival time / Tiempo de sobrevivencia \_\_\_\_\_ min.  
Survival time =  $(\log_{10} \text{ labeled population} - 2) \times \text{labeled D-value}$

Kill time / Tiempo de muerte \_\_\_\_\_ min.  
Kill time =  $(\log_{10} \text{ labeled population} + 4) \times \text{labeled D-value}$

D - value / Valor D (132 °C) \_\_\_\_\_ sec./seg.

Survival time / Tiempo de sobrevivencia \_\_\_\_\_ min.  
Survival time =  $(\log_{10} \text{ labeled population} - 2) \times \text{labeled D-value}$

Kill time / Tiempo de muerte \_\_\_\_\_ min.  
Kill time =  $(\log_{10} \text{ labeled population} + 4) \times \text{labeled D-value}$

D - value / Valor D (135 °C) \_\_\_\_\_ sec./seg.

Survival time / Tiempo de sobrevivencia \_\_\_\_\_ sec./seg.  
Survival time =  $(\log_{10} \text{ labeled population} - 2) \times \text{labeled D-value}$

Kill time / Tiempo de muerte \_\_\_\_\_ min.  
Kill time =  $(\log_{10} \text{ labeled population} + 4) \times \text{labeled D-value}$

Z-value / Valor Z \_\_\_\_\_ °C

Parameters determined at time of manufacture according to ISO 11138-1:2017, ISO 11138-3:2017 and IRAM 37102:1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant.  
ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto  
Director Técnico  
Technical Director

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## EN Biological Indicators For Steam sterilization

**Composition**  
Bionova® BT23 Spore Ampoules consists of a glass ampoule containing *Geobacillus stearothermophilus* ATCC 7953 spore population dispersed in a culture medium especially designed to change its color when microbial growth occurs.

**Product description**  
Bionova® BT23 Self-Contained spore ampoules has been designed for monitoring Steam sterilization processes at temperatures between 121°C and 135°C. If the sterilization process was not effective, the indicator medium will turn to yellow after incubation between 55-62 °C, thus indicating the presence of living *Geobacillus stearothermophilus* spores.

If the sterilization process was effective, the indicator medium will remain the original color after incubation. The final readout should be performed after 48 hours of incubation between 55-62 °C.

### Precautions

Do not use Bionova® BT23 Biological Indicators for monitoring Ethylene Oxide, Dry Heat, Chemical Vapors, Radiation or any other sterilization process different from Steam.  
Do not re-use biological indicators. Must be discarded after use.  
Protect from shocks, falls and sudden temperature changes.  
Check the color of the indicator before starting a sterilization process. It must be purple. Discard it otherwise. The negative control vial can be used to compare the color.  
Do not use an indicator after the expiration date indicated in the package.  
The indicator must be processed by a trained healthcare professional.  
Do not re-use the sterilizer until biological indicator test results are negative (processed indicator culture medium remains its original color).

### Storage

The ampoules should be stored away from light, refrigerated at temperatures between 2-8 °C and 30-80% relative humidity.  
Storage at room temperature (up to 25 °C) is possible for a period no longer than 15 days. **IMPORTANT:** If the product is not stored at recommended temperatures, the stability of the vial content might be affected.  
Do not freeze.  
Do not store close to other sterilization agents or chemical products.

### Instructions for use

1. Take the Bionova® BT23 Ampoules out of the refrigerator and leave them stabilize at room temperature up to ten minutes.
2. Place the Bionova® BT23 Ampoules inside a container with the liquid being sterilized or in a container with liquid similar to that containing the liquid to be sterilized. This container should be the same size and have the same volume as the one containing the liquid being sterilized.
3. Sterilize as usual.
4. Once the sterilization process has finished, open the sterilizer door, wait five minutes and remove the biological indicator from the container. **WARNING:** Wear safety glasses and gloves when removing the Bionova® BT23 Biological Indicator from the sterilized container.
5. Incubate the processed biological indicator, a negative control and a non-sterilized biological indicator as a positive control, during 48 hours at 55-62 °C. The positive control ensures that the incubation conditions were met. The negative control, not being exposed to the sterilization process, can be used as a reference color to compare positive and negative results. Conversely, when used with the ampoules during the process, it ensures that any color change produced immediately after sterilization is the result of the process and not due to germination and/or spore growth.  
Readout should be performed at convenient intervals of 12 hours.  
It is recommended the use of Bionova® IC10/20 Dual Incubator or similar devices that assure an efficient temperature control. **WARNING:** Traditional culture stoves have thermal profiles that differ from the above recommended incubators, thus their use might affect or delay the Bionova® BT23 Biological Indicators result. A color change to yellow of the processed biological indicator medium means a sterilization process failure has occurred. If after 48 hours there is no color change in the processed indicators, a final negative result is obtained, meaning the sterilization process was effective. The positive control indicator should turn to yellow for the result to be valid and the negative control should remain the same color.  
Record the results and discard the positive indicators immediately as it is shown below.

### Disposal

Discard Bionova® BT23 Biological Indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be sterilized in gravity air displacement steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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