Spore Ampoules

Quality certification Certificado de calidad Bionova® BT31

Dry heat sterilization / Esterilización por calor seco
Bacillus atrophaeus ATCC® 9372

LOT

Population / Población _____ CFU / UFC

D-Value / Valor D (160 °C) _____ min.

Survival time / Tiempo de sobrevivida _____ min.
Survival time = not less than (log₁₀ labeled population - 2) x labeled D-value

Kill time / Tiempo de muerte _____ min.
Kill time = not more than (log₁₀ labeled population + 4) x labeled D-value

Z-value / Valor Z _____ °C

Parameters determined at time of manufacture ISO 11138-1:2017, ISO 11138-4:2017 and IRAM 37102-1:1999. The shown values are reproducible only under the same conditions under which they were determined.

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ISO 11138-1 and ISO 11138-4 compliant
ATCC® is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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Explanation of Symbols

-  Product designed for use in dry heat sterilization processes.
-  Batch number.

-  Manufacture date.
-  Expiration date.
-  Manufacturer.

EN Biological Indicators

For dry heat sterilization

Composition

Bionova® BT31 Biological Indicator system consists of a glass ampoule containing a population of *Bacillus atrophaeus* ATCC® 9372 spores embedded in a solid matrix.

Product description

Bionova® BT31 Biological Indicator is designed for the control of dry heat sterilization processes. It is compatible with extremely high temperatures generated during deprogenation cycles and sterilization processes in ovens and dry heat tunnels. The contents of the ampoule should be incubated in the Bionova® MC1030-2 Culture Medium, Bionova® MC1020-2 Culture Medium or a culture medium suitable for the growth of *Bacillus atrophaeus* ATCC® 9372. If the sterilization process was unsuccessful, Bionova® indicator medium will change to yellow after incubation at 37±2 °C, thus indicating the presence of living spores. If the sterilization process was successful, the indicator medium will remain blue (MC1030-2) or pink/red (MC1020-2), after 72-hour incubation at 37±2 °C.

Precautions

Do not use the Bionova® BT31 Biological Indicator to control steam, ethylene oxide, chemical vapors, radiation, or other sterilization processes. Do not reuse. Protect from shocks, falls, sudden temperature changes, etc.

Instructions for use

- Place Bionova® BT31 ampoule next to the material to be sterilized, in a suitable package, and according to the recommended sterilization practices. Put the package in those areas that you consider a priori most inaccessible to the sterilizing agent (dry heat). Generally, a challenging area is the center of the load or near the door.
- Sterilize in the usual way.
- Once the sterilization process finishes, open the sterilizer door, wait for the necessary time for the material to cool down, and remove the biological indicator from the package.

- Let the biological indicator cool down until it reaches room temperature.
- Under sterile conditions (laminar flow or sterile environment):
 - Break the vial at the indicated mark (below the blue dot on the vial neck).
 - Pour the contents of the ampoule into the Bionova® MC1030-2 Culture Medium, Bionova® MC1020-2 Culture Medium or any compatible medium to reveal the result.
 - Mix by inverting, making sure to close the culture medium cap tightly.
- IMPORTANT:** Use gloves and a mask when transferring the solid matrix to the culture medium tube. Ensure that the ampoule contents do not come into contact with the plastic tube rim or any surface.
- Incubate at 37±2 °C.
- NOTE:** The use of devices that ensure efficient temperature control is recommended.
- IMPORTANT:** Use a biological indicator not subjected to the sterilization process as a positive control, at least once per day, when a sterilization cycle is performed. The incubation of a positive control biological indicator is important to check the viability of the spores in it and to ensure that incubation conditions are adequate.
- Incubate the processed biological indicator and the indicator used as a positive control for a maximum of 72 hours at 37±2 °C. Make observations conveniently every 10 hours. The color change to yellow of the Bionova® growth indicator medium indicates a failure in the sterilization process. If after 72 hours no color change is observed in the processed indicators, the result is negative (the sterilization process was efficient). The color of the indicator used as a positive control must change from blue (MC1030-2) or pink/red (MC1020-2) to yellow for results to be valid.
- Record the results and discard them immediately, as it is indicated below.

WARNING: Do not use the load until the biological indicator result is negative (the medium containing the processed indicator does not manifest a color change).

Storage

Store away from light and at a temperature between 10-30 °C, relative humidity between 30-80 %. Do not freeze. Do not store near sterilizing agents or other chemical products.

Disposal

Discard biological indicators according to your country's healthcare and safety regulations. Culture medium positive tubes can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, at 132 °C for 15 minutes or at 134 °C for 10 minutes; or in a dynamic air removal Steam sterilizer at 132 °C for 4 minutes or at 135 °C for 3 minutes.

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