BT40

Rev. 12 / 10.2020

Spore Strips

BT40/2, BT40/3, BT40/4, BT40/5, BT40/6, BT40/7, BT40/8

Quality certification Certificado de calidad Bionova® BT40/

Ethylene Oxide and Dry Heat Sterilization
Esterilización por Óxido de Etileno y Calor Seco
Bacillus atrophaeus ATCC 9372

LOT

Population / Población _____ CFU / UFC

EO

D-value / Valor D _____ min.
(54°C, 60% RH, 600 mg EO / l)Survival time / Tiempo sobrevivida _____ min.
Survival time = (log₁₀ labeled population - 2) x labeled D-valueKill time / Tiempo de muerte _____ min.
Kill time = (log₁₀ labeled population + 4) x labeled D-value

DRY

D-value / Valor D (160 °C) _____ min.

Survival time / Tiempo sobrevivida _____ min.
Survival time = (log₁₀ labeled population - 2) x labeled D-valueKill time / Tiempo de muerte _____ min.
Kill time = (log₁₀ labeled population + 4) x labeled D-value

Z-value / Valor Z _____ °C

Parameters determined at time of manufacture according to ISO 11138-1: 2017, ISO 11138-2: 2017, ISO 11138-4: 2017 and IRAM 37102: 1999 (Parts 1 and 2) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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EN Spore Strips

For Ethylene Oxide and Dry Heat sterilization

Composition

Each Bionova® BT40 envelope consists of a paper strip soaked with *Bacillus atrophaeus* ATCC 9372 spore population.

Product Description

The envelope containing Bionova® BT40 Spore Strip is specifically designed for the monitoring of Ethylene Oxide (EO) sterilization processes or for Dry Heat sterilization processes. Spore strips should be incubated in Bionova® MC1030-2, Bionova® MC1020-2 Growth Medium tube or in appropriate growth medium for *Bacillus atrophaeus* ATCC 9372 growth. If the sterilization process was not successful, Bionova® MC1030-2 or Bionova® MC1020-2 growth medium will turn to yellow after incubation at 37±2 °C, thus indicating the existence of living *Bacillus atrophaeus* spores on the strip.

If the sterilization process was successful, MC1030-2 or MC1020-2 Indicating Medium will remain its original color. A final readout

should therefore be performed after 48 hours of growth medium incubation at 37±2 °C.

Precautions

Do not use Bionova® BT40 Spore Strip envelopes to control Steam sterilization cycles, radiation or other sterilization processes different from EO and/or Dry Heat sterilization.
Do not reuse spore strip envelopes.

Storage

Best stored in original package under the next conditions: Temperature between 10-30 °C, 30-80 % relative humidity. Do not freeze. Do not store these biological indicators near sterilizing agents or other chemical products.

Instructions for use

- Identify Bionova® BT40 envelope by writing on its back the sterilizer number (in case of having more than one), load number and processing date.
- Pack it in appropriate package along with material to be sterilized, according to recommended sterilization practices. Place the envelope in those areas you a priori consider more inaccessible for sterilizing agent. A typical problematic area is the load center.
- Sterilize as usual.
- After the sterilization process has finished, you can:
 - Remove the Bionova® BT40 envelope from the package for processing and incubation.
 - For Ethylene Oxide sterilization, you can aerate the package containing the biological indicator with the rest of the load or remove the biological indicator from the sterilization package for incubation in MC1030-2, MC1020-2 culture medium or any other appropriate growth medium prior to material aeration.
- Tear up the envelope at one end very carefully under sterility conditions (eg. laminar flow cabinet) and transfer spore strip with a sterile clamp to MC1030-2, MC1020-2 Growth Medium tube, or to any other appropriate growth medium. **IMPORTANT!** Use latex gloves and chin strap when transferring spore strips from the envelope to growth medium tube. **WARNING!** Avoid spore strip contact with any kind of surface including the outer part of growth medium tube and of Bionova® BT40 envelope.

IMPORTANT: Use a non-sterilized spore strip as a positive control at least once per day, when a sterilization cycle is run. Both, processed spore strip and the strip used as positive control must belong to the same batch. Positive control guarantees that the incubation was carried out under appropriate conditions.
- Incubate processed spore strips along with the strip used as positive control during a maximum of 48 hours at 37±2 °C. It is advisable to make observations every 10 hours. Color change of MC1030-2 or MC1020-2 Growth medium to yellow means a failure on sterilization process. If after 48 hours no colour change is visible on MC1030-2 or MC1020-2 Growth Medium, the result is negative (i.e. the sterilization process was successful). MC1030-2 or MC1020-2 Indicating Medium used for incubating the spores strips which have not been sterilized (positive control) should turn to yellow in order for the results to be valid. Record the positive ones and dismiss them immediately as indicated below.

Shelf life

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date.

Disposal

Discard Growth Medium tubes already used to analyse Bionova® BT40 Spore Strips growth according to your country's healthcare and safety regulations. Positive MC1030-2 or MC1020-2 Growth Medium tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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Terragene S.A.