Spore Strips

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Parameters determined at time of manufacture by calculations adapted from ISO 11138-1: 2017 e IRAM 37102-1: 1999 standards. The showed values are reproducible only under the same conditions under which it were determined.

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
ISO and USP Compliant.
ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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EN Spore Strips

For Liquid Peracetic Acid High Level Disinfection processes

Composition

Each Bionova® envelope contains a polyethylene fibers strip soaked with *Geobacillus stearothermophilus* ATCC 7953 spore population.

Product description

The envelope containing Bionova® Spore Strip is specifically designed to control Liquid Peracetic Acid High Level Disinfection processes. Spore strips should be incubated in MC1020-2 Bionova® Growth Medium tube or other medium appropriate for *Geobacillus stearothermophilus* ATCC 7953 growth.

If disinfection process was not successful, MC1020-2 Growth Medium will turn to yellow after incubation at 60±2 °C, thus indicating the existence of living *Geobacillus stearothermophilus* spores in spore strips.

If disinfection process was effective, MC1020-2 Growth Indicating medium will remain its original color. Final reading should be done after a 24-hour incubation at 60±2 °C of the growth medium containing the spore strip.

Precautions

Do not use spore strip envelopes to control EO, Dry Heat, Radiation or any sterilization or disinfection process other than the one for which the indicator was designed.

Do not reuse spore strip envelopes.

Do not reuse the equipment until spore strip growth result turns negative (MC1020-2 Growth Medium containing processed spore strip remains its original colour).

Storage

Store in a dark place and at temperatures between 10-30 °C, 30-80 % relative humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

Directions for use

1. With the reusable silicone strip provided along with the product, locate the envelope with spore strip inside the High Disinfection equipment, accompanying the load to be disinfected with PAA. **WARNING:** The indicator must be completely submerged in the washing solution to work properly.

2. Disinfect as usual.

3. After high disinfection process has finished, remove Bionova® envelope for processing and incubation.

4. Open the envelope at the end marked with "PEEL OFF" very carefully and in sterility conditions (eg. laminar flow cabinet), and transfer the Spore Strip using a sterile clamp to MC1020-2 Growth Medium tube, or any other appropriate growth medium.

IMPORTANT: Use latex gloves and chin strap when transferring spore strips from the envelope to the growth medium tube. **WARNING:** Avoid spore strip contact with any kind of surface including the outer part of Growth Medium tube and Bionova® envelope.

IMPORTANT: Use a spore strip which has not been subjected to disinfection process as a positive control whenever you incubate a processed strip. Positive control guarantees that the incubation was carried out under appropriate conditions.

5. Incubate processed spore strips along with the strip used as positive control during a maximum of 24 hours at 60±2 °C. It is advisable to make observations every 10 hours. Growth medium indicator colour turning to yellow means a failure in disinfection process. If after 24 hours no colour change is visible on MC1020-2 Growth Medium, the result is negative (i.e. the disinfection process was successful). MC1020-2 Growth Medium used for incubating the spores strips which have not been disinfected (positive control) should turn to yellow in order for the results to be valid. Record the positive samples and dispose them immediately as indicated below.

Shelf life

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date.

Disposal

Dispose the growth medium tubes used to evaluate Bionova® Spore Strips growth according to your county's healthcare and safety regulations. Positive MC1020-2 Growth Medium tubes can be sterilized in gravity air displacement Steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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