BT50 Biological Indicator

Rev. 17 / 11.2020

Spore Strips

BT50/4 | BT50/5 | BT50/6 | BT50/7

Quality certification Certificado de calidad Bionova® BT50/X

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Parameters determined at time of manufacture according to ISO 11138-1: 2017, ISO 11138-3: 2017, ISO 11138-5:2017 and IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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EN Spore Strip For Steam and Formaldehyde sterilization

Composition
Each Bionova® BT50 envelope consists of a paper strip soaked with *Geobacillus stearothermophilus* ATCC 7953 spore population.

Product description
The envelope containing Bionova® BT50 Spore Strip is specifically designed to control Steam 121 °C-135 °C (Steam) and Formaldehyde (Form) sterilization processes. Spore strips should be incubated in Bionova® MC20-2 or MC1020-2 Growth Medium tube or in appropriate Growth medium for *Geobacillus stearothermophilus* ATCC 7953 growth.
If the sterilization process was not effective, MC20-2 or MC1020-2 Growth Medium will turn to yellow after incubation between 55-62°C, thus indicating the existence of living *Geobacillus*

stearothermophilus spores on the strips.
If the sterilization process was effective, MC20-2 or MC1020-2 Indicating Medium will remain its original color. A final readout should therefore be performed after 24 hours of growth medium incubation between 55-62°C for Steam sterilization processes control, or 48 hours between 55-62°C for Formaldehyde sterilization processes control.

Precautions
Do not use Bionova® BT50 Spore Strip envelopes to control EO, Dry Heat sterilization, Radiation or other sterilization processes different from those for which were designed.
Do not reuse spore strip envelopes.
Do not reuse sterilizer until spore strip growth result is negative (MC20-2 or MC1020-2 Growth Medium containing processed spore strip remains its original color).

Storage
Best stored in original box under the next conditions: Temperature between 10-30 °C, 30-80 % relative humidity.
Do not freeze.
Do not store these biological indicators near sterilizing agents or other chemical products.

Instructions for use
1. Identify Bionova® BT50 envelope by writing on its back the sterilizer number (in case of having more than one), load number and processing date.
2. Pack it in an appropriate package along with material to be sterilized, according to recommended sterilization practices. Place the envelope in those areas you a priori consider more inaccessible for sterilizing agent. A typical problematic area is the load center.
3. Sterilize as usual.
4. After sterilization process has finished, remove Bionova® BT50 envelope from the package for processing and incubation.
5. Tear up the envelope at one end very carefully under sterility conditions (eg. laminar flow cabinet) and transfer the spore strip with a sterile clamp to MC20-2 or MC1020-2 Growth Medium tube or any other appropriate growth medium. **WARNING:** Avoid spore strip contact with any kind of surface including the outer part of Growth Medium tube and of Bionova® BT50 envelope.
IMPORTANT: Use latex gloves and a mask when transferring spore strips from the envelope to growth medium tube.
Use a spore strip which has not been subjected to the sterilization process as a positive control each time a processed strip is incubated. Both, processed spore strip and the strip used as positive control must belong to the same batch. Positive control guarantees that the incubation was carried out under appropriate conditions.
6. Incubate processed spore strips along with the strip used as a positive control between 55-62°C during a maximum of 24 hours for Steam or 48 hours for Formaldehyde. It is advisable to make observations every 10 hours.
Color change of MC20-2 or MC1020-2 indicating medium to yellow means a failure of the sterilization process. If after 24 hours (Steam) or 48 hours (Form) no color change is visible on MC20-2 or MC1020-2 growth medium, the result is negative (i.e. the sterilization process was effective). MC20-2 or MC1020-2 Growth Medium used for incubating the spores strips which have not been sterilized (positive control) should turn to yellow for the results to be valid. Record the positive ones and dismiss them immediately as indicated below.

Shelf life
Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date.

Disposal
Discard Growth Medium tubes after use to analyse Bionova® BT50 Spore Strips growth according to your country's healthcare and safety regulations. Positive MC20-2 or MC1020-2 Growth Medium tubes can be sterilized in gravity air displacement Steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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