BT70

Rev. 8 / 11.2020

Spore Strips

BT70/5, BT70/6, BT70/7

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Population / Población _____ CFU/UFC

D-value / Valor D _____ kGy

Parameters determined at time of manufacture according to ISO 11138-1:2017 and IRAM 37102-1:1999 standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant
ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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EN Spore Strips

For Gamma Radiation sterilization

Composition

Each Bionova® BT70 envelope consists of a paper strip soaked with a *Bacillus pumilus* ATCC 27142 spore population.

Product description

The envelope containing Bionova® BT70 Spore Strip is specifically designed to control Gamma Radiation sterilization processes. Spore strips should be incubated in Bionova® MC1020-2 Growth Medium tube or in appropriate growth medium for *Bacillus pumilus* growth. If the sterilization process was not successful, MC1020-2 Growth Medium will turn to yellow after incubation at 37±2 °C, thus indicating the existence of living *Bacillus pumilus* spores on the strips. If the sterilization process was successful, MC1020-2 Indicating Medium will remain its original color. A final readout should therefore be performed after 48 hours of growth medium incubation at 37±2 °C.

Precautions

Do not use Bionova® BT70 Spore Strip envelopes to control Steam sterilization cycles, Chemical Vapors, Ethylene Oxide or other sterilization processes different from those for which they were designed. Do not reuse spore strip envelopes. Do not reuse the sterilizer until the spore strip growth result is negative (MC1020-2 Growth Medium containing processed spore strips remains its original color).

Storage

Store in a dark place and at temperatures between 10-30 °C, 30-80 % relative humidity. Do not freeze. Do not store near sterilizing agents or other chemical products.

Directions for use

1. Identify Bionova® BT70 envelope by writing on its back the sterilizer number (in case of having more than one), load number and processing date.
2. Pack it in appropriate package along with material to be sterilized, according to recommended sterilization practices. Place the envelope in those areas you consider a priori more inaccessible for the sterilizing agent. A typical problematic area is the load center.
3. Sterilize as usual.

4. After sterilization process has finished remove Bionova® BT70 envelope from the sterilization package used for processing and incubation.

5. Tear up the envelope at one end very carefully under sterility conditions (eg. laminar flow cabinet) and transfer the spore strip with a sterile clamp to MC1020-2 Growth Medium tube, or to any other appropriate growth medium. **IMPORTANT:** Use latex gloves and chin strap when transferring spore strips from the envelope to MC1020-2 Growth Medium tube. **WARNING:** Avoid spore strip contact with any kind of surface including the outer part of growth medium tube and of Bionova® BT70 envelope.

IMPORTANT: Use a spore strip which has not been subjected to the sterilization process as a positive control each time a processed strip is incubated. Both, processed spore strip and the strip used as positive control must belong to the same batch. Positive control guarantees that the incubation was carried out under appropriate conditions.

6. Incubate processed spore strips along with the strip used as positive control during a maximum of 48 hours at 37±2 °C. It is advisable to make observations every 10 hours. Color change of MC1020-2 Growth Medium to yellow means failure of the sterilization process. If after incubation, no color change is visible on MC1020-2 Growth Medium, the result is negative (i.e. the sterilization process was successful). MC1020-2 Indicating Medium used for incubating the spores strips which have not been sterilized (positive control) should turn to yellow in order for the results to be valid.

Record the positive results and discard them as indicated below.

Shelf life

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date.

Disposal

Discard Growth Medium tubes after used to analyse Bionova® BT70 spore strips growth according to your country's healthcare and safety regulations. Positive MC1020-2 Growth Medium tubes can be sterilized in gravity air displacement Steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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