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Composición

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Composition

Each tube contains a population of *Bacillus atrophaeus* ATCC 9372 spores soaked on a carrier. It also has a blue growth indicator medium contained in the glass ampoule.

Product description

Bionova® BT80 Biological Indicator has been specifically designed for monitoring of Moist Heat decontamination processes, such as microwave and radiofrequency irradiation processes intended for medical waste disposal treatment. If the decontamination process was not successful, the indicator medium would turn from blue to yellow after incubation at 37±2 °C, thus indicating the presence of living *Bacillus atrophaeus* spores. If the decontamination process was successful, indicator medium would remain blue after incubation. The final readout should be carried out after 48 hours of incubation at 37±2 °C.

WARNING: Do not use Bionova® BT80 Biological Indicators for monitoring Steam sterilization cycles, Dry Heat, Chemical Vapor, Ethylene Oxide, Radiation or any sterilization processes other than Moist Heat decontamination. Do not re-use biological indicators.

Storage

Store it in a dark place at temperatures between 10-30 °C, 30-80 % Relative Humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Instructions for use

- Identify the Bionova® BT80 Indicator Tube by writing the machine number (in case there is more than one), load number, and processing date on the indicator label.
 - Place the biological indicator along with the items to be decontaminated in an appropriated package according to recommended decontamination practices. Place this package in those areas which a priori you consider most inaccessible for the decontamination agent (Moist Heat). Generally a problematic area is the center of the load. Proceed as usual.
 - After the decontamination process has finished, allow the load to cool down to room temperature and remove the biological indicator from the load for its incubation.
 - Crush the glass ampoule contained in the biological indicator and incubate at 37±2 °C.
- IMPORTANT:** Use a non-processed biological indicator as positive control each time a processed indicator is incubated. The positive control ensures that correct incubation conditions were met.
- Incubate the processed biological indicator and the indicator used as positive control for a maximum of 48 hours at 37±2 °C. Reading should be performed at convenient intervals of 10 hours.
 - A color change to yellow of the growth indicator medium means a decontamination process failure has occurred. If after 48 hours there is no color change in the processed indicators, a final negative result is obtained, thus the decontamination process was successful. The positive control indicator should show a color change to yellow for the results to be valid. Record the positive biological indicators and discard them immediately as shown below.
- WARNING:** Do not reuse the decontamination machine until the biological indicator test results are negative (process indicator remain the original color).

Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

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Survival time = not less than (log₁₀ labeled population - 2) x labeled D-value

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Kill time = not more than (log₁₀ labeled population + 4) x labeled D-value

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Parameters determined at time of manufacture according to ISO 11138-1:2006 and IRAM 37102-1:1999 standards. The shown values are reproducible only under the same conditions under which they were determined.

ISO and USP Compliant.

ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrian J. Rovetto
Director Técnico
Technical Director

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