Spore Coupon

BT93/5 | BT93/6

Quality certification Certificado de calidad Bionova® BT93/

VH202

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Parameters determined at time of manufacture according to ISO 11138-1: 2017 and IRAM 37102-1: 1999 standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant
ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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EN Spore Coupon

For Plasma or Vaporized Hydrogen Peroxide sterilization

Composition

Each Bionova® BT93 envelope contains a stainless steel coupon soaked with *Geobacillus stearothermophilus* ATCC 7953 spore population.

Product description

The envelope containing Bionova® BT93 Spore Coupon is specifically designed to control Plasma or Vaporized Hydrogen Peroxide sterilization processes. Spore coupons should be incubated in MC1020-2 Bionova® Growth Medium tube or other medium appropriate for *Geobacillus stearothermophilus* ATCC 7953 growth.

If sterilization process was not effective, MC1020-2 Growth Medium will turn to yellow after incubation between 55-62 °C, thus indicating the existence of living *Geobacillus stearothermophilus* spores in the spore coupon.

If sterilization process was effective, MC1020-2 Growth Indicating medium will remain its original color. Final readout should be done after a 24-hour incubation between 55-62 °C of the growth medium containing the spore coupon.

Precautions

Do not use spore coupon envelopes to control EO, Dry Heat, Radiation or any sterilization process other than the one for which the indicator was designed.

Do not re-use spore coupon envelopes.

Do not re-use sterilizer until spore coupon growth result turns negative (MC1020-2 Growth Medium containing processed spore coupon remains its original color).

Storage

Store at temperatures between 10-30 °C, 30-80 % relative humidity. Do not freeze. Do not store near sterilizing agents or other chemical products.

Instructions for use

1. Pack the envelope in an appropriate package along with items for sterilization, according to recommended sterilization practices. Place it in those areas you consider a priori more inaccessible to the sterilizing agent. The load center is generally a problematic area.

2. Sterilize as usual.

3. After sterilization process has finished, remove Bionova® BT93 envelope from the package for processing and incubation.

4. Tear the envelope at one end very carefully in sterility conditions (eg. laminar flow cabinet) and transfer the spore coupon using a sterile clamp to MC1020-2 Growth Medium

tube, or any other appropriate growth medium. **WARNING:** Avoid spore coupon contact with any kind of surface including the outer part of growth medium tube and Bionova® BT93 envelope. **IMPORTANT:** Use latex gloves and a mask when transferring spore coupons from the envelope to the Growth Medium tube.

Use a spore coupon which has not been subjected to sterilization process as a positive control whenever you incubate a processed coupon. Both, processed spore coupon and the coupon used as positive control must belong to the same batch. Positive control guarantees that the incubation was carried out under appropriate conditions.

5. Incubate processed spore coupon along with the coupon used as positive control during a maximum of 24 hours between 55-62 °C. It is advisable to make observations every 10 hours. Color change of MC1020-2 to yellow means a failure in sterilization process. If after 24 hours no color change is visible on MC1020-2 Growth Medium, the result is negative (i.e. the sterilization process was effective. MC1020-2 growth medium used for incubating the spore coupon which have not been sterilized (positive control) should turn to yellow for the results to be valid.

Record the positive results and dispose biological indicators immediately as indicated below.

Disposal

Discard the growth medium tubes used to evaluate Bionova® BT93 spore coupon growth according to health regulations of your country. Positive MC1020-2 Growth Medium tubes can be sterilized in gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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