
BT98 Biological Indicator

Rev. 3 / 07.2023

BIONOVA[®]

Hyper Rapid Readout Fluorescence System

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Quality certification Certificado de calidad Bionova® BT98

Plasma or Vaporized Hydrogen Peroxide sterilization / Esterilización por Plasma o Vapor de Peróxido de Hidrógeno
Geobacillus stearothermophilus ATCC® 7953

LOT

Heat shock population / Población _____ CFU / UFC
D - value / Valor D _____ sec./seg.
(2.0 mg/L VH₂O₂, 50 °C)

Survival time / Tiempo de Sobrevida _____ min.
Survival time = (log₁₀ labeled population - 2) / labeled D-value

Kill time / Tiempo de muerte _____ min.
Kill time = (log₁₀ labeled population + 4) x labeled D-value

System’s sensitivity / Sensibilidad del Sistema:

System’s sensitivity is determined as the difference between the number of positive indicators after 7-day incubation and false negative indicators (negative by fluorescence readout and visually positive) in relation to the number of 7-day positive indicators. / La sensibilidad del sistema se determina como la diferencia entre el número de indicadores positivos después de 7 días de incubación y los indicadores falsos negativos (negativos por lectura de fluorescencia y positivos visualmente) en relación al número de indicadores positivos luego de 7 días de incubación.


Sensibility / = (N° Positives at 7 days) - (N° False negatives) X 100 = 9% Sensitivity
(N° Positives at 7 False)

Parameters determined at time of manufacture according to ISO 11138-1: 2017 and IRAM 37102-1: 1999 standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant.

ATCC® is a registered trademark of American Type Culture Collection.
Terragene® and *Bionova*® are registered trademarks of Terragene S.A.
STERRAD® is a registered trademark of Advanced Sterilization Products.
V-PRO® is a registered trademark of Stenis Corporation.

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Explanation of Symbols

Product designed for use with Vaporized Hydrogen Peroxide sterilization cycles
 Manufacture Date

Batch number
 Expiration Date

Manufacturer

Model	Vaporized Hydrogen Peroxide
	Cycles
BT98	Sterrad® 100S Sterilization System
	Sterrad® NX Sterilization System (Standard and Advanced Cycles)
	Sterrad® NX with ALLClear® Technology Sterilization System (Standard and Advanced Cycles)
	Sterrad® 100NX Sterilization System (Standard, Flex, Express and Duo Cycles)
	Sterrad® 100NX with ALLClear® Technology Sterilization System (Standard, Flex, Express and Duo Cycles)
	V-Pro® S2 Low Temperature Sterilization System (Fast, Non Lumen, Lumen and Flexible Cycles)
V-Pro® maX 2 Low Temperature Sterilization System (Fast-Non Lumen, Non Lumen, Lumen and Flexible Cycles)	
Fluorescence Read Time	pH Color Change
5 minutes	5 days

EN Biological Indicators

For plasma or vaporized hydrogen peroxide sterilization

Indicators for Use

United States
Terragene® Bionova® Hyper Biological Indicator (BT98) is a self-contained biological indicator (SCBI) inoculated with a minimum of 10⁸ viable *Geobacillus stearothermophilus* bacterial spores and is intended for monitoring the efficacy of plasma or vaporized hydrogen peroxide sterilization processes. BT98 has hyper rapid readout at 5 minutes at 60 °C. See intended use systems for reference.

Outside the United States

Terragene® Bionova® BT98 fluorescence hyper rapid readout biological indicators have been designed for quick and easy monitoring of plasma or vaporized hydrogen peroxide sterilization cycles.

Device description

Terragene® Bionova® BT98 fluorescence hyper rapid readout Biological Indicators are single-use self-contained biological indicators (SCBIs) that consist of a polypropylene tube, a spore carrier and a glass ampoule with a culture medium, enclosed with a colored cap. Each tube contains a population of *Geobacillus stearothermophilus* ATCC® 7953 spores inoculated on a spore carrier, a plastic cap with holes and a barrier permeable to plasma or vaporized hydrogen peroxide. Each BT98 has a process indicator on the label that changes from purple to green when exposed to hydrogen peroxide.

Precautions
Do not use the Terragene® Bionova® BT98 SCBI to control EO, dry heat, formaldehyde or other sterilization processes different from those specified. Do not reuse the SCBI.
WARNING: Place one or more SCBIs in hard-to-reach sterilizing areas to ensure all areas of the chamber are sterilized. Please evaluate all load configurations to ensure ALL hard-to-reach areas have been identified, and place a SCBI in each of those locations.
WARNING: Do not reuse the sterilizer until the SCBI test result is negative.

Instructions for use

1.Indentify the Terragene® Bionova® BT98 SCBI by writing the sterilizer number (in case of having more than one), load number and processing date on the label.

2.Pack the SCBI along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place the package in those areas which are considered most inaccessible for the sterilizing agent (e.g., the center of the load and areas near the door).

3.Sterilize as usual.

4.After the sterilization process has finished, open the sterilizer door, wait five minutes and remove the SCBI from the package.
CAUTION: Wear safety glasses and gloves when removing the Terragene® Bionova® BT98 SCBI from the sterilized package.
WARNING: Do not crush or handle the SCBI excessively, since this might cause the glass ampoule to burst. Allow the SCBI to cool down until it reaches room temperature.

5.Check the process indicator on the SCBI label. A color change to green indicates that the SCBI has been exposed to hydrogen peroxide.
IMPORTANT: This color change does not evidence the process’s effectiveness in achieving sterility. If the color of the process indicator has not changed, check the sterilization process.

6.Press the lid to seal the tube. Crush the ampoule contained in the SCBI with an individual ampoule crusher or with the ampoule crusher placed in the rear of the Bionova® Hyper Auto-reader incubator (BHY). Then shake the tube vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the SCBI in the incubator.
IMPORTANT: Use a non-sterilized SCBI as a positive control at least once per day when performing a sterilization cycle. The positive control ensures that the correct incubation conditions have been followed, that the medium promotes rapid growth, that spore viability has not been compromised by improper storage temperature, humidity or proximity to chemicals, and that the Terragene® Bionova® Auto-reader incubators are functioning properly. Both the positive control indicator and the processed indicator should belong to the same batch.

7.Incubate the processed indicator and the positive control indicator in the Bionova® Hyper Auto-reader incubator for a maximum of 5 minutes at (60 ± 2) °C or hyper rapid readout.
NOTE: The waiting time between sterilization and incubation should not exceed 1 hour. Fluorescence detection by the reader (excitation 340-380 nm/emission 455-465 nm) means failure in the sterilization process. If no fluorescence is detected after 5 minutes of incubation, the result is negative. The positive control must give positive fluorescence readout. It is good practice to incubate a positive control for visual color change. Record positive results and discard the SCBIs immediately, as it is indicated below.

Hyper rapid readout: 5 minutes

The hyper rapid readout must be carried out in the Bionova® Hyper Auto-reader incubator. Fluorescence is emitted after the reader stimulates the spore carrier with UV light. Final readout of negative results is readily available after 5 minutes of incubation. Fluorescence readout is an indirect measure of the viability of *Geobacillus stearothermophilus* spores after the sterilization process (positive results). Furthermore, a failure in the sterilization process can also become evident by culture medium color change.

Visual confirmation: 5 days

Optionally, you can perform a visual color change confirmation after a 5-day incubation. If the sterilization process has not been successful, culture medium will turn to yellow during incubation at (60 ± 2) °C, indicating the presence of living spores. If sterilization was successful, culture medium will remain purple after incubation. The positive control must show a color change from purple to yellow for results to be valid.
NOTE: If a 5-day readout is performed, a humidified environment will be required to prevent the medium from drying out.

Readout time: 7 days

A 7-day readout is optional and may be routinely performed. This is an initial validation of the 5-minute readout, where fluorescence results are compared to the 7-day visual readouts.
NOTE: If a 7-day readout is performed, a humidified environment will be required to prevent the medium from drying out.

Monitoring frequency

Follow the policies and procedures with the monitoring frequency specified by the professional associations and/or standards corresponding to your country. As recommended practice and to ensure optimal patient safety, Terragene S.A. recommends that each sterilization load is monitored with the appropriate biological indicator.

Storage
Store in a dark place at temperature between 10-30 °C and 30-80 % relative humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Shelf-life
Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date. The chemical process indicator on the SCBI label has an expiration date of 2 years when used as part of the SCBI. Endpoint stability reaction: the chemical process indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal
Discard biological indicators after use according to your country’s healthcare and safety regulations. The positive biological indicators can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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