Spore Coupon

BTC90-2/4 , BTC90-2/5 , BTC90-2/6

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Parameters determined at time of manufacture according to ISO 11138-1: 2017 and IRAM 37102-1: 1999 standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant
ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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EN Spore Coupon

For Hydrogen Peroxide sterilization

Composition
Each Bionova® Coupon is soaked with a *Geobacillus stearothermophilus* ATCC 7953 spore population.

Product description
Bionova® Spore Coupons have been specifically designed to be used in Process Challenge Devices to control Plasma or Vaporized Hydrogen Peroxide sterilization processes. Alternatively, spore coupons can also be used for load control in Plasma or Vaporized sterilization processes. Once the sterilization process has finished, spore coupons should be incubated in a Bionova® MC20-2 o MC1020-2 culture medium tube or in an appropriate culture medium for *Geobacillus stearothermophilus* growth. If the sterilization process has not been successful, MC20-2 o MC1020-2 culture medium will turn to yellow after incubation at 60±2 °C, thus indicating the existence of living spores on the coupon. If the sterilization process was efficient, MC20-2 o MC1020-2 growth medium will remain its original color. A final readout should be carried out after a 24-hour incubation at 60±2 °C.

Precautions
Do not use spore coupons to control Steam, Chemical Vapors, Dry Heat, Ethylene Oxide, Radiation cycles or any other sterilization processes. Do not re-use spore coupons.

Storage
Best stored in the original box under the following conditions: temperature between 10-30 °C, 30-80 % relative humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Instructions for use
1-A. Remove the screw cap from the Helix PCD device. Only for metal capsule Helix devices, remove the bolt attached to the cap.
1-B. Place the spore coupon in the lumen of the capsule.
2. Close the lid of the capsule by applying moderate pressure.
NOTE: Terragene® suggests using the system Chemdye® Helix-PCD designed for small sterilizers (CDPCD2X025). The Helix-PCD system should be placed in a completely empty chamber. Alternatively, the spore coupon might be placed in an appropriate package along with the load to be sterilized in accordance with recommended sterilization practices. Place the coupon in those areas which a priori you consider the most inaccessible for the sterilizing agent. A problematic area is the center of the load.
3. For Hydrogen Peroxide processes, start the selected program and sterilize as usual.
4. After the sterilization process has finished, remove the coupons

from the package for processing and incubation.
5. Carefully transfer the spore coupons with a sterile clamp, under sterility conditions, to the MC20-2 o MC1020-2 culture medium tube, or to any other appropriate culture medium. **IMPORTANT:** Use latex gloves and a chin strap when transferring spore coupons to the culture medium tube. **WARNING:** Avoid contact of spore coupons with any kind of surface including the outer part of the culture medium tube.
6. Incubate spore coupons at 60±2 °C.
IMPORTANT: Use an unprocessed spore coupon as a positive control to ensure that incubation conditions were met.
7. Incubate processed spore coupons along with the positive control up to 24 hours at 60±2 °C. It is advisable to make observations every 8 hours. Color change of MC20-2 o MC1020-2 growth medium to yellow shows failure in the sterilization process. If after incubation, there is no color change in the growth medium, the result is negative (the sterilization process was successful). MC20-2 o MC1020-2 culture medium used for incubating the spores coupons which have not been sterilized (positive control) should turn to yellow for the results to be valid. Record the positive results and discard them immediately as indicated below.
ATTENTION: Do not re-use the sterilizer until spore coupons growth results are negative (MC20-2 o MC1020-2 culture medium containing processed spore coupons remains its original color).
WARNING: For incubation in culture media different from those specified by the manufacturer, incubate for a maximum of 7 days making observations every 24 hours.

Disposal
Discard Bionova® culture medium tubes after use according to your country's healthcare and safety regulations. Positive MC20-2 o MC1020-2 growth medium tubes can be sterilized in gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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