## Discos con esporas Spore discs Discos com esporos

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For Ethylene Oxide and Dry Heat sterilization

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
Parameters determined at time of manufacture according to EN ISO 11138-1:2006, EN ISO 11138-2:2009, EN ISO 11138-4:2006 and IRAM 37102:1999 (Parts 1 and 2) standards. The showed values are reproducible only under the same conditions under which they were determined.

Lic. Adrián J. Rovetto  
Director Técnico  
Quality Assurance Director

Industria Argentina - Made in Argentina

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## Spore discs BTD40-8/6 For Ethylene Oxide and Dry Heat sterilization

English

### Composition

Each Bionova® paper disc is soaked with *Bacillus atrophaeus* ATCC 9372 spore population.

### Product description

The Bionova® spore disc is specifically designed to control ethylene oxide (EO) and dry heat (Dry) sterilization processes. Spore disc should be incubated in Bionova® MC1020 or MC1030 growth medium tube or in growth medium appropriate for the corresponding *Bacillus atrophaeus* ATCC 9372 growth.

If sterilization process was not successful, MC1020 or MC1030 growth medium will turn to yellow after incubation at 37±2 °C, thus indicating the existence of living spores on the disc.

If sterilization process was correct, MC1020 or MC1030 indicating medium will keep its original color. A final reading should therefore be done after 48 hours at 37±2 °C.

### Warning!

Do not use spore disc to control steam, plasma, chemical steams, radiation cycles or other sterilization processes.

Do not re-use spore discs.

### Storage

Store in a dark place and at temperatures between 15-30 °C, 35-60% relative humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

### Directions for use

1. Pack the discs in appropriate package along with the material for sterilization, according to recommended sterilization practices. Place the packet in those areas you a priori consider more inaccessible for sterilizing agent. A typical problematic area is the load center.

2. Sterilize as usual.

3. After sterilization process has finished, remove the discs from the package for processing and incubation.

4. Transfer spore discs with a sterile clamp, very carefully under sterility conditions MC1020 or MC1030 growth medium tube, or to any other appropriate growth medium. **IMPORTANT!** Use latex gloves and chinstrap when transferring spore discs to growth medium tube. **WARNING!** Avoid spore discs contact with any kind of surface including the outer part of growth medium tube.

5. Incubate spore discs at 37±2 °C.

**IMPORTANT:** Use a spore disc which has not been under sterilization process as a positive control each time a processed disc is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.

6. Incubate processed spore discs along with the strip used as positive control. It is advisable to make observations every 10 hours. Color change of MC1020 or MC1030 growth medium to yellow shows a failure on sterilization process. If after incubation, no color change is visible on the growth media, the result is negative (i.e. the sterilization process was successful). MC1020 or MC1030 indicating medium used for incubating the spores discs which have not been sterilized (positive control) should turn to yellow in order for the results to be valid. Record the positive ones and discharge them immediately as indicated below.

**WARNING!** Do not re-use sterilizer until spore discs growth result turns negative (MC1020 or MC1030 growth medium containing processed spore discs remains in its original color).

### Disposal

Dispose of the growth medium tubes already used to analyse spore discs growth according to health regulations in your country. Positive MC1020 or MC1030 growth medium tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.

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