Discos con esporas Spore discs

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For Steam, Ethylene Oxide and Dry Heat sterilization

Certificado de calidad Quality certification Bionova® Código/Code: **BTD60-4/--**

Esterilización por Vapor, Óxido de Etileno y Calor Seco
Steam, Ethylene Oxide and Dry Heat sterilization
A) *Bacillus atrophaeus* ATCC 9372
B) *Geobacillus stearothermophilus* ATCC 7953

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Parameters determined at time of manufacture according to UNE-EN ISO 11138-1: 2007, UNE-EN ISO 11138-2: 2009, UNE-EN ISO 11138-3: 2009, UNE-EN ISO 11138-4: 2007 and IRAM 37102: 1999 (Parts 1, 2 and 3) standards. The values shown are reproducible only under the same conditions which they were determined.

ISO and USP Compliant
ATCC is registered trademark of American Type Culture Collection

Lic. Adrián J. Róvotto
Director Técnico
Quality Assurance Director

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Spore discs

For Steam, Ethylene Oxide and Dry Heat sterilization
BTD60-4/5, BTD60-4/6

English

Composition

Each Bionova® paper disc is soaked with *Geobacillus stearothermophilus* ATCC 7953 and *Bacillus atrophaeus* ATCC 9372 spore population.

Product description

The Bionova® spore disc is specifically designed to control steam, ethylene oxide (EO) and dry heat sterilization processes. Spore discs should be incubated in an appropriate growth medium for the corresponding *Geobacillus stearothermophilus* or *Bacillus atrophaeus* growth. It is suggested the utilization of the culture medium Bionova® MC20 for the growth of *Geobacillus stearothermophilus* ATCC 7953 at 60 °C, MC1030 for the growth of *Bacillus atrophaeus* ATCC 9372 at 37 °C and/or MC1020 for growth of both spores, selective according to the incubation temperature (*Geobacillus stearothermophilus* 60 °C/*Bacillus atrophaeus* 37 °C).

If the sterilization process was not successful, MC20, MC1030 or MC1020 growth medium will turn to yellow after incubation at the corresponding temperature, thus indicating the existence of living spores on the discs.

If the sterilization process was correct, MC20, MC1030 or MC1020 indicating medium will keep its original color. A final reading should therefore be done after 24 hours of growth medium incubation at 60±2 °C (MC20 and/or MC1020) for steam sterilization control processes, or 48 hours at 37±2 °C (MC1030 and/or MC1020) for ethylene oxide and dry heat sterilization control processes.

Warning!

Do not use spore discs to control formaldehyde or any sterilization processes other than steam, EO and/or dry heat.
Do not re-use the spore discs.

Storage

Store in a dark place and at temperatures between 10-30 °C, 30-80% relative humidity.
Do not freeze.
Do not store near sterilizing agents or other chemical products.

Directions for use

1. Pack the Bionova® discs in appropriate package along with items for sterilization, according to recommended sterilization practices. Place the disc in those areas you a priori consider more inaccessible for sterilizing agent. A typical problematic area is the center of the load.

2. Sterilize as usual.

3. After sterilization process has finished, remove the Bionova® discs from the package for processing and incubation.

4. Transfer spore discs with a sterile clamp, very carefully under sterility conditions to MC20, MC1030 or MC1020 growth medium tube, or to any other appropriate growth medium. **IMPORTANT:** Use latex gloves and chinstrap when transferring spore discs to growth medium tube. **WARNING!** Avoid spore disc contact with any kind of surface including the outer part of growth medium tube.

5. Incubate spore discs at the corresponding temperatures: 60±2 °C or 37±2 °C.

IMPORTANT: Use a spore disc which has not been under sterilization process as a positive control each time a processed disc is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.


6. Incubate processed spore discs along with the disc used as positive control. It is advisable to make observations every 10 hours.

Color change of MC20, MC1030 or MC1020 growth medium to yellow shows a failure on sterilization process. If after incubation, no color change is visible on the growth media, the result is negative (i.e. the sterilization process was successful). The positive control indicator medium MC20, MC1030 or MC1020 should show a yellow color change for the results to be valid. Record the positive ones and discharge them immediately as indicated below.

WARNING! Do not re-use sterilizer until spore discs growth result turns negative (MC20, MC1030 or MC1020 growth medium containing processed spore discs remains in its original color).

Disposal

Dispose of the growth medium tubes already used to analyse spore discs growth according to health regulations in your country. Positive MC20, MC1030 or MC1020 growth medium tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.

Industria Argentina - Made in Argentina

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Fabricado por: Terragene S.A. - Güemes 2879 - (2000) Rosario - Santa Fe - Argentina