**Certificado de calidad**  
**Quality certification**  
**Bionova® BTD92-8/6**

VH<sub>2</sub>O<sub>2</sub>

Esterilización con Plasma o Vapor de Peróxido de Hidrógeno.  
Plasma or Vaporized Hydrogen Peroxide sterilization  
*Geobacillus stearothermophilus* ATCC 7953

LOT

Población / \_\_\_\_\_ UFC /  
Population \_\_\_\_\_ CFU

Valor D/ \_\_\_\_\_ segundos/  
D-value \_\_\_\_\_ seconds  
(2 mg/L VH<sub>2</sub>O<sub>2</sub> Vapor, 50°C)

Tiempo sobrevida / \_\_\_\_\_ segundos/  
Survival time \_\_\_\_\_ seconds  
Survival time = (log<sub>10</sub> labeled population - 2) x labeled D-value

Tiempo de muerte / \_\_\_\_\_ segundos/  
Kill time \_\_\_\_\_ seconds  
Kill time = (log<sub>10</sub> labeled population + 4) x labeled D-value

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**Composition**

Each Bionova® polypropylene disc is soaked with *Geobacillus stearothermophilus* ATCC 7953 spore population.

**Product description**

The Bionova® spore disc is specifically designed to control Plasma or Vaporized Hydrogen Peroxide sterilization processes. Spore discs should be incubated in Bionova® MC20 or MC1020 Growth Medium tube or in appropriate growth medium for *Geobacillus stearothermophilus* growth.

If the sterilization process was not successful, MC20 or MC1020 Growth Medium will turn to yellow after incubation at 60±2 °C, thus indicating the existence of living *Geobacillus stearothermophilus* spores in the discs.

If the sterilization process was successful, MC20 or MC1020 Growth Indicating Medium will remain its original color. Final readout should therefore be performed after 24 hours of growth medium incubation at 60±2 °C.

**Warning!**

Do not use spore discs to control EO, Dry Heat, Radiation, Steam or any other sterilization processes.  
Do not reuse the spore discs.

**Storage**

Best stored in original box under the next conditions: Temperatures between 10-30 °C, 30-80 % Relative Humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

**Directions for use**

1. Pack the discs in an appropriate package along with material to be sterilized according to recommended sterilization practices. Place it in those areas you consider a priori more inaccessible for sterilizing agent. The load center is generally a problematic area.

2. Sterilize as usual.

3. After sterilization process has finished, remove Bionova® discs from the package used for processing and incubation.

4. In sterility conditions transfer the spore discs using a sterile clamp, to MC20 or MC1020 growth medium tube, or to any other appropriate growth medium. **IMPORTANT:** Use latex gloves and chin strap when transferring spore discs to growth medium tube. **WARNING!** Avoid spore disc contact with any kind of surface including the outer part of growth medium tube.

**IMPORTANT:** Use a spore disc which has not been under sterilization process as a positive control each time a processed disc is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.

5. Incubate processed spore discs along with the disc used as positive control during a maximum of 24 hours at 60±2 °C. It is advisable to make observations every 10 hours. Color change of MC20 or MC1020 Growth Medium to yellow means a failure on sterilization process. If after 24 hours no color change is visible on MC20 or MC1020 Growth Medium, the result is negative (i.e. the sterilization process was successful). MC20 and MC1020 Indicating Medium used for incubating the spores discs which have not been sterilized (positive control) should turn to yellow in order for the results to be valid. Record the positive samples and dismiss them immediately as indicated below.

**WARNING!** Do not reuse sterilizer until spore discs growth result turns negative (MC20 or MC1020 growth medium containing processed spore disc remains its original color).

**Disposal**

Discard the growth medium tubes used to evaluate the spore growth of Bionova® Spore Discs according to your country's healthcare and safety regulations. Positive MC20 or MC1020 Growth Medium tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted Steam sterilizers.

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ATCC is a registered trademark of American Type Culture Collection