## Mini spore strips

BTM40.2-10 | BTM40G

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
Z-value / Valor Z \_\_\_\_\_ °C

Parameters determined at time of manufacture according to ISO 11138-1: 2017, ISO 11138-2: 2017, ISO 11138-4: 2017 and IRAM 37102: 1999 (Parts 1 and 2) standards. The values shown are reproducible only under the same conditions which they were determined.

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ISO and USP Compliant.  
ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rowetto  
Director Técnico  
Technical Director

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## EN Mini Spore strips

For Ethylene Oxide and Dry Heat sterilization

## Composition

Each Bionova® Paper Strip is soaked with *Bacillus atrophaeus* ATCC 9372 spore population.

## Product description

Bionova® Spore Strip is designed to control Ethylene Oxide (EO) and Dry Heat (Dry) sterilization processes. Spore strips should be incubated in Bionova® MC1030-2 or MC1020-2 Growth Medium tube or other appropriate growth medium for *Bacillus atrophaeus* growth.If sterilization process was not successful, MC1030-2 or MC1020-2 growth medium will turn to yellow after incubation at 37±2 °C, thus indicating the existence of living *Bacillus atrophaeus* spores on the strip.

If sterilization process was effective, MC1030-2 or MC1020-2 indicator medium will remain its original color. A final readout should therefore be performed after 48 hours of growth medium incubation at 37±2 °C.

## Precautions

Do not use spore strips to control Steam, Radiation cycles or other sterilization processes.

Do not re-use spore strips.

Do not re-use sterilizer until spore strip growth result is negative (MC1030-2 or MC1020-2 growth medium containing processed spore strip remains in its original color).

## Storage

Best stored in the original packaging under the following

conditions: temperatures between 10-30 °C, 30-80 % relative humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

## Directions for use

1. Pack the strip in appropriate package along with items for sterilization, according to recommended sterilization practices. Place the packet in those areas you consider a priori more inaccessible to the sterilizing agent. A typical problematic area is the center of the load.

2. Sterilize as usual.

3. After sterilization process has finished, remove the strip from the package for processing and incubation.

4. Transfer spore strip using a sterile clamp, very carefully under sterility conditions to MC1030-2 or MC1020-2 growth medium tube or to any other appropriate growth medium.

**IMPORTANT:** Use latex gloves and chinstrap when transferring spore strips to growth medium tube. **WARNING:** Avoid spore strip contact with any kind of surface including the outer part of growth medium tube.**IMPORTANT:** Use a spore strip which has not been subjected to the sterilization process as a positive control each time a processed strip is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.

5. Incubate processed spore strips along with the strip used as positive control during a maximum of 48 hours at 37±2 °C. It is advisable to make observations every 10 hours.

Color change of MC1030-2 or MC1020-2 growth medium to yellow means a failure on sterilization process. If after 48-hour incubation, no color change is observed on the growth medium, the result is negative (i.e. the sterilization process was successful). MC1030-2 or MC1020-2 indicator medium used for incubating the spores strip which has not been sterilized (positive control) should turn to yellow in order for the results to be valid. Record the positive ones and discard them immediately as indicated below.

## Shelf life

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date.

## Disposal

Discard the growth medium tubes used to evaluate the spores growth of spore strips according to your country's healthcare and safety regulations. Positive MC1030-2 or MC1020-2 growth medium tubes can be sterilized in gravity air displacement Steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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