## Mini spore strips

BTM50.2-10 | BTM50G

**Quality certification**  
**Certificado de calidad**  
**Bionova® BTM50.2-10 | BTM50G**

 Steam sterilization / Esterilización por Vapor  
*Geobacillus stearothermophilus* ATCC 7953

STEAM

LOT

Population / Población \_\_\_\_\_ CFU / UFC

D-Value / Valor D (121 °C) \_\_\_\_\_ min.

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Survival time = (log<sub>10</sub> labeled population - 2) x labeled D-value

Kill time / Tiempo de muerte \_\_\_\_\_ min.

Kill time = (log<sub>10</sub> labeled population + 4) x labeled D-value

Z-value / Valor Z \_\_\_\_\_ °C

Parameters determined at time of manufacture according to ISO 11138-1: 2017, ISO 11138-3: 2017 e IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions which they were determined.

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
 ISO and USP Compliant  
 ATCC is a registered trademark of American Type Culture Collection

 Lic. Adrián J. Rovetto  
 Director Técnico  
 Technical Director

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For Steam sterilization

**Composition**

 Each Bionova® paper Strip is soaked with *Geobacillus stearothermophilus* ATCC 7953 spore population.

**Product description**

 The Bionova® Spore Strip is specifically designed to control Steam sterilization processes at 121-135 °C. Spore Strips should be incubated in Bionova® MC20-2 or MC1020-2 growth medium tube or other appropriate growth medium for *Geobacillus stearothermophilus* growth.

 If the sterilization process was not successful, MC20-2 or MC1020-2 growth medium will turn to yellow after incubation at 60±2 °C, thus indicating the existence of living *Geobacillus stearothermophilus* spores on the strips.

If sterilization process was successful, MC20-2 or MC1020-2 indicator medium will remain its original color. A final readout should be therefore performed after 24 hours of growth medium incubation at 60±2 °C.

**Precautions**

Do not use spore mini strips to control Radiation, EO, Dry Heat sterilization cycles or other sterilization processes.

Do not re-use the spore strips.

Do not reuse sterilizer until spore strip growth result is negative (MC20-2, MC1020-2 growth medium containing processed spore strips remains its original color).

**Storage**

Best stored in the original box under the following conditions: temperatures between 10-30 °C, 30-80 % relative humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

**Instructions for use**

1. Pack the spore strip in an appropriate package along with material to be sterilized, according to recommended sterilization practices. Place the strip in those areas you consider a priori more inaccessible to the sterilizing agent. A typical problematic area is the center of the load.

2. Sterilize as usual.

3. After sterilization process has finished, remove Bionova® strip from the package used for processing and incubation.

 4. Transfer the spore strip using a sterile clamp, in sterility conditions to MC20-2, MC1020-2 growth medium tube or to any other appropriate growth medium. **IMPORTANT:** Use latex gloves and chin strap when transferring spore strip to growth medium tube. **WARNING:** Avoid spore strip contact with any kind of surface including the outer part of growth medium tube. **IMPORTANT:** Use a spore strip which has not been subjected to the sterilization process as a positive control each time a processed strip is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.

5. Incubate processed spore strip along with the strip used as positive control during a maximum of 24 hours at 60±2 °C. It is advisable to make observations every 10 hours. Color change of MC20-2 or MC1020-2 growth medium to yellow means a failure on the sterilization process. If after incubation, no color change is visible on MC20-2, MC1020-2 growth medium, the result is negative (i.e. the sterilization process was successful). MC20-2, MC1020-2 indicator medium used for incubating the spores strip which have not been sterilized (positive control) should turn to yellow in order for the results to be valid. Record the positive ones and discard them immediately as indicated below.

**Shelf life**

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date.

**Disposal**

Discard the growth medium tubes used to evaluate the spore growth of Bionova® spore strips according to your country's healthcare and safety regulations. Positive MC20-2 or MC1020-2 growth medium tubes can be sterilized in gravity air displacement Steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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