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For Steam, Ethylene Oxide and Dry Heat sterilization

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
Parameters determined at time of manufacture according to ISO 11138 (Parts 1, 2, 3 and 4) and IRAM 37102 (Parts 1, 2 and 3) standards. The showed values are reproducible only under the same conditions under which they were determined.

Lic. Adrián J. Rovetto  
Director Técnico  
Quality Assurance Director

TERRAGENE

Industria Argentina - Made in Argentina

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# Mini Spore Strips For Steam, Ethylene Oxide and Dry Heat sterilization BTM60.2-10/5, BTM60.2-10/6

English

### Composition

Each Bionova® paper strip is soaked with *Geobacillus stearothermophilus* ATCC 7953 and *Bacillus atrophaeus* ATCC 9372 spore population.

### Product description

The Bionova® spore strip is specifically designed to control steam (Steam), ethylene oxide (EO) and dry heat (Dry) sterilization processes. Spore strips should be incubated in Bionova® MC1020, MC20 or MC1030 growth medium tube or in growth medium appropriate for the corresponding *Geobacillus stearothermophilus* or *Bacillus atrophaeus* growth.

If sterilization process was not successful, MC1020, MC20 or MC1030 growth medium will turn to yellow after incubation at 56±2 °C or 37±2 °C respectively, thus indicating the existence of living spores on the strip.

If sterilization process was correct, MC1020, MC20 or MC1030 indicating medium will keep its original color. A final reading should therefore be done after 24 hours of growth medium incubation at 56±2 °C (MC1020, MC20) for steam sterilization control processes, or 48 hours at 37±2 °C (MC1020, MC1030) for ethylene oxide and dry heat sterilization control processes.

### Warning!

Do not use spore strips to control radiation cycles or other sterilization processes.  
Do not re-use spore strips.

### Storage

Store in a dark place and at temperatures between 15 - 30 °C, 35 - 60 % relative humidity.  
Do not freeze.  
Do not store near sterilizing agents or other chemical products.

### Directions for use

- Pack the strips in appropriate package along with items for sterilization, according to recommended sterilization practices. Place the packet in those areas you a priori consider more inaccessible for sterilizing agent. A typical problematic area is the load center.
- Sterilize as usual.
- After sterilization process has finished, remove the strips from the package for processing and incubation.
- Transfer spore strips with a sterile clamp, very carefully under sterility conditions to MC1020, MC20 or MC1030 growth medium tube, or to any other appropriate growth medium. **IMPORTANT!** Use latex gloves and chin strap when transferring spore strips to growth medium tube. **WARNING!** Avoid spore strips contact with any kind of surface including the outer part of growth medium tube.
- Incubate spore strips at the corresponding temperatures: 56±2 °C or 37±2 °C. **IMPORTANT:** Use a spore strip which has not been under sterilization process as a positive control each time a processed strip is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.
- Incubate processed spore strips along with the strip used as positive control. It is advisable to make observations every 10 hours. Color change of MC1020, MC20 or MC1030 growth medium to yellow shows a failure on sterilization process. If after incubation, no color change is visible on the growth media, the result is negative (i.e. the sterilization process was successful). MC1020, MC20 or MC1030 indicating medium used for incubating the spores strips which have not been sterilized (positive control) should turn to yellow in order for the results to be valid. Record the positive ones and discharge them immediately as indicated below. **WARNING!** Do not re-use sterilizer until spore strips growth result turns negative (MC1020, MC20 or MC1030 growth medium containing processed spore strips remains in its original color).

### Disposal

Dispose of the growth medium tubes already used to analyse spore strips growth according to health regulations in your country. Positive MC1020, MC20 or MC1030 growth medium tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for at least 20 minutes or at 132 °C for at least 15 minutes, or at 134 °C for at least 10 minutes in vacuum-assisted steam sterilizers.

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