### Quality certification Certificado de calidad Bionova® BTM92.2-10 | BTM92G

VH202

Plasma or Vaporized Hydrogen Peroxide sterilization/  
Esterilización con Plasma o Vapor de Peróxido de Hidrógeno  
*Geobacillus stearothermophilus* ATCC 7953

LOT

Population / Población \_\_\_\_\_ CFU/ UFC

D - value / Valor D \_\_\_\_\_ sec/seg  
(2 mg/L VHzO<sub>2</sub>, 50 °C)Survival time / Tiempo sobrevivida \_\_\_\_\_ sec/seg  
Survival time = (log10 labeled population - 2) x labeled D-valueKill time / Tiempo de muerte \_\_\_\_\_ sec/seg  
Kill time = (log10 labeled population + 4) x labeled D-value

Parameters determined at time of manufacture according to ISO 11138-1: 2017 e IRAM 37102-1: 1999 standards. The values shown are reproducible only under the same conditions which they were determined.

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
ISO and USP Compliant  
ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rivetto  
Director Técnico  
Technical Director

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### EN Mini Spore strips

For Plasma or Vaporized Hydrogen Peroxide sterilization

#### Composition

Each Bionova® mini strip is soaked with *Geobacillus stearothermophilus* ATCC 7953 spore population.

#### Product description

The Bionova® Spore Strip is specifically designed to control Plasma or Vaporized Hydrogen Peroxide sterilization processes. Spore Strips should be incubated in Bionova® MC20-2 or MC1020-2 growth medium tube or other appropriate growth medium for *Geobacillus stearothermophilus* growth. If the sterilization process was not successful, MC20-2 or MC1020-2 growth medium will turn to yellow after incubation at 60±2 °C, thus indicating the existence of living *Geobacillus stearothermophilus* spores on the strips.

If sterilization process was successful, MC20-2 or MC1020-2 indicator medium will remain its original color. A final readout should be therefore performed after 24 hours of growth medium incubation at 60±2 °C.

#### Precautions

Do not use spore strips to control Steam, EO, Dry Heat, Radiation or any sterilization processes other than Plasma.

Do not reuse the spore strips.

Do not reuse sterilizer until spore strip growth result is negative (MC20-2, MC1020-2 growth medium containing processed spore strips remains its original color).

#### Storage

Best stored in the original box under the following conditions: temperatures between 10-30 °C, 30-80 % relative humidity.

Do not freeze.

Do not store near sterilizing agents or other chemical products.

#### Directions for use

1. Pack the strips in an appropriate package along with material to be sterilized, according to recommended sterilization practices. Place the strips in those areas you consider a priori more inaccessible to the sterilizing agent. A typical problematic area is the center of the load.
2. Sterilize as usual.
3. After sterilization process has finished, remove Bionova® strips from the package for processing and incubation.

4. Transfer the spore strips using a sterile clamp, in sterility conditions to MC20-2, MC1020-2 growth medium tube or to any other appropriate growth medium. **IMPORTANT:** Use latex gloves and chin strap when transferring spore strips to growth medium tube. **WARNING:** Avoid spore strips contact with any kind of surface including the outer part of growth medium tube. **IMPORTANT:** Use a spore strip which has not been subjected to the sterilization process as a positive control each time a processed strip is incubated. Both, processed spore strip and the strip used as positive control must belong to the same batch. Positive control guarantees that the incubation was carried out under appropriate conditions.

5. Incubate processed spore strip along with the strip used as positive control during a maximum of 24 hours at 60±2 °C. It is advisable to make observations every 10 hours. Color change of MC20-2 or MC1020-2 growth medium to yellow means a failure on the sterilization process. If after incubation, no color change is visible on MC20-2 or MC1020-2 growth medium, the result is negative (i.e. the sterilization process was successful). MC20-2 or MC1020-2 indicator medium used for incubating the spores strips which have not been sterilized (positive control) should turn to yellow in order for the results to be valid. Record the positive ones and discard them immediately as indicated below.

#### Shelf life

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date.

#### Disposal

Discard the growth medium tubes used to evaluate the spore growth of Bionova® BTM92 spore strips according to your country's healthcare and safety regulations. Positive MC20-2 or MC1020-2 growth medium tubes can be sterilized in gravity air displacement Steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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