## Culture Medium

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ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto  
Director Técnico  
Technical Director

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#### Explanation of Symbols

- Batch number.
- Manufacture Date.
- Expiration Date.
- Manufacturer.

## EN Culture Medium

### Composition

Each container has 2 ml of growth indicator medium for *Bacillus atrophaeus* ATCC 9372, *Bacillus subtilis* ATCC 35021, *Bacillus pumilus* ATCC 27142 and *Geobacillus stearothermophilus* ATCC 7953.

### Product description

Bionova® MC1020-2 Culture Medium is designed for *Bacillus atrophaeus* ATCC 9372, *Bacillus subtilis* ATCC 35021, *Bacillus pumilus* ATCC 27142 and *Geobacillus stearothermophilus* ATCC 7953 spores growth.

If the sterilization process was not effective, the culture medium will turn to yellow during incubation, indicating the presence of living spores. The incubation conditions will depend on the microorganism and the sterilizing agent:

24 hours between 55-62 °C for spores of *Geobacillus stearothermophilus* ATCC 7953 subjected to Steam or VH<sub>2</sub>O<sub>2</sub> processes and 48 hours for Formaldehyde processes; 48 hours between 35-39 °C for *Bacillus atrophaeus* ATCC 9372 spores subjected to Dry Heat or Ethylene Oxide processes; for *Bacillus subtilis* ATCC 35021 spores subjected to low temperature steam processes; and for *Bacillus pumilus* ATCC 27142 spores subjected to Gamma Radiation processes.

If the sterilization process was effective, the culture medium will remain its original color. Final readout should be carried out after 48 or 24 hours of incubation between 35-39 °C or 55-62 °C, as appropriate.

### Precautions

Bionova® MC1020-2 Culture Medium should not be used for other bacterial species growth different from *Bacillus atrophaeus* ATCC 9372, *Bacillus subtilis* ATCC 35021, *Bacillus pumilus* ATCC 27142 or *Geobacillus stearothermophilus* ATCC 7953.

Do not reuse.

Do not use the sterilizer until the spore carrier growth result is negative (MC1020-2 Culture Medium containing a processed spore carrier remains the original color after the 48 or 24 hours of incubation).

### Storage

Store it in a dark place at temperatures between 4-20 °C, 30-80 % relative humidity.

Do not freeze.

### Instructions for use

1. Identify Bionova® MC1020-2 tube by writing on the label the

sterilizer number (in case of having more than one), load number and processing date.

2. Open the envelope containing the spore carrier by cutting its corner off, under sterile conditions (for example: laminar flow) and transfer using a sterile clamp to MC1020-2 growth medium tube.

**WARNING:** Avoid contact with the carrier and any surface, including the outside of the culture medium tube.

**IMPORTANT:** Use latex gloves and a mask when transferring the spore carrier from the envelope to the MC1020-2 Culture Medium tube.

Every time a processed carrier is incubated, a non-sterilized spore carrier should be used as a positive control. Both, the processed carrier and the carrier used as positive control must belong to the same batch. A positive control assures proper incubation conditions.

3. Incubate the spore carrier inside the MC1020-2 Culture Medium at the corresponding temperatures between 35-39 °C or 55-62 °C. A color change to yellow of the growth indicator medium shows sterilization process has failed. A final negative result is obtained if there is no color change of the MC1020-2 Culture Medium after 48 or 24 hours, meaning sterilization process has been effective. The color of the MC1020-2 Indicator Medium used to incubate the non-sterilized spore carrier (positive control) should turn to yellow for the results to be valid.

Positive results should be recorded and discarded immediately, as explained below.

### Disposal

Discard the culture medium tubes used for analyzing the spores growth according to your country's healthcare and safety regulations. MC1020-2 Culture Medium positive tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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