## Culture Medium

### Quality certification Certificado de calidad Bionova® MC1030-2

For the growth of *Bacillus atrophaeus* ATCC 9372 and *Bacillus subtilis* ATCC 35021  
Para el crecimiento de *Bacillus atrophaeus* ATCC 9372 y *Bacillus subtilis* ATCC 35021

LOT

#### Color transition / Transición de color

Turns to yellow / Vira al amarillo

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ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rowetto  
Director Técnico  
Technical Director

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#### Explanation of Symbols

- Batch number.
- Manufacture Date.
- Expiration Date.
- Manufacturer.

## EN Culture Medium

### Composition

Each container has 2 ml of growth indicator medium for *Bacillus atrophaeus* ATCC 9372 and *Bacillus subtilis* ATCC 35021.

### Product description

Bionova® MC1030-2 Culture Medium is designed for *Bacillus atrophaeus* ATCC 9372 and *Bacillus subtilis* ATCC 35021 growth. If the sterilization process was not successful, indicator medium will change to yellow after incubation at (37 ± 2) °C for 48 hours for *Bacillus atrophaeus* ATCC 9372 spores subjected to Dry Heat and Ethylene Oxide processes, or *Bacillus subtilis* ATCC 35021 spores subjected to Steam at low temperature, thus showing the presence of living *Bacillus atrophaeus* ATCC 9372 or *Bacillus subtilis* ATCC 35021 spores, respectively. If the sterilization process was effective the indicator medium will remain its original color. Final readout should be carried out after 48 hours of incubation at (37 ± 2) °C.

### Precautions

Bionova® MC1030-2 Culture Medium should not be used for growth of other bacterial species different from *Bacillus atrophaeus* ATCC 9372 or *Bacillus subtilis* ATCC 35021. Do not reuse.

### Storage

Store it in a dark place at temperatures between 4-20 °C, 30-80 % relative humidity. Do not freeze.

### Instructions for use

- Identify Bionova® MC1030-2 Indicator tube by writing the sterilizer number (if there is more than one), load number and processing date.
- Open the envelope containing the spores carrier on the corner, under sterile conditions (for example: laminar flow) and transfer using a sterile forcep.  
**IMPORTANT:** Use latex gloves and a mask when transferring the spores carrier from the envelope to the MC1030-2 Culture Medium tube.  
**WARNING:** Avoid contact with the carrier and any surface, including the outside of the culture medium tube.  
**IMPORTANT:** Every time a processed carrier is incubated, a non-sterilized spore carrier should be used as a positive control. A positive control assures proper incubation conditions.
- Incubate the spore carrier inside the MC1030-2 Culture

Medium at (37 ± 2) °C. A color change to yellow of the growth indicator medium shows sterilization process has failed. If after 48 hours of incubation there is no color change of the MC1030-2 Culture Medium, the result is negative (the sterilization process has been successful). The color of the MC1030-2 Indicator Medium used to incubate the non-sterilized spore carrier (positive control) should turn to yellow for the results to be valid.

Positive results should be recorded and discarded immediately, as explained below.

**WARNING:** Do not use the sterilizer until the growth result from the spore carrier is negative (MC1030-2 Culture Medium containing a processed spores carrier remains its original color after 48 hours of incubation).

### Disposal

Discard the culture medium tubes used for analyzing the spores growth according to your country's healthcare and safety regulations. MC1030-2 Culture Medium positive tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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