Culture Medium

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For the growth of *Geobacillus stearothermophilus*
ATCC 7953

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24 hs of incubation between 55-62 °C (spores subject to Steam/VH2O2)

48 hs of incubation between 55-62 °C (spores subject to Formaldehyde)

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ISO and USP Compliant

ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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Explanation of Symbols

- Batch number.
- Manufacture Date.
- Expiration Date.
- Manufacturer.

EN Culture Medium

Composition

Each container has 2 ml of Growth Indicator medium for
Geobacillus stearothermophilus ATCC 7953.

Product description

Bionova® MC20-2 Culture Medium is designed for *Geobacillus*
stearothermophilus ATCC 7953 growth.

If the sterilization process was not effective, the indicator
medium will turn to yellow after incubation between 55-62 °C
for 24 hours for spores subject to Steam/VH2O2 sterilization
processes and 48 hours for Formaldehyde processes
(maximum), thus showing the presence of living *Geobacillus*
stearothermophilus ATCC 7953 spores.

If the sterilization process was effective the indicator medium
will remain its original color. Final readout should be carried
out at 24 or 48 hours of incubation between 55-62 °C,
depending on the sterilization process.

Precautions

Bionova® MC20-2 Culture Medium should not be used for
growth of *Bacillus atrophaeus* or other bacterial species
different from *Geobacillus stearothermophilus* ATCC 7953.
Do not reuse.

Do not reuse the sterilizer until the spore carrier growth result
is negative (MC20-2 Culture Medium containing a processed
spores carrier remains the original color).

Storage

Store it in a dark place at temperatures between 4-20 °C, 30-80 %
relative humidity.
Do not freeze.

Instructions for use

1. Identify Bionova® MC20-2 tube by writing on the label the
sterilizer number (in case of having more than one), load
number and processing date.
2. Open the envelope containing the spore carrier under sterile
conditions (for example: laminar flow) and transfer using a
sterile clamp to MC20-2 growth medium tube.

WARNING: Avoid contact between carrier and any surface,
including the outside of the culture medium tube.

IMPORTANT: Use latex gloves and a mask when transferring
the spore carrier from the envelope to the MC20-2 Culture
Medium tube. Every time a processed carrier is incubated, a
non-sterilized spore carrier should be used as a positive
control. Both, the processed carrier and the carrier used as
positive control must belong to the same batch. A positive

control assures proper incubation conditions.

3. Incubate the spore carrier inside the MC20-2 Culture
Medium between 55-62 °C. Color change of the growth
indicator medium to yellow shows the sterilization process
has failed. A final negative result is obtained if there is no color
change of the MC20-2 Culture Medium at 24 hours of
incubation (spores subject to Steam/VH2O2 processes) or 48
hours of incubation (spores subject to Formaldehyde
processes), meaning the sterilization process has been
effective. The color of the MC20-2 Indicator Medium used to
incubate the unprocessed spore carrier (positive control)
should turn to yellow for the result to be valid.

Positive results should be recorded and discarded
immediately, as explained below.

Disposal

Discard the culture medium tubes used for analyzing the spores
growth, according to your country's healthcare and safety
regulations. MC20-2 Culture Medium positive tubes can be
sterilized in gravity air displacement steam sterilizer at 121 °C
for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes;
or in a dynamic air removal steam sterilizer at 132 °C for 4
minutes or 135 °C for 3 minutes.

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