Process Challenge Device (PCD)

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Parameters determined at time of manufacture according to ISO 11138-1:2017, ISO 11138-2:2017 and IRAM 37102: 1999 (Parts 1 and 2) standards. The shown values are reproducible only under the same conditions under which they were determined.

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
ISO and USP Compliant.
ATCC is a registered trademark of American Type Culture Collection.

Terragene® and Bionova® are registered trademarks of Terragene S.A.

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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EN Process Challenge Device

For a quick and easy monitoring of Ethylene Oxide sterilization processes

Product Description

The Bionova® PCD110 Process Challenge Device (PCD) has been designed for quick and easy monitoring of Ethylene Oxide sterilization processes.

A PCD is a device that simulates a load to be sterilized and constitutes a challenge to the sterilization process. It is used to routinely evaluate the process performance by detecting inadequate Ethylene Oxide penetration. It also allows for periodic monitoring and validation of the sterilizer (after a repair, installation, or relocation).

Bionova® PCD110 consists of a pre-assembled disposable device designed to be as resistant as the test device described in ANSI/AAMI Standards. It contains within it a BT110 Self-Contained Biological Indicator (SCBI) with a Chemical Process Indicator printed on its label. The latter is a Type 1 indicator, according to ISO 11140-1: 2014, which turns from brown to green when exposed to Ethylene Oxide. On its exterior, the Bionova® PCD110 process challenge device also contains a Process Indicator that turns orange after exposure.

Each PCD is designed to challenge the route and arrival of the sterilizing agent to the population of *Bacillus atrophaeus* ATCC 9372 spores contained within the SCBI BT110.

Indications for use

WARNING: Do not use the Bionova® PCD110 to control Steam, Dry Heat, Formaldehyde, or other sterilization processes other than those indicated. Do not open the device before sterilization. Do not reuse the device.

Instructions for use

- Place the device in the area that you consider a priori most inaccessible to the sterilizing agent. Generally, a problematic area is the center of the load or near the sterilizer door.
- Sterilize as usual.
- After the sterilization process is finished, open the sterilizer door, wait 5 minutes, and remove the process challenge device. **CAUTION:** Wear

gloves and safety goggles when removing the Bionova® PCD110.

4. Verify that the Process Indicator on the device changed to orange.

5. Remove the Biological Indicator from the test device. Dispose of the remaining components of the challenge device according to your institution's policies.

6. Verify that the Process Indicator printed on the SCBI label changed to green. The color change confirms that the indicator was exposed to Ethylene Oxide. **IMPORTANT:** This color change does not evidence that the process was sufficient to achieve sterility. If the Process Indicator did not change color, it is necessary to check the sterilization process.

7. Press the cap of the indicator to seal the tube. Break the ampoule contained in the SCBI with an individual ampoule crusher or with the ampoule crusher provided in the Terragene® Bionova® incubator. Then, shake the indicator vigorously with movements similar to those used to lower the temperature in a mercury thermometer, until the culture medium drops and completely soaks the spore carrier. **IMPORTANT:** Use a non-sterilized SCBI as a positive control at least once per day, when using the same indicator batch. The positive control ensures that incubation conditions were adequate and that spore viability was not altered due to inadequate storage temperature, humidity, or proximity to chemicals. Also, the ability of the medium to promote rapid growth and proper functioning of the Terragene® Bionova® Auto-Reader Incubator is evaluated. The processed indicator and the positive control must belong to the same batch.

8. Incubate both the processed indicator and the indicator used as a positive control in the appropriate Terragene® Bionova® auto-reader incubator for up to 4 hours at (37 ± 2) °C for a rapid readout. **NOTE:** The waiting time between sterilization and incubation should not exceed 7 days. The detection of fluorescence by the auto-reader indicates a failure in the sterilization process. If after 4 hours no fluorescence is detected, the result is negative, and, therefore, the sterilization process was effective. For the result to be valid, the indicator used as a positive control must be detected as such in the auto-reader.

If the relative humidity of the sterilizer chamber is unknown or lower than 35 %, continue incubating the biological indicator for a pH color change of the medium, since low relative humidity may increase the fluorescence reading time beyond 4 hours.

Record the results and discard the biological indicators immediately, as indicated below.

Rapid readout: 4 hours

The rapid readout should be carried out in Terragene® Bionova® Auto-Reader incubators. Fluorescence is emitted after stimulating the spore carrier with UV light. The final readout of the negative results is available after 4 hours of incubation. The fluorescence readout is an indirect determination of the viability of *Bacillus atrophaeus* spores after the sterilization process (positive result). On the other hand, a failure in the sterilization process can also be evidenced by the color change of the culture medium after 48 hours of incubation. Due to the high sensitivity of the results obtained by fluorescence at 4 hours, conventional incubation of the SCBI Terragene® Bionova® BT110 for evidence of color change of the medium does not represent an additional advantage.

Visual confirmation: 48 hours

Optionally, it is possible to perform a color-change visual confirmation after 48 hours of incubation. If the sterilization process was unsuccessful, the culture medium will change to yellow during incubation at (37 ± 2) °C, indicating the presence of live spores. If the sterilization was successful, the culture medium will remain blue. A definitive negative result is obtained after 48 hours of incubation. The positive control must show a color change from blue to yellow for the results to be valid.

7-day reading

An optional 7-day visual confirmation can be performed periodically. It is an initial validation of the 4-hour rapid readout. The results obtained by fluorescence measurements are compared with the visual confirmation (color change) at 7 days. **NOTE:** To perform this confirmation, do it in a humidified environment to prevent the medium from drying out.

Monitoring frequency

Follow the policies and procedures with the monitoring frequency specified by your country's professional associations and/or standards. As a recommended practice, and to provide optimal patient safety, Terragene® recommends that each sterilization load should be monitored with the appropriate biological indicator.

Storage conditions

Store preferably away from light, at a temperature between 10-30 °C and relative humidity between 30-80 %. Do not freeze. Do not store biological indicators near sterilizing agents or other chemicals.

Shelf life

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored under recommended conditions. Do not use indicators after their expiration date. The Process Indicator on the SCBI label has an expiration date of 2 years when used as part of the SCBI.

Stability of the endpoint reaction: The endpoint of the Process Indicator will remain unchanged for 6 months when the SCBI is stored under the conditions stated above.

Disposal

Discard biological indicators according to your country's sanitary regulations. Positive biological indicators can be sterilized in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, at 132 °C for 15 minutes or 134 °C for 10 minutes, or in a steam sterilizer with dynamic air removal at 132 °C for 4 minutes, or 135 °C for 3 minutes.

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