Process Challenge Device (PCD)

Quality certification

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STEAM

LOT

BT20 SCBI

Geobacillus stearothermophilus ATCC 7953

LOT

Heat shock population / Población _____ CFU / UFC

D-value / Valor D (121 °C) (*) _____ min.

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Survival time = (log₁₀ labeled population - 2) x labeled D-value

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Kill time = (log₁₀ labeled population + 4) x labeled D-value

Z-value / Valor Z _____ °C

Parameters determined at time of manufacture according to ISO 11138-1: 2017, ISO 11138-3: 2017 and IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant
ATCC is a registered trademark of American Type Culture Collection.
Terragene® and Bionova® are registered trademarks of Terragene S.A.

* Determined by the fraction negative method

PCDBI-2-RC Chemical Integrator

LOT

Stated Values (SV)

Temperature	Time
121 °C	_____ min.
128 °C	_____ min.
135 °C	_____ min.

SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

According to ISO 11140-1:2014 standard

Lic. Adrián J. Rovetto
Technical Director

SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

According to ISO 11140-1:2014 standard

Lic. Adrián J. Rovetto
Technical Director

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Explanation of Symbols

	Product designed for use in Steam sterilization cycles.		Batch number.
	CE mark.		Manufacture Date.
	Authorized representative in the European Community.		Expiration Date.
			Manufacturer.

ATTENTION: If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto** Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes. This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through **Bionova®** Software. Use **Trazanto** scanner to interpret the results (see instructions for use in the manual of the device).

ATTENTION: Do not write on the indicator before scanning.

EN

Process Challenge Device

For a rapid and easy monitoring of Steam sterilization processes

Indications for use

Bionova® PCD20-2 Steam Process Challenge Device has been designed for quick and easy monitoring of dynamic air removal and gravity air displacement Steam sterilization processes at 132/135 °C ≥ 4 minutes and at 121 °C ≥ 30 minutes.

Device description

Bionova® PCD20-2 Process Challenge Device (PCD) Test Pack has been designed to simulate a pack to be sterilized and to pose a challenge to the sterilization process. It is used to evaluate the effective performance of the process by detecting inadequate air removal and steam penetration. It also allows release of routine loads, especially implants. Furthermore, it allows to perform the routine monitoring and periodic validation of the sterilizers (after repair, installation, relocation).

Bionova® PCD20-2 consists of a disposable pre-assembled package as outlined in ANSI/AAMI ST79 which contains a BT20 Self-Contained Biological Indicator (SCBI) and a chemical integrator indicator (Type 5 according to ISO 11140-1:2014 standard) printed on a PCDBI-2-RC Record Card that gives instant visible indication that sterilizing conditions have been reached.

Each pack consists of a stack of porous cards holding a self-contained biological indicator tube that contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores soaked on a carrier as well as growth indicator medium contained in a glass ampoule. Each SCBI has a process indicator (Type 1 according to ISO 11140-1:2014 standard) on label that changes from pink to brown when exposed to steam. The chemical integrator printed on Record Card changes from purple to green when sterilization conditions were reached while process indicator (Type 1 according to ISO 11140-1:2014 standard) on PCD box changes from light blue to grey when exposed to steam.

Precautions

WARNING: Do not use PCD for monitoring Ethylene Oxide, Dry Heat, Formaldehyde or any sterilization process other than Steam. Do not reuse biological indicators.

WARNING: Place one or more PCDs in sterilizing hard-to-reach areas to ensure all areas of the chamber are sterilized. Evaluate all load configurations to ensure ALL hard-to-reach areas have been identified, and place a PCD in each of those locations.

WARNING: Do not reuse the sterilizer until the biological indicator test result is negative.

Instructions for use

1. Place the pack inside a normally loaded steam autoclave, in those areas which are considered most inaccessible for the sterilizing agent (e.g., the center of the load and areas near the door).

2. Run the sterilization cycle.

3. After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the test pack. **NOTE:** The color of the box may vary from the original after undergoing the sterilization cycle. This does not represent a problem regarding the operation or quality of the product.

4. Check that the process indicator printed on box has changed color from light blue to grey. Open the test pack, wait 5 minutes and remove the SCBI. Allow it to cool down to room temperature. **PRECAUTION:** Wear safety glasses and gloves when removing the biological indicator from the sterilized test pack. **WARNING:** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.

5. Check the chemical integrator on the Record Card for correct exposure. Color change to green confirms that the inside of the pack has been exposed to correct sterilization conditions. For chemical integrator reference color, please refer to Result Reference Guide. Otherwise, check the sterilization process.

6. Check the process indicator printed on SCBI's label. A color change to brown confirms that the biological indicator has been exposed to steam. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility.

7. Identify the Bionova® BT20 SCBI by writing the sterilizer number (in case of having more than one), load number, and processing date on the label. Fill out the required information on the Record Card.

8. Crush the glass ampoule contained in the SCBI with an individual ampoule crusher or with the ampoule crusher placed within the incubator's incubation area. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control in order to ensure that correct incubation conditions were met; capability of culture medium to promote rapid growth; no alteration of spore's viability due to improper storage temperature, humidity or proximity to chemicals and also correct functioning of Terragene® Bionova® Auto-reader Incubator. Both, the positive control indicator and the processed indicator should belong to the same batch.

9. Incubate the processed SCBI and the indicator used as positive control in Terragene® Bionova® auto-reader incubators for 24 hours between 55-62 °C. Color change of the culture medium from purple to yellow means that a sterilization process failure has occurred. If after 24-hour incubation there is no color change in the processed indicators, a final negative result is obtained (the sterilization process was effective). The positive control indicator should show culture medium color change to yellow for the results to be valid.

10. Record the SCBI and integrator results and adhere the self-adhesive Record Card or, alternatively, only the area containing the integrator indicator. **WARNING:** Do not use the sterilizer until the biological indicator test results are negative.

11. Discard the pack and the indicators immediately.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Storage

Store in a dark place, at temperatures between 10-30 °C, 30-80 % relative humidity. Do not freeze. Do not store near sterilizing agents or other chemical products.

Shelf life

PCD20-2 has an expiration date of 2 years from the date of manufacture, given by SCBI that carries, when stored at recommended conditions. Do not use after expiration date. Chemical Integrators and Process Indicators have an expiration date of 2 years when used in/on PCD.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicators can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes

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