PCD222-1

Rev. 10 / 03.2020

Process Challenge Device (PCD)

Quality certification

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LOT _____

BT222 SCBI

Geobacillus stearothermophilus ATCC 7953

LOT _____

Heat shock population / Población _____ CFU / UFC

D-value / Valor D (121 °C) (*) _____ min.

Survival time / Tiempo de sobrevida _____ min.
Survival time = (log₁₀ labeled population - 2) x labeled D-value

Kill time / Tiempo de muerte _____ min.
Kill time = (log₁₀ labeled population + 4) x labeled D-value

D-value / Valor D (132 °C) (*) _____ sec./seg.

Survival time / Tiempo de sobrevida _____ min.
Survival time = (log₁₀ labeled population - 2) x labeled D-value

Kill time / Tiempo de muerte _____ min.
Kill time = (log₁₀ labeled population + 4) x labeled D-value

D-value / Valor D (135 °C) (*) _____ sec./seg.

Survival time / Tiempo de sobrevida _____ sec./seg.
Survival time = (log₁₀ labeled population - 2) x labeled D-value

Kill time / Tiempo de muerte _____ min.
Kill time = (log₁₀ labeled population + 4) x labeled D-value

Z-value / Valor Z _____ °C

Parameters determined at time of manufacture according to ISO 11138-1:2017, ISO 11138-3:2017 and IRAM 37102:1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant
ATCC is a registered trademark of American Type Culture Collection.
Terragene® and Bionova® are registered trademarks of Terragene S.A

* Determined by the fraction negative method

PCDBI-1-RC Chemical Integrator

LOT _____

Stated Values (SV)

Temperature	Time
121 °C	_____ min.
128 °C	_____ min.
135 °C	_____ min.

SV = stated value. Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

According to ISO 11140-1:2014 standard

Lic. Adrián J. Rovetto
Technical Director

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Explanation of Symbols

Product is designed for use with Steam sterilization cycles.	Manufacture Date.
Batch number	Expiration Date.
	Manufacturer.

EN Process Challenge Device

For a rapid and easy monitoring of Steam sterilization processes

Composition

Each pack consists of a stack of porous cards holding a self-contained biological indicator (SCBI) tube that contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores soaked on a carrier as well as growth indicator medium contained in a glass ampoule. The pack also contains a Record Card with a self-adhesive Type 5 Chemical integrator indicator and, on the outside, a Type 1 process indicator.

Product description

PCD is a device designed to simulate a load to be sterilized, and to pose a challenge to the sterilization process. These devices are used to evaluate the effective performance of the process. The Steam Process Challenge Device has been designed for a rapid and easy monitoring of dynamic air removal Steam sterilization processes at 132/135 °C \geq 4 minutes and gravity air displacement processes at 121 °C \geq 30 minutes. It detects the inadequate air removal and

Steam penetration in Steam sterilizers. This device allows release of routine loads, specially implants. Furthermore, it allows to perform the routine monitoring and periodic validation of the sterilizers (after repair, installation, relocation).

The device consists of a disposable pre-assembled package as outlined in ANSI/AAMI ST79, and consists of a Self-Contained Biological Indicator (SCBI), a Type 5 chemical integrator indicator printed on the Record Card. The chemical integrator on the Record Card indicates that sterilizing conditions were met.

Warning!

Do not use PCD for monitoring Ethylene Oxide, Dry Heat, Formaldehyde or any sterilization process other than Steam.
Do not re-use the biological indicators.

Storage

Store in a dark place, at 10-30 °C temperature, 30-80 % Relative Humidity. Do not freeze.
Do not store biological indicators near sterilizing agents or other chemical products.

Instructions for use

- Place the pack inside a normally loaded steam autoclave. Place it in areas which you consider a priori most inaccessible for the sterilizing agent (high temperature steam). Generally, the center of the load and areas on the bottom shelf near the door and over the drain are problematic.
- Run the sterilization cycle:
132/135 °C dynamic-air-removal steam sterilization processes \geq 4 minutes
121 °C gravity displacement steam sterilization processes \geq 30 minutes
- After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the test pack. **NOTE:** The color of the box design may vary from the original color after undergoing the sterilization cycle. This does not represent a problem regarding the operation or quality of the product.
- Check that the external process indicator color on the outside of the test pack has changed. Open the test pack, wait 5 minutes and remove the SCBI. Allow it to cool down to room temperature. **PRECAUTION:** Wear safety glasses and gloves when removing the biological indicator from the sterilized test pack. **WARNING:** Do not crush or handle the biological indicator excessively, since this might cause the glass ampule to burst.
- Check the chemical integrator on the Record Card for correct exposure. Color change to the reference confirms that the inside of the pack has been exposed to sufficient sterilization conditions. Otherwise, check the sterilization process.
- Check the process indicator on the label of the biological indicator. A color change to brown confirms that the biological indicator has been exposed to steam. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility.
- Identify the Bionova® BT222 SCBI by writing the sterilizer number (in case of having more than one sterilizer), load number, and processing date on the label. Fill out the required information on the Record Card.
- Crush the glass ampule contained in the SCBI. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spores carrier entirely. Incubate at (60 \pm 2) °C in a compatible Bionova® autoreactor/incubator. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control whenever processed biological indicator is included. The positive control ensures that appropriate incubation conditions are met. Both the positive control indicator and the processed indicator should belong to the same batch. The indicator used as positive control must yield a positive result.
- Incubate the processed biological indicator and the indicator used as positive control for 1 hour at (60 \pm 2) °C to get the final fluorescence result. A 48-hour readout is optional to confirm the 1-hour result by visual color change. A positive fluorescence result (or growth medium color change during the 48-hour optional incubation) of the sterilized biological indicator means a sterilization process failure has occurred. If a negative fluorescence result is obtained (or growth medium remains the original color after the incubation period) the sterilization process was satisfactory.
- Wait for the final results. Fill out the SCBI results and adhere the self-adhesive whole Record Card or alternatively, only the area containing the integrator indicator. **WARNING:** Do not use the sterilizer until the biological indicator test results are negative.
- Discard the test pack and SCBI immediately.

Disposal

Dispose used biological indicators according to your country's healthcare and safety regulations. The positive biological indicator can be sterilized at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes in a gravity air displacement steam sterilizer, or at 132 °C for 4 minutes or 135 °C for 3 minutes in a dynamic air removal steam sterilizer.

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