Process Challenge Device (PCD)

Quality certification

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STEAM

LOT

CE

MD

BT222 SCBI

Geobacillus stearothermophilus ATCC® 7953

LOT

Heat shock population / Población _____ CFU / UFC

D-value / Valor D (121 °C) (*) _____ min.

Survival time / Tiempo de sobrevida _____ min.
Survival time = $(\log_{10} \text{ labeled population} - 2) \times \text{labeled D-value}$

Kill time / Tiempo de muerte _____ min.
Kill time = $(\log_{10} \text{ labeled population} + 4) \times \text{labeled D-value}$

D-value / Valor D (132 °C) (*) _____ sec./seg.

Survival time / Tiempo de sobrevida _____ min.
Survival time = $(\log_{10} \text{ labeled population} - 2) \times \text{labeled D-value}$

Kill time / Tiempo de muerte _____ min.
Kill time = $(\log_{10} \text{ labeled population} + 4) \times \text{labeled D-value}$

D-value / Valor D (135 °C) (*) _____ sec./seg.

Survival time / Tiempo de sobrevida _____ sec./seg.
Survival time = $(\log_{10} \text{ labeled population} - 2) \times \text{labeled D-value}$

Kill time / Tiempo de muerte _____ min.
Kill time = $(\log_{10} \text{ labeled population} + 4) \times \text{labeled D-value}$

Z-value / Valor Z _____ °C

Parameters determined at time of manufacture according to ISO 11138-1: 2017, ISO 11138-3: 2017 and IRAM 37102: 1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant
Terragene® and *Bionova®* are registered trademarks of Terragene S.A.
ATCC® is a registered trademark of American Type Culture Collection.

* Determined by the fraction negative method

IT26-C Chemical Integrator

LOT

Stated Values (SV)

Temperature	Time
121 °C	_____ min.
128 °C	_____ min.
135 °C	_____ min.

SV – stated value: Value of a critical process variable at which the indicator is designed to reach its endpoint as defined by the manufacturer.

According to ISO 11140-1:2014 standard

Lic. Adrián J. Rovetto
Technical Director

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Explanation of Symbols

	Product designed for use with Steam sterilization cycles.		Medical Device.
	CE mark.		Batch number.
	Authorized representative in the European Community.		Manufacture Date.
			Expiration Date.
			Manufacturer.


ATTENTION: If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto** Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes.
This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through **Bionova® Software**. Use **Trazanto** scanner to interpret the results (see instructions for use in the manual of the device).
ATTENTION: Do not write on the indicator before scanning.

Intended Use Table


Model	Autoclave/Steam Cycles		
	Gravity Displacement		
PCD222-C	121 °C	132 °C	135 °C
	30 minutes	15 minutes 25 minutes	10 minutes
	Dynamic Air Removal (Vacuum Assist)		
	121 °C	132 °C	135 °C
	NA	4 minutes	3 minutes
	Fluorescence Read Time	pH Color Change	
1 hour	48 hours		

Result reference guide

Unprocessed

Accept

Reject

EN Process Challenge Device

For a rapid and easy monitoring of Steam sterilization processes

Indications for Use

United States

Terragene® Bionova® PCD222-C provides a defined challenge resistance against the claim cycle shown below and also demonstrated resistance equivalence to the AAMI/ANSI 16 towel pack. The device provides routine monitoring and sterilizer qualification testing steam sterilization processes. See Intended Use Table for reference.

Outside the United States

Bionova® PCD222-C Steam Process Challenge Device has been designed for quick and easy monitoring of dynamic air removal and gravity air displacement Steam sterilization processes at 132/135 °C ≥ 4 minutes and at 121 °C ≥ 30 minutes.

Device Description

Bionova® PCD222-C Process Challenge Device (PCD) Test Pack has been designed to simulate a load to be sterilized and to pose a challenge to the sterilization process. It is used to evaluate the effective performance of the process by detecting inadequate air removal and steam penetration. It also allows release of routine loads, especially implants. Furthermore, it allows to perform the routine monitoring and periodic validation of the sterilizers (after repair, installation, relocation).

Bionova® PCD222-C consists of a disposable pre-assembled package as outlined in ANSI/AAMI ST79 which contains a BT222 Self-Contained Biological Indicator (SCBI), a PCDBI-C-RC Record Card and an IT26-C moving front chemical integrator (Type 5 according to ISO 11140-1:2014 standard) that gives instant visible indication that sterilizing conditions have been reached.

Each pack consists of a stack of porous cards holding a self-contained biological indicator tube that contains a population of *Geobacillus stearothermophilus* ATCC® 7953 spores soaked on a carrier as well as growth indicator medium contained in a glass ampoule. Each SCBI has a process indicator (Type 1 according to ISO 11140-1:2014 standard) on label that changes from pink to brown when exposed to steam. The moving front chemical integrator shows ACCEPT result when sterilization conditions were reached while process indicator (Type 1 according to ISO 11140-1:2014 standard) on PCD box changes from light blue to grey when exposed to steam.

Precautions

WARNING: Do not use PCD for monitoring Ethylene Oxide, Dry Heat, Formaldehyde or any sterilization process other than Steam. Do not reuse biological indicators.

WARNING: Place one or more PCDs in sterilizing hard-to-reach areas to ensure all areas of the chamber are sterilized. Evaluate all load configurations to ensure ALL hard-to-reach areas have been identified, and place a PCD in each of those locations.

WARNING: Do not reuse the sterilizer until the biological indicator test result is negative.

Instructions for Use

- Place the pack inside a normally loaded steam autoclave, in those areas which are considered most inaccessible for the sterilizing agent (e.g., the center of the load and areas near the door).
- Run the sterilization cycle.
- After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the test pack. **NOTE:** The color of the box may vary from the original after undergoing the sterilization cycle. This does not represent a problem regarding the operation or quality of the product.
- Check that the process indicator printed on box has changed color from light blue to grey. Open the test pack, wait 5 minutes and remove the SCBI. Allow it to cool down to room temperature. **PRECAUTION:** Wear safety glasses and gloves when removing the biological indicator from the sterilized test pack. **WARNING:** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.
- Check IR26-C moving front chemical integrator for correct exposure. If the dark bar has reached the ACCEPT zone, this confirms that the inside of the pack has been exposed to correct sterilization conditions.
- Check the process indicator printed on SCBI's label. A color change to brown confirms that the biological indicator has been exposed to steam. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility.
- Identify the Bionova® BT222 SCBI by writing the sterilizer number (in

case of having more than one), load number, and processing date on the label. Fill out the required information on the Record Card.

8. Press the cap to seal the tube. Crush the glass ampoule contained in the SCBI with an individual ampoule crusher or with the ampoule crusher placed within the incubator's incubation area. Then shake the tube down vigorously, until the medium reaches the base of the tube and soaks the spore carrier entirely. Incomplete wetting of the spore carrier may lead to an incorrect fluorescence readout. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control in order to ensure that correct incubation conditions were met; capability of culture medium to promote rapid growth; no alteration of spore's viability due to improper storage temperature, humidity or proximity to chemicals and also correct functioning of Terragene® Bionova® Auto-reader Incubator. Both, the positive control indicator and the processed indicator, should belong to the same batch.

9. Incubate the processed biological indicator and the indicator used as positive control in Terragene® Bionova® Auto-reader Incubators for 1 hour at $(60 \pm 2) ^\circ\text{C}$ to get the final fluorescence result. A 48-hour readout is optional to culture medium color change after the 48-hour incubation) means failure of sterilization process. If a negative fluorescence result is obtained (or culture medium remains the original color after 48-hour incubation), it means that sterilization process was successful.

10. Record the SCBI and integrator results and adhere the self-adhesive Record Card or, alternatively, only the area containing the integrator indicator. **WARNING:** Do not use the sterilizer until the biological indicator test results are negative.

11. Discard the pack and the indicators immediately.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Monitoring frequency

Follow the policies and procedures with the monitoring frequency specified by the professional associations and/or standards corresponding to your country. As recommended practice, and to provide optimal patient safety, Terragene® recommends that each sterilization load be monitored with the appropriate biological indicator.

Storage

Store in a dark place, at temperatures between 10-30 °C, 30-80 % relative humidity. Do not freeze. Do not store near sterilizing agents or other chemical products.

Self Life

PCD222-C has an expiration date of 2 years from the date of manufacture, given by SCBI that carries, when stored at recommended conditions. Do not use after expiration date. Chemical Integrators and Process Indicators have an expiration date of 2 years when used in/on PCD.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicators can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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