Cleaning Indicators

Results Reference Guide

Unprocessed

Failures

- Inadequate chemical agents' concentration.
- Insufficient water pressure.
- Poor water quality, temperature or cycle duration.

Processed and satisfactory

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators. Images are reproduced for illustration purposes only, and do not necessarily match the article described in every detail.

ATTENTION: If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto Automatic Quality Control and Traceability System** for cleaning, disinfection and sterilization processes. This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through **Bionova® Software**. Use **Trazanto** scanner to interpret the results (see instructions for use in the manual of the device). **ATTENTION:** do not write on the indicator before scanning.

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EN Cleaning Indicators

For Endoscopes Cleaning Performance Test for Automated Endoscope Reprocessors

Endoscopes are very complex reusable medical devices. They consist of several channels and tubes which can be 2 meters long, with variable internal diameters. Endoscopes get in contact with human mucous membranes and tissues, consequently, high-level cleaning and disinfection is required after its use, in order to avoid pathogen transmission. Automated Endoscope Reprocessors (AERs) are the preferred method for endoscope washing and disinfection given that they provide a standardized method with minimal human intervention, together with the fact that washing and disinfection cycles are highly reproducible. An optimal result depends on the adequate combination of water pressure and quality, chemical agents' concentration, temperature and cycle duration. Monitoring the efficacy of cleaning processes is crucial in order to ensure compliance with the standards of cleanliness.

Product description

The cleaning indicator Chemdye® CDWE consist of a 38 mm x 6 mm synthetic, self-adhesive, thermostable support printed with a mixture of specially combined natural colored components. The indicator's formulation allows detection of all factors affecting the outcome of the endoscopes cleaning process in Automatic Endoscope Reprocessors. The indicator ink is non-toxic.

Instructions for use

1. Select for the Chemdye® CCDER device of your preference, according to your endoscope and/or reprocessor configuration.
2. Place a Chemdye® CDWE indicator strip inside the capsule, with the indicator ink towards the perforated side.
3. Screw the capsule lid.
4. Connect the Chemdye® CCDER device to a proper AER terminal.
5. Connect the endoscope terminals to proper AER terminals.
6. Run the selected washing program.
7. Once the process is finished, disconnect the Chemdye® CCDER device from the AER.
8. Remove the indicator strip and evaluate the result.
9. Record the results by sticking the self-adhesive strip in your record book.

Results interpretation

Complete elimination of the washing ink indicates that cleaning efficacy was

adequate. If ink residues are observed, the endoscope cannot be considered clean. Failures during a cleaning cycle can be a consequence of inadequate chemical agents' concentration, insufficient water pressure, poor water quality, temperature or cycle duration. All these events diminish cleaning efficacy and it is not recommended to use the endoscope until a satisfactory result is obtained.

NOTE: Endoscope reprocessing monitoring can be complemented with chemical indicators for peracetic acid as disinfectant, the Integron® IT400 and IT401 series. Chemdye® CDWE and CCDER devices allow monitoring endoscope cleaning efficacy but not disinfection level. For disinfection level monitoring use the biological indicator Bionova® BT400, following the instructions for use of this product.

Additional Technical Information

Chemdye® Cleaning Indicator CDWE is manufactured, calibrated and designed following specifications outlined in ISO/TS 15883-5:2005.

CAUTION: Terragene® does not guarantee a proper performance of the washing indicator Chemdye® CDWE if used with Cleaning Challenge Devices of other manufacturers.

Results obtained with the washing indicator Chemdye® CDWE are only valid for the washing cycle under analysis. Conclusions for previous or subsequent cycles should not be drawn.

Terragene® does not guarantee the proper functioning of washing indicator Chemdye® CDWE if used for different washing processes than the ones it was designed for.

Storage

Store the washing indicator Chemdye® CDWE in a dry place, protected from light, at a temperature from 10-30 °C and 30-80% relative humidity. Do not use the indicators after their expiration date. **NOTE:** store the product incorrectly may affect the stability thereof and the results obtained.

Disposal

Discard the indicators after use according to your country's healthcare and safety regulations.

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