# Ultrasonic Cleaning Indicator

## Results reference guide

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### Explanation of Symbols

	CE mark.		Batch number.
	Authorized representative in the European Community.		Manufacture Date.
	Medical Device.		Manufacturer.

## EN Cleaning indicators

For cavitation capacity testing of ultrasonic washers

Control of cleaning and washing process of medical instruments and materials is very important as this result influences the success of subsequent disinfection and / or sterilization processes. Ultrasonic washers are a key point in this process as they allow the removal of soils in areas of difficult access such as joints, crevices and lumens of medical instruments. The cleaning power of ultrasonic washers is due to the cavitation process, where small vacuum bubbles at high temperatures collapse quickly creating liquid flows at high pressure that impact on the instruments removing the residues present in their surface.

### Product description

Chemdye® CDWU indicator was designed to monitor the cavitation process. It allows testing the operation of ultrasonic washing machines with different washing configurations by measuring the generated cavitation energy. CDWU indicator consists of a clear vial, with a reactive blue-colored solution and glass beads immersed in it. When the cavitation energy is adequate, vibration of the glass pearls triggers a color change in the solution, from blue to yellow, through a range of color intermediates.


### Characteristics

CDWU indicator does not require a special holder to be used. Simply, place one vial on the ultrasonic machine zone you wish to monitor and ultrasonic washer cycle is started. CDWU indicator can be used for monitoring ultrasonic cleaning processes in washing machines operating in a temperature range of 18 °C to 70 °C, with detergent or washing solution, and at 35 kHz or higher frequencies.

### Instructions for use

Ultrasonic washers should be checked at least once per week and the results should be recorded. There are two types of tests (see figure) and CDWU indicator can be used in both:  
**-Periodic functional test:** for initial set-up of the washing machine, and at regular intervals of three months or after a repair. The purpose of this test is to verify the uniform distribution of the cavitation process in the washer tank.  
**-Test routine:** is performed weekly or daily, to verify that the cavitation works correctly. This test allows detecting a possible failure in the performance of the ultrasonic washing machine.

### Periodic functional test

### Procedure:

- 1-Prepare the tank with cleaning solution following both the manufacturer's instructions of the ultrasonic washing machine and the manufacturer's instructions of the detergent used.
- 2-Degassify the tank according to indications of use of the detergent manufacturer.
- 3-Adjust the temperature according to indications of use of the detergent manufacturer.
- 4-Take the necessary number of indicators according to the Test you want to perform (Periodic Functional Test or Routine Test) and the size of the washer tank.
- 5-Place them in the basket without load according to the distribution of the Test to be performed.
- 6-Place the basket in the tank and run the desired ultrasonic cleaning program.

- 7- Completed the program, remove the indicators from the machine.
- 8- Analyze the results according to the table and record the data obtained.

### Interpretation of results

Not exposed Cavitation failures Correct cavitation

The color change from blue to yellow is an accepted result and this indicates the presence of cavitation energy suitable for washing. The time in which this color change occurs is a measure of the strength of cavitations energy. Green color result indicates a failure to reach the cavitation energy necessary for a good cleaning.

In both tests, the slower change in color of an indicator (in average) indicates the presence of an area with less cavitation energy or a "weak point" of the tank, while a failure result indicates the presence of a "cold spot" or an area without the cavitation energy necessary for a good wash. In addition to the correct operation of the ultrasonic washer, the cavitation process is influenced by multiple variables such as: presence of dissolved gases in the water, water hardness, temperature and the use of detergents. If there is an unfavourable result, check:

- Working frequency of the ultrasonic washing machine. CDWU indicator can be used in ultrasonic washers with frequencies equal to or greater than 35 kHz. Frequencies greater than 40 kHz cause a faster color change.
  - Presence of bubbles in the washing solution. That occurs because of the dissolved gases presence in the tank solution of the ultrasonic washer. The degassing of the cleaning solution is extremely important to achieve a result a satisfactory result. This process can take several minutes or even hours depending on the temperature and the tank's size. Follow the manufacturer's instructions for proper degassing.
  - Time and temperature of the ultrasonic wash. CDWU indicator is designed to verify cavitation during ultrasonic washing with a minimum duration of 3 minutes (180 seconds) and temperatures between 18-70 °C.
  - General condition of the ultrasonic machine. Any visible deterioration of the washer or unknown sound detected could mean that the machine is not working correctly.
- After checking the above actions, repeat the test. Report any unsatisfactory results for corrective action in accordance with the policies of your Institution. For questions regarding the operation of the washing machine, in order to corroborate possible causes of faults, contact the technical service.

### Precautions

Do not re-use the CDWU indicator.  
 Do not open the indicator before its use, since this will cause its deterioration.  
 Do not use the CDWU indicator for monitoring normal washing using non-ultrasonic cleaners.  
 At 24 hours, the color of the CDWU indicator might vary. It is recommended to read and register the results when performing the test.

### Storage

CDWU indicator should be stored at temperatures between 10-30 °C, specially protected from light and moisture, within its original packaging.

### Disposal

Discard the indicators after use according to your country's healthcare and safety regulations.

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