
Hygiene Monitoring System

Figure 1

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Product description

The first step for a proper cleaning of surgical instruments is to eliminate the blood, body fluids and tissues immediately after use. Biological contaminants such as blood generate high hygienic risk in both patients and health personnel, which can result in the transmission of pathogens affecting large numbers of people.

Chemdye® HemoSwab Hygiene System was designed to verify surgical instruments cleaning by detection of blood residues left behind after improper cleaning. The system is based on the high sensitivity of peroxidase activity of the hemoglobin protein which allows a rapid detection up to 0,1 µg of blood.

This system has a high absorption swab in conjunction with a moisturizer that allows collection of samples from different surfaces with the same efficacy.

A visual readout by a color change indicates detectable levels of hemoglobin.

Indications of use

Use for testing any surface after the cleaning process and for hard-to-reach areas in complex instruments such as serrated edges or box joints.

Characteristics

Chemdye® HemoSwab provides all the equipment required for testing blood residues on surgical instruments or surfaces.

Precautions

Do not open until use. Do not reuse. Wear gloves, surgical mask and other hygiene measures. Do not freeze. Protect from high temperatures. Keep away from light. Oxidizing agents such as hypochlorite (present in some detergents and disinfecting agents) may cause false positives.

Instructions for use

1. Take a swab from the container without touching the absorbent part of any swab.

2. Moisten the swab by immersing for 5 seconds in the moisturizer solution (vial with light blue cap) (Figure 1). Press the swab against the walls of the tube to remove moisturizer's excess.

3. Slip the swab thoroughly onto the surface where you want to take the sample.

NOTE: It is recommended to slide the swab zigzag in one direction and then zigzag in the perpendicular direction. Press the swab and rotate it several times while collecting the sample. To perform comparable cleaning tests, surface sampling should be standardized. So, it is recommended to always take the sample in areas of similar surface size and slide the swab 10 times in each direction.

4. Unscrew the cap of the indicator vial (transparent cap) and transfer the liquid to the reaction activator vial (blue cap) (Figure 2). Close the reaction activator vial and mix by inversion five times (Figure 3).

5. Immerse the sample swab in the activated indicator solution (vial blue cap) (Figure 4).

6. After 30-60 seconds observe if color change occurred on the surface of the swab or in the indicator solution. The color change from colorless to bluish green indicates the presence of blood residues. The more blue and more intense the color developed, the greater the amount of blood present on the surface tested.

7. After using the HemoSwab system reprocess the instrument or surface that was used for blood monitoring.

Special considerations and precautions

-If the product is stored refrigerated let it take room temperature before use.

- If the hydration or immersion procedure was incorrect, the analysis could be affected.

- The results analyzed after the recommended time period are not valid. Small amounts of blood can give a color change that is not stable over time. In addition the activated solution can change to a soft green color after a few hours without this meaning the presence of blood.

- Oxidizing agents such as hypochlorite (present in some detergents and disinfecting agents) can cause false positives. In these cases do not use this system to detect blood residues.

Storage

The Chemdye® HemoSwab system should be stored at temperatures between 2 and 25 °C and keep away from light. Shelf-life is 18 months from the date of manufacture when stored under recommended conditions.

Disposal

Discard the vials after use, according to your country healthcare and safety regulations. Due to the possibility of the vial containing blood, it is recommended that the HemoSwab vials used to detect blood residues be discarded in a biological waste container.