



Indicator Ink

Instructions for use

Use exclusive for professionals e Instituciones Sanitarias.

Explanation of Symbols

Batch number.

Manufacture Date.

Expiration Date.

Manufacturer.

EN Indicator Ink

For washing efficiency tests in washer-disinfectors

Product description

Chemink® IKW4 washing Indicator Ink contains a mixture of colored components (proteins, sugars and colored fatty acids) specially designed to monitor the washing efficiency in thermal disinfection washing machines. IKW4 poses a high challenge to the washing process and it can be removed from metallic surfaces. Washability depends on the presence, type and concentration of detergent; temperature and pressure of the washing solution and cycle time. The indicators efficiency is validated by comparing against the test soils detailed in ISO/TS 15883-5:2005 standard (Annex N). The IKW4 ink could also be used for validating cleaning processes in ultrasonic washing machines operating at frequencies equal to or greater than 35 kHz, between 20 °C and 65 °C, in synergy with detergents suitable for this process.

IMPORTANT: Chemink® IKW4 Indicator Ink should be used as a validation tool for washing processes. Its application, the amount placed and sites where adhered depend on the methodologies applied by the user. It is recommended for use in difficult-to-wash areas of metallic surgical materials.

Instructions for use

- Shake before use.
- Remove the brush from the bottle and apply a thin film of indicating ink on the surface of the product to clean. **IMPORTANT:** apply on those areas that are inaccessible to the washing solution, for example, hinges of metallic surgical instruments, serrated edges or box joints. The thickness of the layer spread affects the result.
- Leave to air dry for 30 minutes.
- After ink is dry, place the material to clean on the tray and subject to the selected washing cycle.

5. Once the cleaning process has finished; remove the material from the washer.

6. Visually observe the presence of remaining indicator ink on the material. Analyze results.

NOTE: if any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Interpretation of results

Systematize the application method and validate it in order to use it as an important tool in washing machines performance analysis.

Complete removal of Chemink® IKW4 Ink indicates that the cleaning efficiency was adequate. If residues of the ink spot are observed, the material cannot be considered clean. Failure during a cleaning cycle can be a consequence of inadequate chemical agent concentration, insufficient water pressure, temperature or process duration. These factors affect the cleaning efficiency, therefore it is not recommended to use the material until a satisfactory result is obtained.

NOTE: for uncertain results, the cleaning process can be validated by using systems for protein detection on surfaces such as Chemdye® PRO1 MICRO or PRO2. In case of an unfavorable outcome, check the indicated steps and repeat the cleaning process.

Precautions

Do not use the product after its expiration date. Do not use Chemink® IKW4 Indicator Ink to monitor other processes different from indicated. In case of spill, clean with absorbent material. Do not wet.

WARNING: Professional or industrial use only. Not for domestic or cosmetic use. In case of contact with eyes or skin, rinse with water. Seek medical attention if irritation occurs.

Storage and product shelf-life

The indicator ink has a shelf-life of 12 months from the date of manufacture when stored at temperatures between 10-30 °C, relative humidity 30-80 %, kept in a dark place and away from sterilizing agents, vapors and/or gases.

Disposal

Discard along with glass waste according to your country's healthcare regulations.

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