Result reference guide

Unprocessed

Failures

One or more disinfection parameters were not achieved. The materials cannot be considered disinfected and must be processed again.

The main failure causes are:

- Insufficient time
- Insufficient temperature

Processed and Satisfactory

Explanation of Symbols

	CE mark.
	Authorized representative in the European Community.
	Medical Device.

	Batch number.
	Manufacture Date.
	Expiration Date.
	Manufacturer.

EN Thermodesinfection and cleaning indicators

For moist heat disinfection processes in washer-disinfectors

Product description

After cleaning, surgical instruments must be disinfected to ensure safe subsequent handling and processing. Thermal disinfection with moist heat, is the most common method for disinfection of medical devices in the hospital setting. Disinfection effects can be defined by the A_0 value, a physical parameter specified in ISO 15883-1, which denotes the inactivation of microorganisms. A particular A_0 value can be achieved with the most diverse temperature/time combinations. For processing surgical instruments, an A_0 value of 600 is specified in ISO 15883-2, corresponding to a hold time of 1 minute at 90 °C, ensuring lethality of bacteria, mycobacteria, fungi and heat-sensitive viruses. If efficacy is also to be ensured against heat-resistant viruses, e.g. hepatitis B, many authorities in the field of disinfection, like the Robert Koch Institute, recommend an A_0 value of 3000, corresponding to a temperature of 90 °C with a hold time of 5 min. A higher safety approach from this institute comprises a cycle at 93 °C, 10 minutes, a program known as the "BGA/RKI" or epidemic program.

Integron® Indicators have been designed to react in moist heat disinfection processes in washer-disinfectors, providing an accurate, convenient method of routine control, and assuring disinfection efficiency. The green indicating ink was developed to turn to purple when temperature and time conditions are met.

Characteristics

Integron® Indicators consist of one-sided-laminated strips of 30 mm x 70 mm printed with green reactive ink, designed with state-of-the-art technology to withstand conditions present during moist heat disinfection processes.

Indicators 100 % free of toxic heavy metals

Indications for use

Use the correspondent Integron® Indicator for monitoring thermal disinfection step in washer-disinfectors according to the following programs:

IT27W-1: 90 °C, 1 minute. A_0 = 600
 IT27W-5: 90 °C, 5 minutes. A_0 = 3.000
 IT27W-10: 93 °C, 10 minutes. BGA/RKI/Epidemic Program. A_0 = 12.000

Precautions

Do not expose this product to Steam, EO, Dry Heat, or any sterilization or disinfection processes other than Moist Heat Disinfection Processes in washer-disinfectors. Do not reuse.

WARNING: The results are not reproducible when immersing the indicators in a thermal bath.

Instructions for use

1. Place the indicator on the Chemdye® CDWAH Holder.
2. Place the system on the tray along with the material.
3. Run the selected disinfection program and once finished remove the indicator.
4. Visually analyze indicators and record results. The indicator must change to the reference color for considering that the desired disinfection conditions were met. If unsatisfactory results are obtained, perform the actions indicated therein and repeat the disinfection process.

Use one Integron® indicator for each cycle.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Interpretation of results

The color change of the indicating ink from green to purple indicates that the thermal disinfection conditions were attained, thus allowing a quick and easy assessment of the success of the process. If no color change is observed, or the color obtained does not reach the reference color, the load cannot be considered disinfected. Failures during a disinfection cycle can be a consequence of inadequate low temperature or cycle duration, which lowers its efficiency. It is not recommended to release the load until a satisfactory result is obtained.

Storage

Indicators should be stored in dry places, away from light, at temperatures between 0 °C and 4 °C, 30-80 % relative humidity.

Do not wet.

Do not store close to sterilizing agents.

Shelf life

Indicators have an expiration date of 5 years from the date of manufacture when stored under recommended conditions. The expiration date as well as the batch number of the product are specified on the package label and in each unit of the product.

Disposal

Discard indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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