
# IT401 Chemical Indicator


Rev.6 / 03.2023

Integron®

## Liquid Peracetic Acid

### Quality Certification

Code \_\_\_\_\_ Concentration \_\_\_\_\_ SV (Time) \_\_\_\_\_  
 \_\_\_\_\_ 2.500 ppm \_\_\_\_\_ min.  
 \_\_\_\_\_ 1.000 ppm \_\_\_\_\_ min.  
 \_\_\_\_\_

Lic. Adriano J. Rovetto  
Director Técnico  
Technical Director

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### Explanation of Symbols

 CE mark.

  Authorized representative in the European Community.

 Batch number.

 Manufacture Date.

 Expiration Date.

 Manufacturer.

### EN Chemical Indicator

For Liquid Peracetic Acid High Level Disinfection processes

#### Product description

Integron® IT401 Indicators are designed to react in High Level Disinfection processes that use Peracetic Acid (PAA) as disinfecting agent, assuring an adequate control of the disinfection process efficiency (temperature, time and PAA concentration). The reactive ink was developed for processes that work at temperatures between 35-43 °C and turns to blue when integration conditions reach the established values stated in the quality certificate.

#### Characteristics

Integron® IT401 Chemical Indicators consist of an indicator strip of 10 mm x 45 mm, printed with a brown reactive ink specific for PAA, contained within an envelope of non absorbent material of 25 mm x 70 mm, which avoids direct contact of the reactive ink with the aqueous medium.

Indicators 100 % free of toxic heavy metals

#### Indications for use

Use Integron® IT401 Chemical Indicators for monitoring Liquid Peracetic High Level Disinfection cycles in equipments working at 35-43 °C.

**ATTENTION:** Processes at higher temperatures could give results with bluish-green to green coloration. The same results could be obtained if exposure time is extended or Peracetic Acid concentration is above the values stated for that batch.

#### Precautions

Do not expose Integron® IT401 to Steam, EO, Dry Heat, or any sterilization or disinfection processes other than Liquid PAA.

#### Instructions of use

1. With the reusable silicon strip provided along with the product, locate the envelope with the chemical test strip inside the equipment, together with the load to be processed with PAA.
2. Process the load according to the established procedures.
3. After completing the process, analyze the results. The indicator must change to the reference color for considering that the desired disinfection conditions were reached.

**NOTE:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

#### Storage and shelf life

- Keep this indicator in dry places, protected from light and at temperatures between 10-30 °C and 30-80 % relative humidity. Do not wet. Do not storage close to sterilizing or disinfection agents.
- Indicators have an expiration date of 5 years from the date of manufacture when stored in the recommended conditions. Both expiration date and batch of the product are specified on box's label and on each unit of product. Do not use indicators after their expiration date.
- After use, the indicator will not change in 6 months, if stored at the temperature and relative humidity previously indicated.

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