
PRO2

Rev. 14 / 05.2021

Hygiene Monitoring System

KPRO2-E69 | KPRO2-E250

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	CE mark.		Manufacture Date.
	Authorized representative in the European Community.		Expiration Date.
	Batch number.		Manufacturer.

EN Hygiene Monitoring System

For detecting residual protein on surfaces

Product description

The first step in properly cleaning surgical instruments is to rinse off all blood, bodily fluids and tissue immediately after use. If cleaning is not carried out adequately, it may render the disinfection stage ineffective. This results patient's exposure to body fluids and tissue contaminants from prior subjects, which might result in the transmission of pathogens and affect large numbers of people.

Chemdye® KPRO2-E69 and KPRO2-E250 Hygiene Systems are designed to check cleanliness of surgical equipment by detecting protein residues after cleaning. The systems have a high absorption swab, allowing the collection of samples from different surfaces with the same efficacy. The systems are compatible with cleaning verification of endoscopes and other reusable instruments with hard to reach internal channels or cannulated instruments. For this, use the special SW250 endoscope swabs, which may be introduced through the biopsy channel of the endoscope to sample contaminants left after the cleaning process.

A visual color change indicates the presence of detectable levels of protein. KPRO2-E69 and KPRO2-E250 have high sensitivity and can detect as little as 1 µg of protein.

Characteristics

KPRO2-E69 and KPRO2-E250 provide all the necessary equipment required for HTM 01-01 / HTM 01-05 / ISO 15883-1 Residual Protein Test.

Precautions

Do not open until its use. Do not touch the swab tip with your hands. Do not reuse.

Wear gloves, surgical mask and other hygiene measures to avoid cross-contamination of protein residues from your skin. Do not freeze. High concentrations of detergents may cause false positive results. Although KPRO2-E69 and KPRO2-E250 have high sensitivity, they are not a microbiological monitoring device, and surfaces may not be considered to be sterile based on the results obtained.

Indications for use

- Take a swab from the foil pouch without making contact with the absorbent part.
- Moisturize the swab by immersing the swab tip into the moisturizer tube for 5 seconds before swabbing. Do not push the swab against the walls of the tube. Do not shake.
- Slide the swab thoroughly over the surface where you want to collect the sample.

A)When using SW250 endoscope swabs insert the swab into the endoscope channel and push it all the way through one time. Then, cut the tip of the swab (approximately 7 cm from the end) with a pair of scissors. Be careful not to cut the distal end. Finally, pull the wire from the distal end to remove the remaining swab segment from the biopsy channel. B)If you work under the ISO 15883 standard, it is recommended to sample approximately 10 cm² of surface. Swab in zigzag in one direction and then do the same perpendicularly. Apply pressure on the swab and rotate it while collecting the sample. For performing comparable cleaning tests, you must standardize surface sampling. In order to achieve this, it is advisable to take the sample always at similar points of the surface, with an invariable area of 10 cm². swab 10 times in each direction. C)If you work under the HTM 01-01 standard, pick one size of the surgical instrument to be tested and swab thoroughly over its surface. Apply pressure on the swab and rotate it while collecting the sample. 4.Unscrew the cap of the vial and immerse the swab tip (or the initial segment of the endoscope swab) in the reactive colorimetric solution, move it for 15 seconds into the solution.

IMPORTANT: Do not push the swab against the walls of the tube, just immerse and move the swab into the solution. If hydration or immersion procedure was incorrect, the analysis might be affected.

5.Discard the swab and screw the cap back on the reagent vial. Color development will be complete after 1 minute and will be stable during 20 minutes.

6.Interpret results visually by using the reference guide. By comparing the color of the test against the reference guide, an estimation of the surface cleanliness can be made:

IMPORTANT: Results analyzed beyond the recommended period are not valid.

7.After using the PRO2 system, reprocess the instrument or the surface used for the monitoring of proteins.

Storage

KPRO2-E69 and KPRO2-E250 systems should be stored at temperatures between 2-30 °C in dark places. Shelf life is 24 months when stored refrigerated and 6 months when stored at room temperature up to 30 °C.

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Terragene S.A.