Results reference guide

Unprocessed

Failures

- Inadequate chemical agents' concentration.
- Insufficient water pressure.
- Poor water quality, temperature or cycle duration.

Processed and satisfactory

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators. Images are reproduced for illustration purposes only, and do not necessarily match the article described in every detail.

ATTENTION: If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto Automatic Quality Control and Traceability System** for cleaning, disinfection and sterilization processes. This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through **Bionova® Software**. Use **Trazanto** scanner to interpret the results (see instructions for use in the manual of the device).

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

Figure 4

Cleaning challenge kit

For monitoring cleaning performance of endoscopes in Automatic Endoscope Reprocessors.

Endoscopes are highly-complex reusable medical devices. They have different channels, up to 2 meters long, with different inner diameters. Since these devices are in contact with human tissues and mucosa, they must undergo cleaning and high-level disinfection processes before being reused, in order to avoid the transmission of pathogens. Automatic Endoscope Reprocessors (AERs) offer the cleaning and disinfection method for excellence, since they require minimal user intervention, and also, have highly reproducible cleaning and disinfection cycles. An optimal result depends on the proper combination of water pressure and quality, effective concentration of chemical agents, temperature and cycle duration. The monitoring of cleaning performance is essential for guaranteeing compliance with cleaning standards.

Product description

Chemdye® Lumena SixFlow LSF2 cleaning challenge kit, which consists of a challenge device and Chemdye® Lumena CDWE cleaning indicators, allows to monitor cleaning efficacy in Automatic Endoscope Reprocessors (AERs). The challenge device consists of six tubes of polytetrafluoroethylene (PTFE), three 4 mm inner diameter and three 2 mm inner diameter, all of them 1.5 m long. Each tube has a special connector in one of their ends for plugging it to the Automatic Endoscope Reprocessor, and a device in the other end that simulates an endoscope's lumen and where two Chemdye® Lumena CDWE cleaning indicators can be placed (one inside and one outside). Each tube is identified at both ends with a color mark that allows to fastly recognize the Automatic Endoscope Reprocessor's output terminal to which it will be connected. Chemdye® Lumena CDWE cleaning indicator has been specially developed for monitoring these processes. It consists of a synthetic, self-adhesive and thermostable support of 38 mm x 6 mm, printed with a test soil made up from a mixture of specially combined colored natural components. Chemdye® Lumena CDWE cleaning indicator complies with ISO/TS 15883-5:2005 standard.

Chemdye® Lumena CDWE cleaning indicator allows the detection of failures that may affect both the internal and external cleaning efficacy of endoscopes in Automatic Endoscope Reprocessors. Both cleaning indicators' test soil must be completely washed (Figure 1). It is recommended to verify both internal cleaning efficacy and external cleaning efficacy, although it can be used one indicator and verify only the inner cleaning efficacy.

Precautions

- Do not reuse cleaning indicators.
- Do not touch the indicator ink with your hands.

Instructions for use

- Cleaning performance test
1. Check that the challenge device is clean, dry and working properly.
 2. Check that your hands are clean and dry.
 3. Unscrew the 6 devices' caps.
 4. Place a Chemdye® Lumena CDWE cleaning indicator in the inner cavity of each device in a way that the test soil is facing upwards (Fig. 2).
 - WARNING:** before placing the indicator in the inner cavity of the device, remove its liner. Secure the position of the indicator in the inner cavity using the tweezers provided, so that it is firmly adhered to the device (Fig. 3).
 5. Screw the 6 caps again, ensuring color matching of caps and devices.
 6. Place a Chemdye® Lumena CDWE cleaning indicator in the external slot of each device in a way that the ink is visible through its window (Fig. 4).
 - WARNING:** before placing the indicator in the external slot, remove its liner. Secure the position of the indicator in the external slot using the tweezers provided, so that it is firmly adhered to the device (Fig. 4).
 7. Connect the 6 PTFE tubes of the challenge device to the Automatic Endoscope Reprocessor through the output terminals.
 8. Run a normal cycle of the Automatic Endoscope Reprocessor and cancel it after the cleaning phase and just before the disinfection phase. If it is not possible to interrupt the cycle, check the result when the complete cycle has finished.
 9. Once the cycle has finished, unplug the connectors and remove the challenge device from the Automatic Endoscope Reprocessor.
 10. Carefully remove the indicators from the devices, by using the tweezers provided.
 11. Examine the indicators, by device color, in order to find remains of test soil by placing them on a white background. Record the results.
 - NOTE:** in case of having Trazanto traceability software, analyze the indicators using the scanner.
 12. If any of the indicators corresponding to a particular device color has test soil remains, record the corresponding color and take the applicable measures in case of failure of the cleaning performance test.

Obstruction test (channel permeability)

1. Check that the challenge device is clean, dry and working properly.
2. Unscrew the 6 devices' caps and reserve them.
3. Screw a grey cap (0.5 mm inner diameter) on each device for obstruction test.
4. Connect the 6 PTFE tubes of the challenge device to the Automatic Endoscope Reprocessor through the output terminals.
5. Run a normal cycle of the Automatic Endoscope Reprocessor.
6. The machine should detect an obstruction in each channel and trigger an alarm. Record the results.
7. Once the cycle has finished, unplug the connectors and remove the challenge device from the Automatic Endoscope Reprocessor.
8. Remove the grey obstruction caps.
9. For those channels where no obstruction is detected, record the channel number and start an investigation.

Interpretation of results

If Chemdye® Lumena CDWE cleaning indicator is completely washed, it means that the cycle conditions are appropriate for guaranteeing a correct endoscope's cleaning (both internally and externally). Therefore, the flexible endoscope can be reprocessed by using the recently validated program.

The failures on a cleaning cycle may result from an inadequate concentration of chemical agents, insufficient water pressure, inadequate water quality, temperature or duration of the cycle, among others. These facts decrease the cleaning efficacy, therefore it is not recommended to reprocess or use the endoscope until a satisfactory result is obtained.

NOTE: Chemdye® Lumena SixFlow LSF2 kit allows the monitoring of endoscopes' cleaning processes, but it does not make it possible to determine if the appropriate disinfection level has been reached. If the equipment does also perform a high-level disinfection by using peracetic acid as disinfecting agent, the monitoring can be completed through the use of Integron® IT400 and IT401 chemical indicators, as well as through Bionova® KBT400 biological indicator.

Additional technical information

Before each use, it must be verified that the connectors and the PTFE tubes are in perfect conditions. In case any connector is damaged and its replacement is justified, please contact a representative and request a spare part. In case the tube is blocked or it has water leaks, the challenge device must be replaced with a new one.

WARNING: Terragene® does not guarantee the well-functioning of Chemdye® Lumena SixFlow LSF2 Kit if other manufacturers' cleaning indicators or challenge devices are used. The results are only valid for the cleaning cycle tested with the particular parameters. Do not draw conclusions for cycles with different parameters. Terragene® does not guarantee the well-functioning of Chemdye® Lumena SixFlow LSF2 kit if it is used for cleaning processes different from the one for which it was originally designed for.

Storage

Store Chemdye® Lumena SixFlow LSF2 kit in a dry place, protected from light, at a temperature between 10-30°C and at a relative humidity between 30-80%. Do not use Chemdye® Lumena CDWE cleaning indicators after their expiration date.

NOTE: an incorrect storage of the product may affect its stability as well as the results obtained.

Disposal
Discard the kit components after their use, according to your country's healthcare and safety regulations.

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