
Challenge Kit for Endoscope Reprocessors

Results reference guide

Unprocessed

Failures

- Inadequate chemical agents' concentration.
- Insufficient water pressure.
- Poor water quality, temperature or cycle duration.

Processed and satisfactory

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators. Images are reproduced for illustration purposes only, and do not necessarily match the article described in every detail.

ATTENTION: If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto Automatic Quality Control and Traceability System** for cleaning, disinfection and sterilization processes. This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through **Bionova® Software**. Use **Trazanto** scanner to interpret the results (see instructions for use in the manual of the device).

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Explanation of Symbols

- MD** Medical Device.
- LOT** Batch number.
- Manufacture Date.
- Expiration Date.
- Manufacturer.

Figure 1


Figure 2


Figure 3


Figure 4


Figure 5

	Chemdye® LUMENIA	Terragene®	Trazanto®
LSF1 R.0			
CAPSULE			
ORANGE			
GREEN			
BROWN			
YELLOW			
BLACK			
WHITE			

Cleaning challenge kit

For monitoring cleaning performance of endoscopes in Automated Endoscope Reprocessors

Endoscopes are highly-complex reusable medical devices. They have different channels, up to 2 meters long, with different inner diameters. Since these devices are in contact with human tissues and mucosa, they must undergo cleaning and high-level disinfection processes before being reused, in order to avoid the transmission of pathogens. Automated Endoscope Reprocessors (AERs) offer the cleaning and disinfection method par excellence, since they require minimal user intervention, and, also, have highly reproducible cleaning and disinfection cycles. An optimal result depends on the proper combination of water pressure and quality, an effective concentration of chemical agents, temperature and cycle duration. The monitoring of cleaning performance is essential for guaranteeing compliance of cleaning standards.

Product description

Chemdye® Lumenia SixFlow LSF1 Test Kit consists of a challenge device (Chemdye® Lumenia SixFlow CDDER61) and cleaning indicators (Chemdye® Lumenia CDWE) for monitoring cleaning performance in Automated Endoscope Reprocessors. The challenge device consists of three polytetrafluoroethylene (PTFE) tubes with an inner diameter of 1 mm, two PTFE tubes with an inner diameter of 2 mm, and one PTFE tube with an inner diameter of 4 mm; all of them 1.5 m long. Each tube has special connectors in one of their ends for plugging it to the Endoscope Reprocessor, and a device in the other end that simulates an endoscope's lumen and where two Chemdye® Lumenia CDWE cleaning indicators can be placed (one inside and one outside). Each of the six devices have different colors that match the channels' diameters to be monitored, and are in line with the color scheme that every Reprocessor Machine has in its output terminals. Each PTFE tube is also identified with a color mark at the opposite end of the one that is connected to the device, allowing a fast recognition of the output terminal to which it will be connected. The lumens are joined and protected through a hose that covers them entirely. Chemdye® Lumenia CDWE cleaning indicator was especially designed for monitoring these processes. It consists of a synthetic, self-adhesive and thermostable support of 38 mm x 6 mm, printed with a test soil made up from a mixture of especially combined colored natural components. Chemdye® Lumenia CDWE cleaning indicator is manufactured, calibrated and designed in compliance with ISO/TS 15883-5:2005 standard. The cleaning indicator's formulation allows the detection of failures that may affect the cleaning performance of endoscopes (both internally and externally) in Automated Endoscope Reprocessors. The test soil must be completely washed in the internal (A) as well as in the external (B) side of the device (Fig. 1). It is recommended to use both indicators, although optionally, you can verify only the internal cleaning.

Precautions

- Do not reuse cleaning indicators.
- Do not touch the indicator ink with your hands.

Instructions for use

Cleaning performance test

1. Check that Chemdye® Lumenia SixFlow LSF1 Device is clean, dry and well-functioning.
2. Check that your hands are clean and dry.
3. Unscrew the six devices' caps.
4. Place a Chemdye® Lumenia CDWE cleaning indicator in the inner cavity of each device in a way that the ink is facing upwards (Fig.2). **WARNING:** Before placing the indicator in the device, remove the self-adhesive paper. Make sure that the indicator is firmly adhered to the device by using the provided tweezers (Fig.3).
5. Screw the six caps again, making sure that the colors of the caps and the devices match.
6. Place a Chemdye® Lumenia CDWE cleaning indicator in each device's external slot in a way that the ink is visible through its window (Fig.4). **WARNING:** Before placing the indicator in the slot, unlock the device (Fig. 4) and remove the self-adhesive paper. Make sure that the indicator is firmly adhered to the slot by using the provided tweezers. Secure the indicator with the device's lock (Fig. 4).
7. Connect the 6 PTFE tubes of the Lumenia SixFlow LSF1 to the Automated Endoscope Reprocessor through the output terminals.
8. Run a normal cleaning cycle and cancel it after the cleaning phase and just before the disinfection phase. If it is not possible to interrupt the cycle, check the results when the complete cycle finishes.
9. Once the cycle is completed, unplug the connectors and remove the Lumenia SixFlow LSF1 from the washing machine.
10. Carefully remove the indicators from the devices, by using the given tweezers.
11. Examine (by color device) all the indicators, in order to find remains of test soil by placing them on a white background. If you have the Trazanto traceability software, place the indicator in the scanner and proceed with the results recording, considering the reference color for each device. To facilitate the simultaneous recording of the indicators, the Trazanto card can be used.

To position the indicators on the card, only the following two configurations are allowed (Fig. 5):

- 6 internal indicators placed in the first column + 6 external indicators placed in the second column.
- 6 internal indicators placed in the first column + empty second column

(no external indicators placed in the second column).

The Trazanto card does NOT recognize the following configurations:

- First recording column incomplete or empty.
 - Second logging column incomplete (use only full or empty).
- The Trazanto log and analysis system will return 'rejected' if at least one of the indicators on the card is unprocessed or failed.
12. If any of the indicators corresponding to a particular device color has test soil remains, take note of the indicators that have failed and their corresponding color, and take the applicable measures in case of failure of the cleaning performance test.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Obstruction test (channel permeability)

1. Check that Chemdye® Lumenia SixFlow LSF1 Device is clean, dry and well-functioning.
2. Unscrew the six devices' caps and reserve them.
3. Screw a grey cap for obstruction test on each device (grey caps with a hole of 0.5 mm diameter).
4. Connect the 6 PTFE tubes of the Lumenia SixFlow LSF1 to the Automated Endoscope Reprocessor through the output terminals.
5. Run a normal cleaning cycle.
6. The machine should detect an obstruction and trigger an alarm. An obstruction should be detected on each lumen. Record the results.
7. Once the cycle is finished, unplug the device and remove it from the machine.
8. Remove the grey obstruction caps.
9. If no obstructions are detected, take note of the channel number corresponding to the end of the lumen connected to the machine and launch an investigation.

Interpretation of results

If the Chemdye® Lumenia CDWE cleaning indicator is completely washed, it means that the cycle conditions are appropriate for guaranteeing a correct endoscope's cleaning (both internally and externally). Therefore, the flexible endoscope can be reprocessed by using the recently validated program.

The failures on a cleaning cycle may result from an inadequate concentration of chemical agents, insufficient water pressure, inadequate water quality, wrong temperature or duration of the cycle, among others. These facts decrease the cleaning performance, therefore it is not recommended to reprocess or use the endoscope until a satisfactory result is obtained.

NOTE: Chemdye® Lumenia SixFlow LSF1 Kit allows the monitoring of endoscopes' cleaning processes, but it does not make it possible to determine if the appropriate disinfection level has been reached. If your equipment does also perform a high-level disinfection by using peracetic acid as disinfecting agent, you can complete the monitoring through the use of Integrone® IT400 and IT401 chemical indicators, as well as through Bionova® BT400 biological indicator.

Additional technical information

You should verify that the connectors and the PTFE tubes are in perfect conditions in each use. In case the connector is damaged and its replacement is justified, please contact a representative and request a spare part. In case the tube is blocked or it has water leaks, you should replace the device with a new one.

WARNING: Terragene® does not guarantee the well-functioning of Chemdye® Lumenia SixFlow LSF1 Kit if other manufacturers' cleaning indicators or challenge devices are used. Chemdye® Lumenia SixFlow LSF1 device's results are only valid for a cleaning cycle tested with these particular parameters. Do not draw conclusions for cycles with different parameters. Terragene® does not guarantee the well-functioning of Chemdye® Lumenia SixFlow LSF1 Kit if it is used for a cleaning process different from the one it was originally designed for.

Storage

Store the Chemdye® Lumenia SixFlow LSF1 Kit in a dry place, protected from light, at a temperature from 10-30 °C and 30-80 % relative humidity. Do not use Chemdye® Lumenia CDWE indicators after their expiration date. **NOTE:** An incorrect storage of the product may affect its stability as well as the obtained results.

Disposal

Discard the kit components after their use, according to your country's healthcare and safety regulations.

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