

Hygiene Monitoring System

Instructions for use


Remove swab and add the moisturizer


Swab the desired surface by applying a strong pressure


Place the swab back into the device


Activate


Shake downwards vigorously until the solution turns to green and reaches the readout cone reference volume


Shake for 15 seconds with the swab inside the readout cone


Slide swab upwards without completely removing


Again, shake downwards vigorously and reach the readout cone reference volume


Readout cone without swab

Incubate immediately and read

Results reference guide

NOTE: Reference colors exhibited in printed prospectus and boxes as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

Explanation of Symbols

Batch number. Expiration Date. Manufacture Date. Manufacturer.

EN Hygiene Monitoring System

Designed for the detection of residual food-related allergen proteins after cleaning processes

Product Description

The Chemdye® PRO1 ALERT Hygiene Monitoring System is designed for the fast and simple cleaning efficiency verification in the Food Industry. This device can be used to detect the presence of residual proteins deposits after cleaning processes. The system is used along with Bionova® MiniPRO Auto-reader incubators for analytical analysis.

The validation of this system reduces the use of multiple and time-consuming specific allergen tests. Its high sensitivity allows to alert the presence of residual proteins from different sources like milk, nuts, gluten, soy, egg, fish and seafood. Certain food-related protein residues can cause affliction to consumers due to intolerances or allergies. In the food industry, retail and gastronomic sector, ineffective cleaning processes can leave food-related residues over the contacted surfaces; this is highly unsafe for consumers with certain food intolerance or allergy. Since the exact amounts of food-related proteins necessary to elicit a hypersensitive reaction are unknown, a precise and sensitive protein detection and quantification method like the Chemdye® PRO1 ALERT is essential for the Food Industry. It reduces the possibilities of food cross-contact, allows the safe labeling of products and guarantees the safety of consumers. Furthermore, surface soil can protect virus and bacteria from disinfection or sterilization processes and allow their survival. In this regard, the pre-cleaning step of instruments or equipment is critical for the favorable outcome of these processes. Proteins are the greatest challenge to the cleaning step due to its

high surface adhesion properties, especially when they undergo a denaturation process, for example in a thermo-disinfection or cooking process. Therefore, residual protein level is an ideal indicator to reduce contamination.

The Chemdye® PRO1 ALERT system includes a high absorption swab for the collection of samples from different types of surfaces with high efficiency. A visual readout of a color change indicates the presence of detectable levels of protein. Chemdye® PRO1 ALERT has a high sensitivity with a limit of detection (LOD) of 0.5 µg and a quantization range from 1.0 to 50 µg of protein.

Characteristics

Chemdye® PRO1 ALERT along with Bionova® MiniPRO Auto-readers provide all the necessary equipment for cleaning efficiency verification in accordance with the different international guidelines for maintenance and hygiene for the food industry.

Indications for use

Use for testing any surface after the cleaning process and for hard-to-reach areas like irregular surfaces or box joints and edges.

Precautions

Do not open until use. Do not reuse. Do not touch the swab or the inside of the device with your hands. Wear gloves, surgical mask and other hygiene measures. Do not freeze.

High concentrations of alkaline detergents may cause false negative results in some cases. Hydrogen peroxide-based disinfectants may cause a color change reaction giving false positive results.

Although Chemdye® PRO1 ALERT has high sensitivity, it is not a microbiological monitoring device, and the surfaces should not be considered sterile according to the results obtained.

Instructions for use

1. Remove the device from the pouch and allow it to reach room temperature (20-25 °C).

2. Remove the swab from the device and carefully apply 2 drops of moisturizer on the swab (Fig. A).

3. Swab thoroughly the surface to be tested. Press the swab and rotate it several times, while collecting the sample. To validate the system, standardize the surface for sampling and swab the surface in zigzag in one direction, and then do the same procedure perpendicularly (Fig. B).

4. Place the swab back into the device (Fig. C).

5. Activate it, by pressing it firmly downwards (Fig. D).

6. Shake vigorously downwards, until the solution turns to green and reaches the reference volume of the readout cone (Fig. E).

7. Shake vigorously downwards for 15 seconds, with the swab submerged in the solution (Fig. F).

8. Slide the swab upwards without removing it from the device (Fig. G).

9. Shake vigorously downwards and collect the solution until reaching the volume level indicated in the reference readout cone (Fig. H and I).

CAUTION: If the device is not shaken properly (so as the solution in the internal compartment falls completely), or if there is a great presence of bubbles in the solution, the results can be invalid.

10. Make sure the swab is outside the readout cone (Fig. I).

11. For a quantitative analysis, use Bionova® MiniPro Auto-readers incubators. Read the instructions for use before using these devices.

IMPORTANT: The incubation must be done without the swab submerged in the solution and immediately after the activation stage.

IMPORTANT: During the activation of PRO1 ALERT Hygiene Monitoring System, the swab tip could become detached and retained in the reading cone. If this happens, a visual estimation can be carried out after incubation for 4 minutes in the Auto-readers using the reference color guide or repeat the test using a new pen.

IMPORTANT: Color estimations should be done within 5 minutes of removing the device out of the incubator. Results analyzed beyond this period are not valid.

12. After using the Chemdye® PRO1 ALERT System, reprocess the surface used for protein monitoring.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Storage

Chemdye® PRO1 ALERT System should be stored at temperatures between 2-30 °C, away from light and at relative humidity (RH) of 30-80 %. Store vertically as indicated on the packaging. Shelf life is 24 months when stored refrigerated (2-8 °C). Maximum shelf life at room temperature is 6 months (up to 30 °C). 1 month of storage at room temperature reduces the expiration date by 4 months.

Disposal

Discard the product after use according to your country's safety regulations.

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