

Once Chemdye® PRO1 ENDO reaches a temperature of 20-25 °C, remove the swab from the packaging and moisten it by immersing it in the moisturizer solution for 5 seconds.


⚠ Avoid touching the absorbent part of the swab or the inside of the device with your hands. Use gloves, a mask, and other hygiene measures.

Insert the swab into the internal channel of the endoscope or cannulated instrument.


Once the swab has traveled through the entire instrument, cut off the tip of the swab containing the absorbent part (approximately 20 cm from the end) using scissors.


Activate the Chemdye® PRO1 ENDO by pressing the handle all the way down. You should hear two clicks confirming that the activation is successful.

⚠ Once activated, proceed immediately with the following steps.


Collect the solution from the inner compartment by shaking vigorously downward until the solution turns green.


⚠ If the device is not shaken properly to achieve complete lowering of the liquid from the inner compartment or large amounts of bubbles are observed in the solution, the results may be invalid.


Remove the handle from the device.

Gently place the swab for 10 seconds in the solution.


Slide the swab back up and remove it from the device.

Place the handle and immediately incubate it in the Bionova® Auto-readers.

⚠ The pen must be immediately incubated for final reading.
⚠ Do not incubate a solution that turns instantly dark purple after shaking. This may indicate excess proteins that could invalidate the quantification.

⚠ Incubation should be performed with the swab out of the solution (reading cone) and immediately after activation.

After completing the reading, pull the swab remaining inside the cannulated instrument from the distal end to fully remove it.

⚠ After using the Chemdye® PRO1 ENDO system, it is recommended to reprocess the cannulated instrument.


Quantitative analysis

You can use the Bionova® IC10/20FR, IC10/20FRLCD and MiniPro auto-reader incubators with automatic reading systems.

- ⚠ Ensure the proper placement of the pen in the incubator; it should be inserted up to the bottom of the position and should not be able to rotate.
- ⚠ Incubation should be performed with the swab out of the solution (reading cone) and immediately after activation.


Qualitative analysis

You can also use the incubator Bionova® IC10/20.
⚠ Visual estimation should be done within 5 minutes of removing the device from the incubator. Results analyzed beyond this period are not valid.

Results reference guide

NOTE: reference colors exhibited in printed prospects and boxes as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

Explanation of Symbols

- MD** Medical Device.
- CE** CE mark.
- EC REP** Authorized representative in the European Community.
- LOT** Batch number.
- MFG** Manufacture date.
- EXP** Expiration date.
- MAN** Manufacturer.

EN Hygiene Monitoring System

For protein detection and quantification on cannulated instruments

Product Description
The first step in proper cleaning of surgical instruments is to remove blood, body fluids, and tissues immediately after use. If cleaning is not performed properly, the disinfection/sterilization steps will not be effective, leading to the exposure of patients to contaminants such as body fluids, tissues from other patients, and microorganisms.
The Chemdye® PRO1 ENDO Hygiene Monitoring System is designed to verify the cleanliness of cannulated instruments by detecting protein residues that may remain after the cleaning process. The system is used in conjunction with Chemdye® SWE high-absorption swabs, specially designed for cannulated instruments. These swabs must be inserted through the internal channel of the cannulae to sample for contaminants that may remain after the cleaning process.
The Chemdye® PRO1 ENDO system provides quantitative results in the range of 1 to 50 µg of protein with a limit of detection (LOD) of 0.5 µg*, when incubated in the Bionova® MiniPro auto-reader.
*Refers to µg BSA (bovine serum albumin).

Precautions
Handle with gloves at all times. Do not open until use. Do not reuse. Do not re-incubate. Do not swab surfaces that have not been subjected to a cleaning process or visibly contaminated areas. Avoid touching the absorbent part of the swab or the inside of the device with your hands. Use

gloves, a mask, and other hygiene measures. Do not freeze. Do not use Chemdye® SWE swabs on any device other than the Chemdye® PRO1 ENDO. High concentrations of alkaline detergents may cause false-negative results in some cases. Oxidizing disinfectants (such as hydrogen peroxide or peracetic acid) may cause a false-positive color change reaction.

Although Chemdye® PRO1 ENDO is highly sensitive, its intended use corresponds to channel monitoring after the cleaning process, so the test result is not evidence that the device is sterile.

Instructions for use

1. Remove the Chemdye® PRO1 ENDO device from the pouch ten minutes before use, or until it reaches a temperature of 20-25 °C.
2. Remove the Chemdye® SWE swab of the diameter that best suits your needs (see table of suggested diameters) from the pouch, making sure not to touch the absorbent part.
3. Moisten the swab by immersing it for 5 seconds in the wetting solution. Do not shake or press the swab against the sides of the tube (Fig. A).
4. Insert and pass the swab through the internal channel of the endoscope or cannulated instrument by pushing it in one go (Fig. B). Then, use scissors to cut the swab approximately 20 cm from the absorbent end. Be careful not to touch or damage the absorbent part of the swab (Fig. C).
5. Activate the Chemdye® PRO1 ENDO pen by pushing the handle until it stops (Fig. D).

CAUTION: once activated, proceed immediately with the following steps.
6. Shake vigorously downwards until the solution turns green and reaches the reference volume of the reading cone (Fig. E).

CAUTION: if the device is not shaken properly to achieve complete lowering of the liquid from the inner compartment or large amounts of bubbles are observed in the solution, the results may be invalid.

7. Remove the handle from the device (Fig. F).
8. Insert the swab into the device until it is fully immersed in the solution within the reading cone. Gently hold the swab in the solution for 10 seconds (Fig. G).
9. Slide the swab back up and remove it from the device (Fig. H). Place the handle (Fig. I).

CAUTION: the pen must be immediately incubated for final reading.

CAUTION: do not incubate a solution that turns instantly dark purple after shaking. This may indicate excess proteins that could invalidate the quantification.

10. For quantitative analysis you can use the Bionova® IC10/20FR, IC10/20FRLCD and MiniPro auto-reader incubators with automatic reading systems. Read the instructions for use of the equipment beforehand. Incubate the Chemdye® PRO1 ENDO device at 60 ± 2 °C, for:

- 4 minutes in Bionova® MiniPro incubator.
- 7 minutes in Bionova® IC10/20FR and IC10/20FRLCD incubators.

CAUTION: ensure the proper placement of the pen in the incubator; it should be inserted up to the bottom of the position and should not be able to rotate.

CAUTION: incubation should be performed with the swab out of the solution (reading cone) and immediately after activation.

11. For qualitative analysis, you can also use the incubator Bionova® IC10/20. Read the instructions for use of the equipment beforehand. After the incubation time is complete, remove the product from the incubator and interpret it visually using the reference color guide.

• 7 minutes in Bionova® IC10/20.

CAUTION: visual estimation should be done within 5 minutes of removing the device from the incubator. Results analyzed beyond this period are not valid.

12. Pull the cable from the distal end to remove the remaining segment of the swab from the channel (Fig. J).

13. After using the Chemdye® PRO1 ENDO system, it is recommended to reprocess the endoscope or cannulated instrument.

NOTE: in the event of any serious incident involving the product, the incident must be reported to Terragene S.A. and to the competent authority in the State where the user is located.

Storage

The Chemdye® PRO1 ENDO system should be stored at temperatures between 2-8 °C, protected from light and relative humidity (RH) of 30-80 %. Store vertically as indicated on the packaging. Its shelf life is 24 months when refrigerated (2-8 °C). It can also be stored at room temperature (25-30 °C); however, its shelf life is reduced to 6 months in this case.

Waste treatment

Discard the product after use according to your country's health regulations.

Table of suggested swab diameters according to the internal diameter of the cannula:

Cannula internal diameter (mm)	Chemdye® SWE swab diameter (mm)
2,4	1,7
2,8	2
3,5	2,7
4,5	3

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