

# Hygiene Monitoring System

## Instructions for use


Remove swab and add the moisturizer


Swab the desired surface by applying a strong pressure


Place the swab back into the device


Activate


Shake downwards vigorously for 15 seconds. The solution turns to green and reaches the readout cone reference volume

Incubate for 10 minutes at room temperature

## Results reference guide

NOTE: Reference colors exhibited in printed prospectus and boxes as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

### Explanation of Symbols

- CE mark
- Authorized representative in the European Community
- Batch number
- Manufacture Date
- Expiration Date
- Manufacturer

## EN Hygiene Monitoring System

Designed for the detection of residual surface proteins

### Product Description

Surface cleaning and disinfection are mandatory processes applied on a regular basis in medical centers, bio-pharmaceutical industries, food industries, and retail and gastronomy sectors. Protein is the soil component that represents the greatest challenge to cleaning processes due to its high surface adhesion, especially after a denaturation process e.g. during thermo-disinfection or cooking steps. Thus, the residual protein level on surfaces and rinsing water is an appropriate indicator of the efficiency of the cleaning processes. The Chemdye® PRO1 RT Hygiene Monitoring System was designed for the fast and simple protein level determination on instrument's surfaces, hard-to-access areas or in the rinsing water of instruments and equipment after cleaning procedures. The system includes a high absorption swab, which allows the collection of samples from different types of surfaces or solutions with high efficiency. A visual readout of a color change indicates the presence of detectable levels of protein in only 10 minutes allowing to take fast corrective decisions in the presence of contaminants. This system is highly sensitive detecting 15 µg of surface protein or 0.3 mg/ml of protein from rinsing water.

### Characteristics

Chemdye® PRO1 RT provides a fast and simple method for residual protein detection. Allows the verification of cleaning processes in any workplace, it requires minimal operator capabilities and do not require any special instrumentation for the result interpretation.

### Indications for use

Use for testing any surface, hard-to-reach areas in complex instruments, such as irregular surfaces or box joints and edges or rinsing water after the cleaning process.

### Precautions

Do not open until use. Do not reuse. Do not touch the swab or the inside of the device with your hands. Wear gloves, surgical mask and other hygiene measures. Do not freeze. High concentrations of alkaline detergents may cause false negative results in some cases. Hydrogen peroxide-based disinfectants may cause a false positive color change reaction. Although Chemdye® PRO1 RT has high sensitivity, it is not a microbiological monitoring device, and the surfaces should not be considered sterile according to the results obtained.

### Instructions for use

- Remove the device from the pouch and allow it to reach room temperature (20-30 °C).
- Remove the swab from the device.
  - For surface protein detection, carefully apply 2 drops of moisturizer on the swab (Fig. A). On the surface to take the sample, swab thoroughly. Press the swab and rotate it several times, while collecting the sample. To validate the system standardize the surface for sampling and swab the

surface in zigzag in one direction, and then do the same procedure perpendicularly (Fig. B).

b. For solutions, submerge the tip of the swab for at least 5 seconds to allow complete liquid absorption by the swab.

3. Place the swab back into the device (Fig. C).

4. Activate it, by pressing it firmly downwards (Fig. D).

5. Shake vigorously downwards for 15 seconds with the swab submerged in the solution (Fig. E). The solution must turn to green initially in the readout cone.

**WARNING:** Once activated, the incubation time starts running!

6. Incubate for 10 minutes at room temperature (Fig. F).

**IMPORTANT:** The incubation must be done with the swab submerged in the solution and immediately after the activation stage.

The activated device must remain in a vertical position over the incubation time to allow contact of the swab with the reactive solution.

7. Interpretation of the results: After 10 minutes interpret visually by using the reference guide available inside the product's box. **IMPORTANT:** The interpretation must be done immediately. For later interpretation, the reaction can be stopped for 15 minutes by incubation in iced water.

**NOTE:** The time required for a perceptible color change can also be used to estimate protein amounts. This allows optimization of monitoring times according to the following table.

Time registered	Surface protein	Protein in solution
1 minute or less	> 100 µg	> 2 mg/ml
2 minutes or less	> 50 µg	> 1 mg/ml

8. After using the Chemdye® PRO1 RT System, reprocess the surface or instrument used for protein monitoring.

**NOTE:** If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

### Storage

Chemdye® PRO1 RT System should be stored at temperatures between 2-30 °C, away from light and at relative humidity (RH) of 30-80 %. Store vertically as indicated on the packaging. The shelf life is of 24 months when stored refrigerated (between 2-8 °C). Maximum shelf life without refrigeration (up to 30 °C) is 6 months. 1 month of not refrigerated storage reduces the expiration date by 4 months.

### Disposal

Discard the product after use according to your country's safety regulations.

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